

Solutions for **Food Safety**

Application Notebook



Pesticides

Shimadzu Pesticide MRM Library Support for LC/MS/MS

The Library contains information that can be used to accelerate method development in LC/MS/MS pesticide analysis

Expanding Capabilities in Multi-Residue Pesticide Analysis Using The LCMS-8060

The expanded capability to help accelerate method development workflows and support increased pesticide monitoring programs

Applying 'MRM Spectrum Mode' and Library Searching for Enhanced Reporting Confidence in Routine Pesticide Residue Analysis

It was applied to quantify and identify 193 pesticides using 1,291 MRM transitions without compromising limits of detection, linearity or repeatability

LC/MS/MS Method Package for Residual Pesticides Ver. 2 [Brochure]

Pesticide Screening Analysis for 646 Pesticides in 10.5 Minutes

LabSolutions Insight Library Screening [Brochure]

Software helps with screening target components

Using the Nexera UC Online SFE-SFC-MS System to Analyze Residual Pesticides in Agricultural Products

An example of using the Nexera UC online SFE-SFC-MS system to analyze residual pesticides in agricultural products

Analysis of Organophosphorus Pesticides Using Nexis GC-2030

An analysis of organophosphorus pesticides using Nexis GC-2030 equipped with the FPD-2030 was reported

Veterinary Drugs

Quantitative Analysis of Veterinary Drugs Using the Shimadzu LCMS-8050 Triple Quadrupole Mass Spectrometer

An example of the high-sensitivity analysis of 89 veterinary drugs in a crude extract of livestock and fishery products

Multi-Residue Veterinary Drug Analysis of >200 Compounds using MRM Spectrum Mode by LC-MS/MS

In this research, use of the MRM Spectrum mode was examined by quantifying and identifying 212 veterinary drugs

Natural toxins

Analysis of Nivalenol, Deoxynivalenol, 3-Acetyldeoxynivalenol and 15-Acetyldeoxynivalenol Using Triple Quadrupole LC/MS/MS (LCMS-8050)

The test methods specified for deoxynivalenol are HPLC for both qualitative and quantitative analysis, and LC/MS for verification testing

Multi-Residue Analysis of 18 Regulated Mycotoxins by LC/MS/MS

A single LC/MS/MS method was developed and limits of quantification were at or below the maximum levels set in the EC/1886/2006 document

Assay of Aflatoxin M₁ in Milk Based on Notification Test Methodology, Using Prominence-i and the RF-20A_{xs} Fluorescence Detector

High sensitive analysis method of AFM₁ in milk by HPLC with fluorescence detection was evaluated

Analysis of Mycotoxins in Grain Using Mycotoxin Screening System

The screening analysis for mycotoxins in grain products using the i-Series Solution Package mycotoxin screening system

Mycotoxin Screening System (i-Series Solution Package) [Brochure]

EU criteria concentrations of 10 mycotoxin components detected with high sensitivity in only 14 Minutes

Analysis of Diarrhetic Shellfish Toxin Using Triple Quadrupole LC/MS/MS (LCMS-8050)

Instrumental analysis of shellfish by LC/MS/MS offers high sensitivity and accuracy, making it a highly effective analytical method

Antimicrobials

Analysis of Residual Antimicrobials in Meat with Antimicrobial Screening System (Part 1)

An example screening analysis targeting 12 widely used quinolones (old quinolones, new quinolones)

Analysis of Residual Antimicrobials in Meat with Antimicrobial Screening System (Part 2)

An example screening analysis of 12 antimicrobial target compounds including sulfanomides

Antimicrobial Screening System (i-Series Solution Package) [Brochure]

Screening of 24 synthetic antimicrobial compounds that remain in meat

Fertilizers

Analysis of Sulfamic Acid in Fertilizers Using LC/MS (LCMS-2020)

Good quantitative results were obtained, confirming the applicability of this method using byproduct compound fertilizer as the actual sample

Analysis of Nitrous Acid and Ammonium Thiocyanate in Fertilizers

An example of simultaneous analysis of the nitrous acid and ammonium thiocyanate content of fertilizer by HPLC

Analysis of Melamine and Its Related Substances in Fertilizers

An example of pretreating and analyzing melamine and its related substances in fertilizer by HPLC

Insecticides

Ultra-Sensitive and Rapid Assay of Neonicotinoids, Fipronil and Some Metabolites in Honey by UHPLC-MS/MS (LCMS-8060)

A method for ultra-sensitive assay of neonicotinoids in honey was set up

Sensitive method for the determination of Fipronil and its metabolite Fipronil Sulfone in egg using QuEChERS sample pretreatment and LC-MS/MS detection [LCMS-8060]

High sensitive LC-MS/MS method with simplified QuEChERS sample pretreatment was developed for Fipronil determination in egg

Allergens

High Sensitivity Analysis of Peanut Allergen in Cumin and Spice Mix (LCMS-8060)

A method for the analysis of peanut allergen Ara h1 in spices and seasonings was successfully developed

Preservatives

Fast and High Sensitivity Analysis of Six Preservatives in Beverages by UHPLC with Photodiode Array Detection

A rapid and high sensitivity UHPLC method for quantitation of six preservatives in beverages was established

Bisphenols

Development of An UHPLC Method for Simultaneous Determination of Thirteen Bisphenols in Milk Samples

An UHPLC system with a high sensitivity fluorescence detector was adopted to develop a fast and high sensitivity method to meet the requirements of regulations

Application News

No. C135

Liquid Chromatography Mass Spectrometry

Shimadzu Pesticide MRM Library Support for LC/MS/MS

 David R. Baker, Alan Barnes, Neil Loftus
 Shimadzu Corporation, UK

■ Abstract

To help expand capabilities in LC/MS/MS pesticide monitoring programs we have created the Shimadzu Pesticide MRM Library. The Library has been created with 766 certified reference standards and has been verified for use with Shimadzu LCMS-8050 and 8060 systems.

The Library contains information that can be used to accelerate method development in LC/MS/MS pesticide analysis including;

An average of 8 MRM transitions for each reference standard (with optimized collision energies) are registered in the database including positive and negative ionization mode. In total, more than 6,000 MRM transitions are part of the Library.

Meta-data for each library entry such as CAS#, formula, activity, mono-isotopic mass and adduct masses, rank of MRM transitions, synonyms, InChI, InChIKey, compound names translation (Japanese and Chinese) and links to websites offering further information (alanwood.net, PAN pesticide database, Chemical Book, ChemSpider). The meta-data is intended not only to set up new methods quickly but to help search for compound properties.

Key words; Pesticide MRM Library,
766 compound library

■ Using the Shimadzu Pesticide MRM Library

Expanding pesticide monitoring programmes (or creating focused methods) can be quickly set up using the Library data base (Table 1) and create fully optimized MRM methods for LC/MS/MS analysis.

Users select the target pesticides and corresponding transitions from the Library and simply copy the list into a Shimadzu LabSolutions analytical method. The method will include optimized MRM transitions. Once the acquisition method is created users can start to acquire data for screening or quantitative LC/MS/MS analysis.

Table 1 The Shimadzu Pesticide MRM Library supports a list of over 766 compounds. Designed to build extended LC/MS/MS methods quickly and to review pesticide information easily.

Library entries	
Compound information	Compound Name Synonyms Japanese name Chinese name CAS Chemical Formula Mono-isotopic mass Theoretical m/z ([M+H] ⁺ , [M+Na] ⁺ , [M+K] ⁺ , [M+NH ₄] ⁺ , [M-H] ⁻) Activity InChI InChIKey
MS/MS parameters	Ionization mode Q1 (m/z) Q3 (m/z) Q1 Pre Bias CE Q3 Pre Bias
Web links	Alanwood.net PAN Pesticide Database Chemical Book ChemSpider

	Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
1	(E)-Fenpyroximate	134098-61-6	C24H27N3O4	421.2002	422.2075	420.1929	ESI+	6
2	(E)-Ferimzone	89269-64-7	C15H18N4	254.1531	255.1604	253.1458	ESI+	2
3	(Z)-Fenpyroximate	149054-53-5	C24H27N3O4	421.2002	422.2075	420.1929	ESI+	2
4	(Z)-Ferimzone	89269-64-7	C15H18N4	254.1531	255.1604	253.1458	ESI+	6
5	1-(3, 4-Dichlorophenyl)-3-methylurea	3567-62-2	C8H8Cl2N2O	218.0014	219.0087	216.9941	ESI+	19
6	1-(3, 4-Dichlorophenyl)urea	2327-02-8	C7H6Cl2N2O	203.9857	204.9930	202.9784	ESI+	17
7	1-(4-Isopropylphenyl)-3-methylurea	34123-57-4	C11H16N2O	192.1263	193.1336	191.1190	ESI+	6
8	1-(4-Isopropylphenyl)urea	56046-17-4	C10H14N2O	178.1106	179.1179	177.1033	ESI+	6
9	1-naphthaleneacetamide	86-86-2	C12H11NO	185.0841	186.0914	184.0768	ESI+	4
10	1-Naphthaleneacetic Acid	86-87-3	C12H10O2	186.0681	187.0754	185.0608	ESI-	1
11	2, 4, 5-T	93-76-5	C8H5Cl3O3	253.9304	254.9377	252.9231	ESI-	7
12	2, 4, 6-Tribromophenol	118-79-6	C6H3Br3O	327.7734	328.7807	326.7661	ESI+	10
13	2, 4, 6-Trichlorophenol	88-06-2	C6H3Cl3O	195.9249	196.9322	194.9176	ESI-	3
14	2, 4-D (2, 4-PA)	94-75-7	C8H6Cl2O3	219.9694	220.9767	218.9621	ESI-	7
15	2, 4-DB	94-82-6	C10H10Cl2O3	248.0007	249.0080	246.9934	ESI-	5
16	2, 4-dimethylaniline	95-68-1	C8H11N	121.0891	122.0964	120.0818	ESI+	5
17	2, 6-Dichlorobenzamide	2008-58-4	C7H5Cl2NO	188.9748	189.9821	187.9675	ESI+	13
18	2-Naphthoxy acetic acid	120-23-0	C12H10O3	202.0630	203.0703	201.0557	ESI-	2
19	2-Phenylphenol	90-43-7	C12H10O	170.0732	171.0805	169.0659	ESI-	2
20	3-(3-Indolyl)-propionic acid	830-96-6	C11H11NO2	189.0790	190.0863	188.0717	ESI+	6
21	3, 4, 5-Trimethacarb	2686-99-9	C11H15NO2	193.1103	194.1176	192.1030	ESI+	12
22	3-Indolyl-acetic acid	87-51-4	C10H9NO2	175.0633	176.0706	174.0560	ESI+	12
23	3-Methylphosphinopropionic acid	15090-23-0	C4H9O4P	152.0238	153.0311	151.0165	ESI+	12
24	4-(3-Indolyl)-butyric acid	133-32-4	C12H13NO2	203.0946	204.1019	202.0873	ESI+	14
25	4-Chlorophenoxyacetic acid	122-88-3	C8H7ClO3	186.0084	187.0157	185.0011	ESI-	4
26	6-chloro-3-phenylpyridazin-4-ol	40020-01-7	C10H7ClN2O	206.0247	207.0320	205.0174	ESI+	6
27	6-Furfurylaminopurine	525-79-1	C10H9N5O	215.0807	216.0880	214.0734	ESI+	9
28	Acephate	30560-19-1	C4H10NO3PS	183.0119	184.0192	182.0046	ESI+	6
29	Acequinocyl	57960-19-7	C24H32O4	384.2301	385.2374	383.2228	ESI+	6
30	Acetamidiprid	135410-20-7	C10H11ClN4	222.0672	223.0745	221.0599	ESI+	10
31	Acibenzolar-S-methyl	135158-54-2	C8H6N2OS2	209.9922	210.9995	208.9849	ESI+	6
32	Acifluorfen	50594-66-6	C14H7ClF3NO5	360.9965	362.0038	359.9892	ESI-	12
33	Acionifen	74070-46-5	C12H9ClN2O3	264.0302	265.0375	263.0229	ESI+	2
34	Acrinathrin	101007-06-1	C26H21F6NO5	541.1324	542.1397	540.1251	ESI+	12
35	Alachlor	15972-60-8	C14H20ClNO2	269.1183	270.1256	268.1110	ESI+	12
36	Alanycarb	83130-01-2	C17H25N3O4S2	399.1286	400.1359	398.1213	ESI+	6
37	Aldicarb	116-06-3	C7H14N2O2S	190.0776	191.0849	189.0703	ESI+	5
38	Aldicarb-sulfone (Aldoxycarb)	1646-88-4	C7H14N2O4S	222.0674	223.0747	221.0601	ESI+	5
39	Aldicarb-sulfoxide	1646-87-3	C7H14N2O3S	206.0725	207.0798	205.0652	ESI+	8
40	Allethrin	584-79-2	C19H26O3	302.1882	303.1955	301.1809	ESI+	12
41	Allidochlor	93-71-0	C8H12ClNO	173.0607	174.0680	172.0534	ESI+	12
42	Ametoctradin	865318-97-4	C15H25N5	275.2110	276.2183	274.2037	ESI+	6
43	Ametryn	834-12-8	C9H17N5S	227.1205	228.1278	226.1132	ESI+	6
44	Amidosulfuron	120923-37-7	C9H15N5O7S2	369.0413	370.0486	368.0340	ESI+	8
45	Aminocarb	2032-59-9	C11H16N2O2	208.1212	209.1285	207.1139	ESI+	6
46	Aminopyralid	150114-71-9	C6H4Cl2N2O2	205.9650	206.9723	204.9577	ESI+	7
47	Amisulbrom	348635-87-0	C13H13BrFN5O4S2	464.9576	465.9649	463.9503	ESI+	10
48	Amitraz	33089-61-1	C19H23N3	293.1892	294.1965	292.1819	ESI+	2
49	Amitrole	61-82-5	C2H4N4	84.0436	85.0509	83.0363	ESI+	5
50	AMPA	1066-51-9	CH6NO3P	111.0085	112.0158	110.0012	ESI-	3
51	Ancymidol	12771-68-5	C15H16N2O2	256.1212	257.1285	255.1139	ESI+	6
52	Anilazine	101-05-3	C9H5Cl3N4	273.9580	274.9653	272.9507	ESI+	12
53	Anilofos	64249-01-0	C13H19ClNO3PS2	367.0232	368.0305	366.0159	ESI+	12
54	Aramite	140-57-8	C15H23ClO4S	334.1006	335.1079	333.0933	ESI+	12
55	Asulam	3337-71-1	C8H10N2O4S	230.0361	231.0434	229.0288	ESI+	9
56	Atraton	1610-17-9	C9H17N5O	211.1433	212.1506	210.1360	ESI+	6
57	Atrazine	1912-24-9	C8H14ClN5	215.0938	216.1011	214.0865	ESI+	8
58	Atrazine-2-hydroxy	2163-68-0	C8H15N5O	197.1277	198.1350	196.1204	ESI+	6
59	Atrazine-desethyl	6190-65-4	C6H10ClN5	187.0625	188.0698	186.0552	ESI+	9
60	Atrazine-desethyl-2-hydroxy	19988-24-0	C6H11N5O	169.0964	170.1037	168.0891	ESI+	5
61	Atrazine-desisopropyl	1007-28-9	C5H8ClN5	173.0468	174.0541	172.0395	ESI+	10
62	Avermectin B1a	65195-55-3	C48H72O14	872.4922	873.4995	871.4849	ESI+	4
63	Avermectin B1b	65195-56-4	C47H70O14	858.4766	859.4839	857.4693	ESI+	3
64	Azaconazole	60207-31-0	C12H11Cl2N3O2	299.0228	300.0301	298.0155	ESI+	8
65	Azadirachtin	11141-17-6	C35H44O16	720.2629	721.2702	719.2556	ESI+	8
66	Azamethiphos	35575-96-3	C9H10ClN2O5PS	323.9737	324.9810	322.9664	ESI+	11
67	Azimsulfuron	120162-55-2	C13H16N10O5S	424.1026	425.1099	423.0953	ESI+	5
68	Azinphos-ethyl	2642-71-9	C12H16N3O3PS2	345.0371	346.0444	344.0298	ESI+	5
69	Azinphos-methyl	86-50-0	C10H12N3O3PS2	317.0058	318.0131	315.9985	ESI+	6
70	Aziprotryne	4658-28-0	C7H11N7S	225.0797	226.0870	224.0724	ESI+	4
71	Azobenzene	103-33-3	C12H10N2	182.0844	183.0917	181.0771	ESI+	2
72	Azoxystrobin	131860-33-8	C22H17N3O5	403.1168	404.1241	402.1095	ESI+	5
73	Barban	101-27-9	C11H9Cl2NO2	257.0010	258.0083	255.9937	ESI+	11
74	Beflubutamid	113614-08-7	C18H17F4NO2	355.1195	356.1268	354.1122	ESI+	10
75	Benalaxyl	71626-11-4	C20H23NO3	325.1678	326.1751	324.1605	ESI+	6

	Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
76	Benazolin	3813-05-6	C9H6ClNO3S	242.9757	243.9830	241.9684	ESI+	6
77	Benazolin-ethyl	25059-80-7	C11H10ClNO3S	271.0070	272.0143	269.9997	ESI+	18
78	Bendiocarb	22781-23-3	C11H13NO4	223.0845	224.0918	222.0772	ESI+	6
79	Benfuracarb	82560-54-1	C20H30N2O5S	410.1875	411.1948	409.1802	ESI+	5
80	Benfuresate	68505-69-1	C12H16O4S	256.0769	257.0842	255.0696	ESI+	2
81	Benodanil	15310-01-7	C13H10INO	322.9807	323.9880	321.9734	ESI+	6
82	Benoxacor	98730-04-2	C11H11Cl2NO2	259.0167	260.0240	258.0094	ESI+	17
83	Bensulfuron-methyl	83055-99-6	C16H18N4O7S	410.0896	411.0969	409.0823	ESI+	6
84	Bensulide	741-58-2	C14H24NO4PS3	397.0605	398.0678	396.0532	ESI+	9
85	Bentazone	25057-89-0	C10H12N2O3S	240.0569	241.0642	239.0496	ESI-	5
86	Benthiavalicarb-isopropyl	177406-68-7	C18H24FN3O3S	381.1522	382.1595	380.1449	ESI+	5
87	Benthiazole	21564-17-0	C9H6N2S3	237.9693	238.9766	236.9620	ESI+	6
88	Benzanilide	93-98-1	C13H11NO	197.0841	198.0914	196.0768	ESI+	4
89	Benzo fenap	82692-44-2	C22H20Cl2N2O3	430.0851	431.0924	429.0778	ESI+	2
90	Benzoximate	29104-30-1	C18H18ClNO5	363.0874	364.0947	362.0801	ESI+	12
91	Benzoylprop-ethyl	22212-55-1	C18H17Cl2NO3	365.0585	366.0658	364.0512	ESI+	6
92	Benzthiazuron	1929-88-0	C9H9N3OS	207.0466	208.0539	206.0393	ESI+	9
93	Benzyl dimethyl dodecyl ammonium	139-07-1	C21H37N	303.2926	304.2999	302.2853	ESI+	4
94	Benzyl dimethyl hexadecyl ammonium	122-18-9	C25H45N	359.3552	360.3625	358.3479	ESI+	3
95	Benzyl dimethyl tetradecyl ammonium	139-08-2	C23H41N	331.3239	332.3312	330.3166	ESI+	3
96	Bifenazate	149877-41-8	C17H20N2O3	300.1474	301.1547	299.1401	ESI+	6
97	Bifenox	42576-02-3	C14H9Cl2NO5	340.9858	341.9931	339.9785	ESI+	8
98	Bifenthrin	82657-04-3	C23H22ClF3O2	422.1260	423.1333	421.1187	ESI+	5
99	Bioresmethrin	28434-01-7	C22H26O3	338.1882	339.1955	337.1809	ESI+	6
100	Bispyribac-sodium	125401-92-5	C19H17N4NaO8	452.0944	453.1017	451.0871	ESI+	8
101	Bitertanol	55179-31-2	C20H23N3O2	337.1790	338.1863	336.1717	ESI+	6
102	Bixafen	581809-46-3	C18H12Cl2F3N3O	413.0310	414.0383	412.0237	ESI+	12
103	Boscalid	188425-85-6	C18H12Cl2N2O	342.0327	343.0400	341.0254	ESI+	12
104	Brodifacoum	56073-10-0	C31H23BrO3	522.0831	523.0904	521.0758	ESI+	12
105	Bromacil	314-40-9	C9H13BrN2O2	260.0160	261.0233	259.0087	ESI+	9
106	Bromadiolone	28772-56-7	C30H23BrO4	526.0780	527.0853	525.0707	ESI-	12
107	Bromfeninfos	33399-00-7	C12H14BrCl2O4P	401.9190	402.9263	400.9117	ESI+	17
108	Bromobutide	74712-19-9	C15H22BrNO	311.0885	312.0958	310.0812	ESI+	10
109	Bromophos-ethyl	4824-78-6	C10H12BrCl2O3PS	391.8805	392.8878	390.8732	ESI+	3
110	Bromophos-methyl	2104-96-3	C8H8BrCl2O3PS	363.8492	364.8565	362.8419	ESI+	6
111	Bromoxynil	1689-84-5	C7H3Br2NO	274.8581	275.8654	273.8508	ESI-	11
112	Bromuconazole	116255-48-2	C13H12BrCl2N3O	374.9541	375.9614	373.9468	ESI+	11
113	Bupirimate	41483-43-6	C13H24N4O3S	316.1569	317.1642	315.1496	ESI+	6
114	Buprofezin	69327-76-0	C16H23N3O5	305.1562	306.1635	304.1489	ESI+	6
115	Butachlor	23184-66-9	C17H26ClNO2	311.1652	312.1725	310.1579	ESI+	12
116	Butafenacil	134605-64-4	C20H18ClF3N2O6	474.0805	475.0878	473.0732	ESI+	10
117	Butamifos	36335-67-8	C13H21N2O4PS	332.0960	333.1033	331.0887	ESI+	12
118	Butocarboxim	34681-10-2	C7H14N2O2S	190.0776	191.0849	189.0703	ESI+	3
119	Butocarboxim-sulfone	34681-23-7	C7H14N2O4S	222.0674	223.0747	221.0601	ESI+	14
120	Butocarboxim-sulfoxide	34681-24-8	C7H14N2O3S	206.0725	207.0798	205.0652	ESI+	6
121	Butralin	33629-47-9	C14H21N3O4	295.1532	296.1605	294.1459	ESI+	6
122	Buturon	3766-60-7	C12H13ClN2O	236.0716	237.0789	235.0643	ESI+	9
123	Butylate	2008-41-5	C11H23NO5	217.1500	218.1573	216.1427	ESI+	3
124	Cadusafos	95465-99-9	C10H23O2PS2	270.0877	271.0950	269.0804	ESI+	5
125	Cafenstrole	125306-83-4	C16H22N4O3S	350.1413	351.1486	349.1340	ESI+	3
126	Captafol	2425-06-1	C10H9Cl4NO2S	346.9108	347.9181	345.9035	ESI+	1
127	Carbaryl (NAC)	63-25-2	C12H11NO2	201.0790	202.0863	200.0717	ESI+	6
128	Carbendazim	10605-21-7	C9H9N3O2	191.0695	192.0768	190.0622	ESI+	5
129	Carbetamide	16118-49-3	C12H16N2O3	236.1161	237.1234	235.1088	ESI+	6
130	Carbofuran	1563-66-2	C12H15NO3	221.1052	222.1125	220.0979	ESI+	6
131	Carbofuran-3-hydroxy (3-Hydroxycarbofuran)	16655-82-6	C12H15NO4	237.1001	238.1074	236.0928	ESI+	12
132	Carbofuran-3-keto	16709-30-1	C12H13NO4	235.0845	236.0918	234.0772	ESI+	12
133	Carbophenothion	786-19-6	C11H16ClO2PS3	341.9739	342.9812	340.9666	ESI+	9
134	Carbosulfan	55285-14-8	C20H32N2O3S	380.2134	381.2207	379.2061	ESI+	6
135	Carboxin	5234-68-4	C12H13NO2S	235.0667	236.0740	234.0594	ESI+	6
136	Carfentrazone-ethyl	128639-02-1	C15H14Cl2F3N3O3	411.0364	412.0437	410.0291	ESI+	5
137	Carpropamid	104030-54-8	C15H18Cl3NO	333.0454	334.0527	332.0381	ESI+	18
138	Cartap	15263-53-3	C7H15N3O2S2	237.0606	238.0679	236.0533	ESI+	3
139	Chinomethionat	2439-01-2	C10H6N2OS2	233.9922	234.9995	232.9849	ESI+	6
140	Chloramphenicol	56-75-7	C11H12Cl2N2O5	322.0123	323.0196	321.0050	ESI-	17
141	Chlorantraniliprole	500008-45-7	C18H14BrCl2N5O2	480.9708	481.9781	479.9635	ESI+	28
142	Chlorbromuron	13360-45-7	C9H10BrClN2O2	291.9614	292.9687	290.9541	ESI+	12
143	Chlorbufam	1967-16-4	C11H10ClNO2	223.0400	224.0473	222.0327	ESI+	4
144	Chlordimeform	6164-98-3	C10H13ClN2	196.0767	197.0840	195.0694	ESI+	12
145	Chlorfenvinphos	470-90-6	C12H14Cl3O4P	357.9695	358.9768	356.9622	ESI+	12
146	Chlorflazuron	71422-67-8	C20H9Cl3F5N3O3	538.9630	539.9703	537.9557	ESI+	17
147	Chloridazon	1698-60-8	C10H8ClN3O	221.0356	222.0429	220.0283	ESI+	11
148	Chlorimuron-ethyl	90982-32-4	C15H15ClN4O6S	414.0401	415.0474	413.0328	ESI+	12
149	Chlormequat-chloride	999-81-5	C5H13Cl2N	157.0425	158.0498	156.0352	ESI+	6
150	Chlorophacinone	3691-35-8	C23H15ClO3	374.0710	375.0783	373.0637	ESI-	15

	Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
151	Chlorotoluron	15545-48-9	C10H13ClN2O	212.0716	213.0789	211.0643	ESI+	8
152	Chloroxuron	1982-47-4	C15H15ClN2O2	290.0822	291.0895	289.0749	ESI+	12
153	Chloroxynil	1891-95-8	C7H3Cl2NO	186.9592	187.9665	185.9519	ESI-	6
154	Chlorpropham	101-21-3	C10H12ClNO2	213.0557	214.0630	212.0484	ESI+	2
155	Chlorpyrifos	2921-88-2	C9H11Cl3NO3PS	348.9263	349.9336	347.9190	ESI+	16
156	Chlorpyrifos-methyl	5598-13-0	C7H7Cl3NO3PS	320.8950	321.9023	319.8877	ESI+	12
157	Chlorpyrifos-oxon	5598-15-2	C9H11Cl3NO4P	332.9491	333.9564	331.9418	ESI+	24
158	Chlorsulfuron	64902-72-3	C12H12ClN5O4S	357.0299	358.0372	356.0226	ESI+	18
159	Chlorthiamid	1918-13-4	C7H5Cl2NS	204.9520	205.9593	203.9447	ESI+	16
160	Chromafenozide	143807-66-3	C24H30N2O3	394.2256	395.2329	393.2183	ESI+	6
161	Cinidon-ethyl	142891-20-1	C19H17Cl2NO4	393.0535	394.0608	392.0462	ESI+	24
162	Cinosulfuron	94593-91-6	C15H19N5O7S	413.1005	414.1078	412.0932	ESI+	6
163	Clethodim	99129-21-2	C17H26ClNO3S	359.1322	360.1395	358.1249	ESI+	10
164	Climbazole	38083-17-9	C15H17ClN2O2	292.0979	293.1052	291.0906	ESI+	9
165	Clodinafop (free acid)	114420-56-3	C14H11ClFNO4	311.0361	312.0434	310.0288	ESI+	8
166	Clodinafop-propargyl	105512-06-9	C17H13ClFNO4	349.0517	350.0590	348.0444	ESI+	12
167	Clofentazine	74115-24-5	C14H8Cl2N4	302.0126	303.0199	301.0053	ESI+	10
168	Clomazone	81777-89-1	C12H14ClNO2	239.0713	240.0786	238.0640	ESI+	8
169	Clomeprop	84496-56-0	C16H15Cl2NO2	323.0480	324.0553	322.0407	ESI+	21
170	Cloprop	101-10-0	C9H9ClO3	200.0240	201.0313	199.0167	ESI-	2
171	Clopyralid	1702-17-6	C6H3Cl2NO2	190.9541	191.9614	189.9468	ESI-	2
172	Cloquintocet-mexyl	99607-70-2	C18H22ClNO3	335.1288	336.1361	334.1215	ESI+	9
173	Cloransulam-methyl	147150-35-4	C15H13ClFN5O5S	429.0310	430.0383	428.0237	ESI+	12
174	Clothianidin	210880-92-5	C6H8ClN5O2S	249.0087	250.0160	248.0014	ESI+	7
175	Coumachlor	81-82-3	C19H15ClO4	342.0659	343.0732	341.0586	ESI+	18
176	Coumaphos	56-72-4	C14H16ClO5PS	362.0145	363.0218	361.0072	ESI+	12
177	Coumatetralyl	5836-29-3	C19H16O3	292.1099	293.1172	291.1026	ESI+	6
178	Crimidine	535-89-7	C7H10ClN3	171.0563	172.0636	170.0490	ESI+	12
179	Crotoxyphos	7700-17-6	C14H19O6P	314.0919	315.0992	313.0846	ESI+	6
180	Cruformate	299-86-5	C12H19ClNO3P	291.0791	292.0864	290.0718	ESI+	12
181	Cumyluron	99485-76-4	C17H19ClN2O	302.1186	303.1259	301.1113	ESI+	2
182	Cyanazine	21725-46-2	C9H13ClN6	240.0890	241.0963	239.0817	ESI+	6
183	Cyanofenphos	13067-93-1	C15H14N4O2PS	303.0483	304.0556	302.0410	ESI+	6
184	Cyazofamid	120116-88-3	C13H13ClN4O2S	324.0448	325.0521	323.0375	ESI+	5
185	Cyclanilide	113136-77-9	C11H9Cl2NO3	272.9959	274.0032	271.9886	ESI-	20
186	Cycloate	1134-23-2	C11H21NOS	215.1344	216.1417	214.1271	ESI+	5
187	Cycloheximide	66-81-9	C15H23NO4	281.1627	282.1700	280.1554	ESI+	12
188	Cycloprothrin	63935-38-6	C26H21Cl2NO4	481.0848	482.0921	480.0775	ESI+	2
189	Cyclosulfamuron	136849-15-5	C17H19N5O6S	421.1056	422.1129	420.0983	ESI+	6
190	Cycloxydim	101205-02-1	C17H27NO3S	325.1712	326.1785	324.1639	ESI+	10
191	Cycluron	2163-69-1	C11H22N2O	198.1732	199.1805	197.1659	ESI+	5
192	Cyflufenamid	180409-60-3	C20H17F5N2O2	412.1210	413.1283	411.1137	ESI+	6
193	Cyflumetofen	400882-07-7	C24H24F3NO4	447.1657	448.1730	446.1584	ESI+	8
194	Cyhalofop-butyl	122008-85-9	C20H20FNO4	357.1376	358.1449	356.1303	ESI+	3
195	Cymiazole	61676-87-7	C12H14N2S	218.0878	219.0951	217.0805	ESI+	6
196	Cymoxanil	57966-95-7	C7H10N4O3	198.0753	199.0826	197.0680	ESI+	4
197	Cypermethrin	52315-07-8	C22H19Cl2NO3	415.0742	416.0815	414.0669	ESI+	10
198	Cyphenothrin	39515-40-7	C24H25NO3	375.1834	376.1907	374.1761	ESI+	12
199	Cyproconazole	94361-06-5	C15H18ClN3O	291.1138	292.1211	290.1065	ESI+	10
200	Cyprodinil	121552-61-2	C14H15N3	225.1266	226.1339	224.1193	ESI+	6
201	Cyromazine	66215-27-8	C6H10N6	166.0967	167.1040	165.0894	ESI+	6
202	Daimuron (Dymron)	42609-52-9	C17H20N2O	268.1576	269.1649	267.1503	ESI+	6
203	Dalapon	75-99-0	C3H4Cl2O2	141.9588	142.9661	140.9515	ESI-	10
204	Daminozide	1596-84-5	C6H12N2O3	160.0848	161.0921	159.0775	ESI+	6
205	Dazomet	533-74-4	C5H10N2S2	162.0285	163.0358	161.0212	ESI+	6
206	Deet	134-62-3	C12H17NO	191.1310	192.1383	190.1237	ESI+	2
207	Deltamethrin	52918-63-5	C22H19Br2NO3	502.9732	503.9805	501.9659	ESI+	12
208	Demeton-O	298-03-3	C8H19O3PS2	258.0513	259.0586	257.0440	ESI+	2
209	Demeton-S	126-75-0	C8H19O3PS2	258.0513	259.0586	257.0440	ESI+	3
210	Demeton-S-methyl	919-86-8	C6H15O3PS2	230.0200	231.0273	229.0127	ESI+	2
211	Demeton-S-methyl-sulfone	17040-19-6	C6H15O5PS2	262.0099	263.0172	261.0026	ESI+	6
212	Desmedipham	13684-56-5	C16H16N2O4	300.1110	301.1183	299.1037	ESI+	6
213	Desmetryn	1014-69-3	C8H15N5S	213.1048	214.1121	212.0975	ESI+	4
214	Diafenthion	80060-09-9	C23H32N2OS	384.2235	385.2308	383.2162	ESI+	12
215	Dialifos	10311-84-9	C14H17ClNO4PS2	393.0025	394.0098	391.9952	ESI+	12
216	Diallate	2303-16-4	C10H17Cl2NO5	269.0408	270.0481	268.0335	ESI+	12
217	Diazinon	333-41-5	C12H21N2O3PS	304.1010	305.1083	303.0937	ESI+	6
218	Dicamba	1918-00-9	C8H6Cl2O3	219.9694	220.9767	218.9621	ESI-	2
219	Dichlofenthion	97-17-6	C10H13Cl2O3PS	313.9700	314.9773	312.9627	ESI+	8
220	Dichlofluanid	1085-98-9	C9H11Cl2FN2O2S2	331.9623	332.9696	330.9550	ESI+	11
221	Dichlormid	37764-25-3	C8H11Cl2NO	207.0218	208.0291	206.0145	ESI+	19
222	Dichlorprop	120-36-5	C9H8Cl2O3	233.9850	234.9923	232.9777	ESI-	8
223	Dichlorvos	62-73-7	C4H7Cl2O4P	219.9459	220.9532	218.9386	ESI+	17
224	Diclobutrazol	75736-33-3	C15H19Cl2N3O	327.0905	328.0978	326.0832	ESI+	4
225	Diclofop	40843-25-2	C15H12Cl2O4	326.0113	327.0186	325.0040	ESI-	4

	Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
226	Diclofop-methyl	51338-27-3	C16H14Cl2O4	340.0269	341.0342	339.0196	ESI+	12
227	Dicloran	99-30-9	C6H4Cl2N2O2	205.9650	206.9723	204.9577	ESI+	4
228	Diclofulam	145701-21-9	C13H10Cl2FN5O3S	404.9865	405.9938	403.9792	ESI+	9
229	Dicrotophos	141-66-2	C8H16NO5P	237.0766	238.0839	236.0693	ESI+	6
230	Dicyclanil	112636-83-6	C8H10N6	190.0967	191.1040	189.0894	ESI+	6
231	Didecyldimethylammonium	7173-51-5	C22H47N	325.3709	326.3782	324.3636	ESI+	6
232	Diethanolamine	111-42-2	C4H11NO2	105.0790	106.0863	104.0717	ESI+	6
233	Diethofencarb	87130-20-9	C14H21NO4	267.1471	268.1544	266.1398	ESI+	6
234	Difenacoum	56073-07-5	C31H24O3	444.1725	445.1798	443.1652	ESI+	12
235	Difenoconazole	119446-68-3	C19H17Cl2N3O3	405.0647	406.0720	404.0574	ESI+	12
236	Difenoxuron	14214-32-5	C16H18N2O3	286.1317	287.1390	285.1244	ESI+	6
237	Difenzoquat-methyl-sulfate	43222-48-6	C17H16N2	248.1313	249.1386	247.1240	ESI+	6
238	Diffubenzuron	35367-38-5	C14H9ClF2N2O2	310.0321	311.0394	309.0248	ESI+	9
239	Diffufenican	83164-33-4	C19H11F5N2O2	394.0741	395.0814	393.0668	ESI+	12
240	Dimefuron	34205-21-5	C15H19ClN4O3	338.1146	339.1219	337.1073	ESI+	5
241	Dimepiperate	61432-55-1	C15H21NOS	263.1344	264.1417	262.1271	ESI+	6
242	Dimethachlon	24096-53-5	C10H7Cl2NO2	242.9854	243.9927	241.9781	ESI-	2
243	Dimethachlor	50563-36-5	C13H18ClNO2	255.1026	256.1099	254.0953	ESI+	12
244	Dimethametryn	22936-75-0	C11H21N5S	255.1518	256.1591	254.1445	ESI+	6
245	Dimethenamid	87674-68-8	C12H18ClNO2S	275.0747	276.0820	274.0674	ESI+	12
246	Dimethirimol	5221-53-4	C11H19N3O	209.1528	210.1601	208.1455	ESI+	3
247	Dimethoate	60-51-5	C5H12NO3PS2	228.9996	230.0069	227.9923	ESI+	6
248	Dimethomorph	110488-70-5	C21H22ClNO4	387.1237	388.1310	386.1164	ESI+	12
249	Dimetilan	644-64-4	C10H16N4O3	240.1222	241.1295	239.1149	ESI+	6
250	Dimoxystrobin	149961-52-4	C19H22N2O3	326.1630	327.1703	325.1557	ESI+	6
251	Diniconazole	83657-24-3	C15H17Cl2N3O	325.0749	326.0822	324.0676	ESI+	7
252	Dinocap	39300-45-3	C18H24N2O6	364.1634	365.1707	363.1561	ESI+	6
253	Dinoseb	88-85-7	C10H12N2O5	240.0746	241.0819	239.0673	ESI-	4
254	Dinotefuran	165252-70-0	C7H14N4O3	202.1066	203.1139	201.0993	ESI+	6
255	Dinoterb	1420-07-1	C10H12N2O5	240.0746	241.0819	239.0673	ESI-	4
256	Dioxacarb	6988-21-2	C11H13NO4	223.0845	224.0918	222.0772	ESI+	6
257	Dioxathion	78-34-2	C12H26O6P2S4	456.0087	457.0160	455.0014	ESI+	6
258	Diphenamid	957-51-7	C16H17NO	239.1310	240.1383	238.1237	ESI+	6
259	Diphenylamine	122-39-4	C12H11N	169.0891	170.0964	168.0818	ESI+	4
260	Dipropetryn	4147-51-7	C11H21N5S	255.1518	256.1591	254.1445	ESI+	6
261	Diquat	6385-62-2	C12H12N2	184.1000	185.1073	183.0927	ESI+	3
262	Disulfoton	298-04-4	C8H19O2PS3	274.0285	275.0358	273.0212	ESI+	3
263	Disulfoton-sulfone	2497-06-5	C8H19O4PS3	306.0183	307.0256	305.0110	ESI+	6
264	Disulfoton-sulfoxide	2497-07-6	C8H19O3PS3	290.0234	291.0307	289.0161	ESI+	6
265	Ditalimfos	5131-24-8	C12H14NO4PS	299.0381	300.0454	298.0308	ESI+	6
266	Dithianon	3347-22-6	C14H4N2O2S2	295.9714	296.9787	294.9641	ESI-	4
267	Dithiopyr	97886-45-8	C15H16F5NO2S2	401.0543	402.0616	400.0470	ESI+	6
268	Diuron (DCMU)	330-54-1	C9H10Cl2N2O	232.0170	233.0243	231.0097	ESI+	7
269	DMST	66840-71-9	C9H14N2O2S	214.0776	215.0849	213.0703	ESI+	4
270	DNOC	534-52-1	C7H6N2O5	198.0277	199.0350	197.0204	ESI-	6
271	Dodemorph	1593-77-7	C18H35NO	281.2719	282.2792	280.2646	ESI+	6
272	Dodine	2439-10-3	C15H33N3O2	287.2573	288.2646	286.2500	ESI+	6
273	Doramectin	117704-25-3	C50H74O14	898.5079	899.5152	897.5006	ESI+	10
274	Edifenphos	17109-49-8	C14H15O2PS2	310.0251	311.0324	309.0178	ESI+	6
275	Eamectin B1a	119791-41-2	C49H75NO13	885.5238	886.5311	884.5165	ESI+	5
276	Eamectin B1b	137335-79-6	C55H79NO15	871.5082	872.5155	870.5009	ESI+	3
277	Endosulfan-sulfate	1031-07-8	C9H6Cl6O4S	419.8118	420.8191	418.8045	ESI-	3
278	EPN	2104-64-5	C14H14NO4PS	323.0381	324.0454	322.0308	ESI+	6
279	Epoxiconazole	133855-98-8	C17H13ClFN3O	329.0731	330.0804	328.0658	ESI+	9
280	EPTC	759-94-4	C9H19NOS	189.1187	190.1260	188.1114	ESI+	5
281	Esfenvalerate	66230-04-4	C25H22ClNO3	419.1288	420.1361	418.1215	ESI+	2
282	Esprocarb	85785-20-2	C15H23NOS	265.1500	266.1573	264.1427	ESI+	5
283	Etaconazole	60207-93-4	C14H15Cl2N3O2	327.0541	328.0614	326.0468	ESI+	12
284	Ethametsulfuron-methyl	97780-06-8	C15H18N6O6S	410.1009	411.1082	409.0936	ESI+	6
285	Ethephon	16672-87-0	C2H6ClO3P	143.9743	144.9816	142.9670	ESI-	6
286	Ethidimuron	30043-49-3	C7H12N4O3S2	264.0351	265.0424	263.0278	ESI+	11
287	Ethiofencarb	29973-13-5	C11H15NO2S	225.0823	226.0896	224.0750	ESI+	10
288	Ethiofencarb-sulfone	53380-23-7	C11H15NO4S	257.0722	258.0795	256.0649	ESI+	8
289	Ethiofencarb-sulfoxide	53380-22-6	C11H15NO3S	241.0773	242.0846	240.0700	ESI+	4
290	Ethion	563-12-2	C9H22O4P2S4	383.9876	384.9949	382.9803	ESI+	6
291	Ethiprole	181587-01-9	C13H9Cl2F3N4OS	395.9826	396.9899	394.9753	ESI+	30
292	Ethirimol	23947-60-6	C11H19N3O	209.1528	210.1601	208.1455	ESI+	6
293	Ethofumesate	26225-79-6	C13H18O5S	286.0875	287.0948	285.0802	ESI+	11
294	Ethoprophos	13194-48-4	C8H19O2PS2	242.0564	243.0637	241.0491	ESI+	6
295	Ethoxyquin	91-53-2	C14H19NO	217.1467	218.1540	216.1394	ESI+	4
296	Ethoxysulfuron	126801-58-9	C15H18N4O7S	398.0896	399.0969	397.0823	ESI+	6
297	Ethylenthionurea	96-45-7	C3H6N2S	102.0252	103.0325	101.0179	ESI+	6
298	Etofenprox	80844-07-1	C25H28O3	376.2038	377.2111	375.1965	ESI+	6
299	Etoazole	153233-91-1	C21H23F2NO2	359.1697	360.1770	358.1624	ESI+	6
300	Etrinfos	38260-54-7	C10H17N2O4PS	292.0647	293.0720	291.0574	ESI+	6

Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
301 Famoxadone	131807-57-3	C22H18N2O4	374.1267	375.1340	373.1194	ESI+	6
302 Famphur	52-85-7	C10H16NO5PS2	325.0208	326.0281	324.0135	ESI+	12
303 Fenamidone	161326-34-7	C17H17N3O5	311.1092	312.1165	310.1019	ESI+	6
304 Fenaminosulf	140-56-7	C8H10N3NaO3S	251.0341	252.0414	250.0268	ESI+	2
305 Fenamiphos	22224-92-6	C13H22NO3PS	303.1058	304.1131	302.0985	ESI+	6
306 Fenamiphos-sulfone	31972-44-8	C13H22NO5PS	335.0956	336.1029	334.0883	ESI+	6
307 Fenamiphos-sulfoxide	31972-43-7	C13H22NO4PS	319.1007	320.1080	318.0934	ESI+	6
308 Fenarimol	60168-88-9	C17H12Cl2N2O	330.0327	331.0400	329.0254	ESI+	12
309 Fenazaquin	120928-09-8	C20H22N2O	306.1732	307.1805	305.1659	ESI+	6
310 Fenazox	495-48-7	C12H10N2O	198.0793	199.0866	197.0720	ESI+	6
311 Fenbuconazole	114369-43-6	C19H17ClN4	336.1142	337.1215	335.1069	ESI+	8
312 Fenbutatin-oxide	13356-08-6	C60H78OSn2	1054.4121	1055.4194	1053.4048	ESI+	11
313 Fenchlorazol-ethyl	103112-35-2	C12H8Cl5N3O2	400.9059	401.9132	399.8986	ESI+	118
314 Fenfuram	24691-80-3	C12H11NO2	201.0790	202.0863	200.0717	ESI+	6
315 Fenhexamid	126833-17-8	C14H17Cl2NO2	301.0636	302.0709	300.0563	ESI+	23
316 Fenitrothion	122-14-5	C9H12NO5PS	277.0174	278.0247	276.0101	ESI+	2
317 Fenobucarb	3766-81-2	C12H17NO2	207.1259	208.1332	206.1186	ESI+	6
318 Fenoprop	93-72-1	C9H7Cl3O3	267.9461	268.9534	266.9388	ESI-	8
319 Fenothiocarb	62850-32-2	C13H19NO2S	253.1136	254.1209	252.1063	ESI+	4
320 Fenoxanil	115852-48-7	C15H18Cl2N2O2	328.0745	329.0818	327.0672	ESI+	29
321 Fenoxaprop	95617-09-7	C16H12ClNO5	333.0404	334.0477	332.0331	ESI+	23
322 Fenoxaprop-ethyl	66441-23-4	C18H16ClNO5	361.0717	362.0790	360.0644	ESI+	12
323 Fenoxaprop-P-ethyl	71283-80-2	C18H16ClNO5	361.0717	362.0790	360.0644	ESI+	12
324 Fenoxycarb	79127-80-3	C17H19NO4	301.1314	302.1387	300.1241	ESI+	6
325 Fenpropathrin	64257-84-7	C22H23NO3	349.1678	350.1751	348.1605	ESI+	11
326 Fenpropidin	67306-00-7	C19H31N	273.2457	274.2530	272.2384	ESI+	6
327 Fenpropimorph	67564-91-4	C20H33NO	303.2562	304.2635	302.2489	ESI+	6
328 Fensulfothion	115-90-2	C11H17O4PS2	308.0306	309.0379	307.0233	ESI+	6
329 Fensulfothion-oxon	6552-21-2	C11H17O5PS	292.0534	293.0607	291.0461	ESI+	6
330 Fensulfothion-oxon-sulfone	6132-17-8	C11H17O6PS	308.0483	309.0556	307.0410	ESI+	4
331 Fensulfothion-sulfone	14255-72-2	C11H17O5PS2	324.0255	325.0328	323.0182	ESI+	6
332 Fenthion	55-38-9	C10H15O3PS2	278.0200	279.0273	277.0127	ESI+	6
333 Fenthion-oxon	6552-12-1	C10H15O4PS	262.0429	263.0502	261.0356	ESI+	6
334 Fenthion-oxon-sulfone	14086-35-2	C10H15O6PS	294.0327	295.0400	293.0254	ESI+	12
335 Fenthion-oxon-sulfoxide	6552-13-2	C10H15O5PS	278.0378	279.0451	277.0305	ESI+	3
336 Fenthion-sulfone	3761-42-0	C10H15O5PS2	310.0099	311.0172	309.0026	ESI+	4
337 Fenthion-sulfoxide	3761-41-9	C10H15O4PS2	294.0149	295.0222	293.0076	ESI+	6
338 Fenuron	101-42-8	C9H12N2O	164.0950	165.1023	163.0877	ESI+	6
339 Fenvalerate	51630-58-1	C25H22ClNO3	419.1288	420.1361	418.1215	ESI+	6
340 Fipronil	120068-37-3	C12H4Cl2F6N4O5	435.9387	436.9460	434.9314	ESI-	12
341 Fipronil-desulfinyl	205650-65-3	C12H4Cl2F6N4	387.9717	388.9790	386.9644	ESI-	12
342 Fipronil-sulfide	120067-83-6	C12H4Cl2F6N4S	419.9438	420.9511	418.9365	ESI-	12
343 Fipronil-sulfone	120068-36-2	C12H4Cl2F6N4O2S	451.9336	452.9409	450.9263	ESI-	12
344 Flamprop-isopropyl	52756-22-6	C19H19ClFNO3	363.1037	364.1110	362.0964	ESI+	10
345 Flamprop-methyl	52756-25-9	C17H15ClFNO3	335.0724	336.0797	334.0651	ESI+	4
346 Flamprop-M-isopropyl	63782-90-1	C19H19ClFNO3	363.1037	364.1110	362.0964	ESI+	10
347 Flazasulfuron	104040-78-0	C13H12F3N5O5S	407.0511	408.0584	406.0438	ESI+	6
348 Flooumafen	90035-08-8	C33H25F3O4	542.1705	543.1778	541.1632	ESI+	12
349 Flonicamid	158062-67-0	C9H6F3N3O	229.0463	230.0536	228.0390	ESI+	8
350 Florasulam	145701-23-1	C12H8F3N5O3S	359.0300	360.0373	358.0227	ESI+	2
351 Fluacrypyrim	229977-93-9	C20H21F3N2O5	426.1403	427.1476	425.1330	ESI+	6
352 Fluazifop	69335-91-7	C15H12F3NO4	327.0718	328.0791	326.0645	ESI+	12
353 Fluazifop-butyl	69806-50-4	C19H20F3NO4	383.1344	384.1417	382.1271	ESI+	6
354 Fluazifop-P (free acid)	83066-88-0	C15H12F3NO4	327.0718	328.0791	326.0645	ESI+	12
355 Fluazifop-P-butyl	79241-46-6	C19H20F3NO4	383.1344	384.1417	382.1271	ESI+	6
356 Fluazinam	79622-59-6	C13H4Cl2F6N4O4	463.9514	464.9587	462.9441	ESI-	12
357 Fluazuron	86811-58-7	C20H10Cl2F5N3O3	505.0019	506.0092	503.9946	ESI+	17
358 Flubendiamide	272451-65-7	C23H22F7IN2O4S	682.0233	683.0306	681.0160	ESI+	5
359 Flucycloxuron	94050-52-9	C25H20ClF2N3O3	483.1161	484.1234	482.1088	ESI+	10
360 Flucythrinate	70124-77-5	C26H23F2NO4	451.1595	452.1668	450.1522	ESI+	4
361 Fludioxonil	131341-86-1	C12H6F2N2O2	248.0397	249.0470	247.0324	ESI-	6
362 Flufenacet	142459-58-3	C14H13F4N3O2S	363.0665	364.0738	362.0592	ESI+	6
363 Flufenoxuron	101463-69-8	C21H11ClF6N2O3	488.0362	489.0435	487.0289	ESI+	8
364 Flumetralin	62924-70-3	C16H12ClF4N3O4	421.0452	422.0525	420.0379	ESI+	3
365 Flumetsulam	98967-40-9	C12H9F2N5O2S	325.0445	326.0518	324.0372	ESI+	2
366 Flumioxazin	103361-09-7	C19H15FN2O4	354.1016	355.1089	353.0943	ESI+	2
367 Fluometuron	2164-17-2	C10H11F3N2O	232.0823	233.0896	231.0750	ESI+	4
368 Fluopicolide	239110-15-7	C14H8Cl3F3N2O	381.9654	382.9727	380.9581	ESI+	11
369 Fluopyram	658066-35-4	C16H11ClF6N2O	396.0464	397.0537	395.0391	ESI+	12
370 Fluoroglycofen-ethyl	77501-90-7	C18H13ClF3NO7	447.0333	448.0406	446.0260	ESI+	12
371 Fluoxastrobin	361377-29-9	C21H16ClFN4O5	458.0793	459.0866	457.0720	ESI+	12
372 Flupyrsulfuron-methyl	144740-54-5	C15H14F3N5O7S	465.0566	466.0639	464.0493	ESI+	12
373 Fluquinconazole	136426-54-5	C16H8Cl2FN5O	375.0090	376.0163	374.0017	ESI+	10
374 Fluridone	59756-60-4	C19H14F3NO	329.1027	330.1100	328.0954	ESI+	4
375 Flurochloridone	61213-25-0	C12H10Cl2F3NO	311.0092	312.0165	310.0019	ESI+	18

Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
376 Fluroxypyr	69377-81-7	C7H5Cl2FN2O3	253.9661	254.9734	252.9588	ESI+	15
377 Fluroxypyr-1-methylheptylester	81406-37-3	C15H21Cl2FN2O3	366.0913	367.0986	365.0840	ESI+	12
378 Flurprimidol	56425-91-3	C15H15F3N2O2	312.1086	313.1159	311.1013	ESI+	6
379 Flurtamone	96525-23-4	C18H14F3NO2	333.0977	334.1050	332.0904	ESI+	6
380 Flusilazole	85509-19-9	C16H15F2N3Si	315.1003	316.1076	314.0930	ESI+	6
381 Fluthiacet-methyl	117337-19-6	C15H15ClFN3O3S2	403.0227	404.0300	402.0154	ESI+	9
382 Flutolanil	66332-96-5	C17H16F3NO2	323.1133	324.1206	322.1060	ESI+	12
383 Flutriafol	76674-21-0	C16H13F2N3O	301.1027	302.1100	300.0954	ESI+	5
384 Fluxapyroxad	907204-31-3	C18H12F5N3O	381.0901	382.0974	380.0828	ESI+	11
385 Fomesafen	72178-02-0	C15H10ClF3N2O6S	437.9900	438.9973	436.9827	ESI+	21
386 Fonofos	944-22-9	C10H15OPS2	246.0302	247.0375	245.0229	ESI+	5
387 Foramsulfuron	173159-57-4	C17H20N6O7S	452.1114	453.1187	451.1041	ESI+	6
388 Forchlorfenuron	68157-60-8	C12H10ClN3O	247.0512	248.0585	246.0439	ESI+	12
389 Fosetyl-aluminium	39148-24-8	C2H7O3P	110.0133	111.0206	109.0060	ESI-	3
390 Fosthiazate	98886-44-3	C9H18NO3PS2	283.0466	284.0539	282.0393	ESI+	6
391 Fuberidazole	3878-19-1	C11H8N2O	184.0637	185.0710	183.0564	ESI+	6
392 Furalaxyl	57646-30-7	C17H19NO4	301.1314	302.1387	300.1241	ESI+	3
393 Furametpyr	123572-88-3	C17H20ClN3O2	333.1244	334.1317	332.1171	ESI+	12
394 Furathiocarb	65907-30-4	C18H26N2O5S	382.1562	383.1635	381.1489	ESI+	6
395 Furmecyclox	60568-05-0	C14H21NO3	251.1521	252.1594	250.1448	ESI+	6
396 Gibberellic acid (Gibberellin)	77-06-5	C19H22O6	346.1416	347.1489	345.1343	ESI-	11
397 Glufosinate	77182-82-2	C5H12NO4P	181.0504	182.0577	180.0431	ESI+	10
398 Glyphosate	1071-83-6	C3H8NO5P	169.0140	170.0213	168.0067	ESI+	8
399 Halofenozide	112226-61-6	C18H19ClN2O2	330.1135	331.1208	329.1062	ESI+	12
400 Halosulfuron-methyl	100784-20-1	C13H15ClN6O7S	434.0411	435.0484	433.0338	ESI+	11
401 Haloxyfop	69806-34-4	C15H11ClF3NO4	361.0329	362.0402	360.0256	ESI+	9
402 Haloxyfop-2-ethoxyethyl	87237-48-7	C19H19ClF3NO5	433.0904	434.0977	432.0831	ESI+	12
403 Haloxyfop-methyl	69806-40-2	C16H13ClF3NO4	375.0485	376.0558	374.0412	ESI+	12
404 Haloxyfop-R-methyl	72619-32-0	C16H13ClF3NO4	375.0485	376.0558	374.0412	ESI+	12
405 Heptenophos	23560-59-0	C9H12ClO4P	250.0162	251.0235	249.0089	ESI+	9
406 Hexaconazole	79983-71-4	C14H17Cl2N3O	313.0749	314.0822	312.0676	ESI+	10
407 Hexaflumuron	86479-06-3	C16H8Cl2F6N2O3	459.9816	460.9889	458.9743	ESI-	12
408 Hexazinone	51235-04-2	C12H20N4O2	252.1586	253.1659	251.1513	ESI+	3
409 Hexythiazox	78587-05-0	C17H21ClN2O2S	352.1012	353.1085	351.0939	ESI+	11
410 Hydramethylnon	67485-29-4	C25H24F6N4	494.1905	495.1978	493.1832	ESI+	12
411 Hymexazol	10004-44-1	C4H5NO2	99.0320	100.0393	98.0247	ESI+	3
412 Imazalil	35554-44-0	C14H14Cl2N2O	296.0483	297.0556	295.0410	ESI+	12
413 Imazamethabenz-methyl	81405-85-8	C16H20N2O3	288.1474	289.1547	287.1401	ESI+	12
414 Imazamox	114311-32-9	C15H19N3O4	305.1376	306.1449	304.1303	ESI+	10
415 Imazapic	104098-48-8	C14H17N3O3	275.1270	276.1343	274.1197	ESI+	11
416 Imazapyr	81334-34-1	C13H15N3O3	261.1113	262.1186	260.1040	ESI+	11
417 Imazaquin	81335-37-7	C13H17N3O3	311.1270	312.1343	310.1197	ESI+	6
418 Imazethapyr	81335-77-5	C15H19N3O3	289.1426	290.1499	288.1353	ESI+	12
419 Imazosulfuron	122548-33-8	C14H13ClN6O5S	412.0357	413.0430	411.0284	ESI+	13
420 Imibenconazole	86598-92-7	C17H13Cl3N4S	409.9927	411.0000	408.9854	ESI+	22
421 Imidacloprid	138261-41-3	C9H10ClN5O2	255.0523	256.0596	254.0450	ESI+	8
422 Indanofan	133220-30-1	C20H17ClO3	340.0866	341.0939	339.0793	ESI+	6
423 Indoxacarb	173584-44-6	C22H17ClF3N3O7	527.0707	528.0780	526.0634	ESI+	12
424 Iodosulfuron-methyl	144550-36-7	C14H14IN5O6S	506.9710	507.9783	505.9637	ESI+	8
425 Ioxynil	1689-83-4	C7H3I2NO	370.8304	371.8377	369.8231	ESI-	4
426 Ipconazole	125225-28-7	C18H24ClN3O	333.1608	334.1681	332.1535	ESI+	5
427 Iprobenfos	26087-47-8	C13H21O3PS	288.0949	289.1022	287.0876	ESI+	3
428 Iprodione	36734-19-7	C13H13Cl2N3O3	329.0334	330.0407	328.0261	ESI+	4
429 Iprovalicarb	140923-17-7	C18H28N2O3	320.2100	321.2173	319.2027	ESI+	6
430 Irgarol 1051	28159-98-0	C11H19N5S	253.1361	254.1434	252.1288	ESI+	6
431 Isazofos	42509-80-8	C9H17ClN3O3PS	313.0417	314.0490	312.0344	ESI+	12
432 Isocarbamid	30979-48-7	C8H15N3O2	185.1164	186.1237	184.1091	ESI+	6
433 Isocarbofos	24353-61-5	C11H16N4O4PS	289.0538	290.0611	288.0465	ESI+	6
434 Isofenphos	25311-71-1	C15H24NO4PS	345.1164	346.1237	344.1091	ESI+	6
435 Isofenphos-methyl	99675-03-3	C14H22NO4PS	331.1007	332.1080	330.0934	ESI+	6
436 Isofenphos-oxon	31120-85-1	C15H24NO5P	329.1392	330.1465	328.1319	ESI+	3
437 Isomethiozin	57052-04-7	C12H20N4O5S	268.1358	269.1431	267.1285	ESI+	6
438 Isonoruron	28805-78-9	C13H22N2O	222.1732	223.1805	221.1659	ESI+	6
439 Isoprocab	2631-40-5	C11H15NO2	193.1103	194.1176	192.1030	ESI+	3
440 Isopropalin	33820-53-0	C15H23N3O4	309.1689	310.1762	308.1616	ESI+	6
441 Isoprothiolane	50512-35-1	C12H18O4S2	290.0647	291.0720	289.0574	ESI+	6
442 Isoproturon	34123-59-6	C12H18N2O	206.1419	207.1492	205.1346	ESI+	6
443 Isopyrazam	881685-58-1	C20H23F2N3O	359.1809	360.1882	358.1736	ESI+	9
444 Isoxaben	82558-50-7	C18H24N2O4	332.1736	333.1809	331.1663	ESI+	6
445 Isoxadifen-ethyl	163520-33-0	C18H17NO3	295.1208	296.1281	294.1135	ESI+	12
446 Isoxaflutole	141112-29-0	C15H12F3NO4S	359.0439	360.0512	358.0366	ESI+	5
447 Isoxathion	18854-01-8	C13H16NO4PS	313.0538	314.0611	312.0465	ESI+	6
448 Ivermectine	70288-86-7	C48H74O14	874.5079	875.5152	873.5006	ESI+	6
449 Karbutilate	4849-32-5	C14H21N3O3	279.1583	280.1656	278.1510	ESI+	16
450 Kasugamycin	6980-18-3	C14H25N3O9	379.1591	380.1664	378.1518	ESI+	3

Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
451 Kresoxim-methyl	143390-89-0	C18H19NO4	313.1314	314.1387	312.1241	ESI+	6
452 Lactofen	77501-63-4	C19H15ClF3NO7	461.0489	462.0562	460.0416	ESI+	12
453 Lambda-Cyhalothrin	91465-08-6	C23H19ClF3NO3	449.1006	450.1079	448.0933	ESI+	4
454 Lenacil	2164-08-1	C13H18N2O2	234.1368	235.1441	233.1295	ESI+	6
455 Linuron	330-55-2	C9H10Cl2N2O2	248.0119	249.0192	247.0046	ESI+	12
456 Lufenuron	103055-07-8	C17H8Cl2F8N2O3	509.9784	510.9857	508.9711	ESI-	10
457 Malaoxon	1634-78-2	C10H19O7PS	314.0589	315.0662	313.0516	ESI+	6
458 Malathion	121-75-5	C10H19O6PS2	330.0361	331.0434	329.0288	ESI+	12
459 Maleic-hydrazide	123-33-1	C4H4N2O2	112.0273	113.0346	111.0200	ESI+	3
460 Mandipropamid	374726-62-2	C23H22ClNO4	411.1237	412.1310	410.1164	ESI+	12
461 MCPA (MCP)	94-74-6	C9H9ClO3	200.0240	201.0313	199.0167	ESI-	3
462 MCPA-butoxyethyl ester	19480-43-4	C15H21ClO4	300.1128	301.1201	299.1055	ESI+	12
463 MCPB	94-81-5	C11H13ClO3	228.0553	229.0626	227.0480	ESI-	3
464 Mecarbam	2595-54-2	C10H20NO5PS2	329.0521	330.0594	328.0448	ESI+	6
465 Mecoprop (MCP)	93-65-2	C10H11ClO3	214.0397	215.0470	213.0324	ESI-	2
466 Mecoprop-P	16484-77-8	C10H11ClO3	214.0397	215.0470	213.0324	ESI-	4
467 Mefenacet	73250-68-7	C16H14N2O2S	298.0776	299.0849	297.0703	ESI+	6
468 Mefenpyr-diethyl	135590-91-9	C16H18Cl2N2O4	372.0644	373.0717	371.0571	ESI+	24
469 Mefluidide	53780-34-0	C11H13F3N2O3S	310.0599	311.0672	309.0526	ESI+	10
470 Mepanipyrim	110235-47-7	C14H13N3	223.1109	224.1182	222.1036	ESI+	6
471 Mephosfolan	950-10-7	C8H16NO3PS2	269.0309	270.0382	268.0236	ESI+	6
472 Mepiquat	24307-26-4	C7H16N	114.1283	115.1356	113.1210	ESI+	6
473 Mepronil	55814-41-0	C17H19NO2	269.1416	270.1489	268.1343	ESI+	5
474 Meptyldinocap	6119-92-2	C18H24N2O6	364.1634	365.1707	363.1561	ESI-	6
475 Mesosulfuron-methyl	208465-21-8	C17H21N5O9S2	503.0781	504.0854	502.0708	ESI+	6
476 Mesotrione	104206-82-8	C14H13NO7S	339.0413	340.0486	338.0340	ESI+	6
477 Metaflumizone	139968-49-3	C24H16F6N4O2	506.1177	507.1250	505.1104	ESI+	6
478 Metalaxyl	57837-19-1	C15H21NO4	279.1471	280.1544	278.1398	ESI+	6
479 Metalaxyl-M	70630-17-0	C15H21NO4	279.1471	280.1544	278.1398	ESI+	6
480 Metamitron	41394-05-2	C10H10N4O	202.0855	203.0928	201.0782	ESI+	3
481 Metazachlor	67129-08-2	C14H16ClN3O	277.0982	278.1055	276.0909	ESI+	6
482 Metconazole	125116-23-6	C17H22ClN3O	319.1451	320.1524	318.1378	ESI+	4
483 Methabenzthiazuron	18691-97-9	C10H11N3OS	221.0623	222.0696	220.0550	ESI+	6
484 Methacrifos	62610-77-9	C7H13O5PS	240.0221	241.0294	239.0148	ESI+	12
485 Methamidophos	10265-92-6	C2H8NO2PS	141.0013	142.0086	139.9940	ESI+	6
486 Methfuroxam	28730-17-8	C14H15NO2	229.1103	230.1176	228.1030	ESI+	4
487 Methidathion	950-37-8	C6H11N2O4PS3	301.9619	302.9692	300.9546	ESI+	7
488 Methiocarb	2032-65-7	C11H15NO2S	225.0823	226.0896	224.0750	ESI+	6
489 Methiocarb-sulfone	2179-25-1	C11H15NO4S	257.0722	258.0795	256.0649	ESI+	9
490 Methiocarb-sulfoxide	2635-10-1	C11H15NO3S	241.0773	242.0846	240.0700	ESI+	6
491 Methomyl	16752-77-5	C5H10N2O2S	162.0463	163.0536	161.0390	ESI+	6
492 Methoprene	40596-69-8	C19H34O3	310.2508	311.2581	309.2435	ESI+	12
493 Methoprotryne	841-06-5	C11H21N5OS	271.1467	272.1540	270.1394	ESI+	6
494 Methoxyfenozide	161050-58-4	C22H28N2O3	368.2100	369.2173	367.2027	ESI+	6
495 Metobromuron	3060-89-7	C9H11BrN2O2	258.0004	259.0077	256.9931	ESI+	12
496 Metolachlor	51218-45-2	C15H22ClNO2	283.1339	284.1412	282.1266	ESI+	12
497 Metolcarb	1129-41-5	C9H11NO2	165.0790	166.0863	164.0717	ESI+	6
498 Metominostrobin	133408-50-1	C16H16N2O3	284.1161	285.1234	283.1088	ESI+	6
499 Metosulam	139528-85-1	C14H13Cl2N5O4S	417.0065	418.0138	415.9992	ESI+	24
500 Metoxuron	19937-59-8	C10H13ClN2O2	228.0666	229.0739	227.0593	ESI+	5
501 Metrafenone	220899-03-6	C19H21BrO5	408.0572	409.0645	407.0499	ESI+	12
502 Metribuzin	21087-64-9	C8H14N4OS	214.0888	215.0961	213.0815	ESI+	5
503 Metsulfuron-methyl	74223-64-6	C14H15N5O6S	381.0743	382.0816	380.0670	ESI+	6
504 Mevinphos	7786-34-7	C7H13O6P	224.0450	225.0523	223.0377	ESI+	5
505 Mexacarbate	315-18-4	C12H18N2O2	222.1368	223.1441	221.1295	ESI+	6
506 Molinate	2212-67-1	C9H17NOS	187.1031	188.1104	186.0958	ESI+	6
507 Monalide	7287-36-7	C13H18ClNO	239.1077	240.1150	238.1004	ESI+	20
508 Monocrotophos	6923-22-4	C7H14NO5P	223.0610	224.0683	222.0537	ESI+	12
509 Monolinuron	1746-81-2	C9H11ClN2O2	214.0509	215.0582	213.0436	ESI+	10
510 Monuron	150-68-5	C9H11ClN2O	198.0560	199.0633	197.0487	ESI+	10
511 Morpholine	110-91-8	C4H9NO	87.0684	88.0757	86.0611	ESI+	6
512 Moxidectin	113507-06-5	C37H53NO8	639.3771	640.3844	638.3698	ESI+	12
513 Myclobutanil	88671-89-0	C15H17ClN4	288.1142	289.1215	287.1069	ESI+	8
514 <i>N</i> -(2, 4-Dimethylphenyl) formamide	60397-77-5	C9H11NO	149.0841	150.0914	148.0768	ESI+	6
515 <i>N</i> -(2, 4-Dimethylphenyl) - <i>N'</i> -methylformamidine	33089-74-6	C10H14N2	162.1157	163.1230	161.1084	ESI+	6
516 <i>N</i> , <i>N'</i> -Diphenylurea	102-07-8	C13H12N2O	212.0950	213.1023	211.0877	ESI+	4
517 Naled	300-76-5	C4H7Br2Cl2O4P	377.7826	378.7899	376.7753	ESI+	6
518 Naproanilide	52570-16-8	C19H17NO2	291.1259	292.1332	290.1186	ESI+	2
519 Napropamide	15299-99-7	C17H21NO2	271.1572	272.1645	270.1499	ESI+	6
520 Naptalam	132-66-1	C18H13NO3	291.0895	292.0968	290.0822	ESI+	6
521 Neburon	555-37-3	C12H16Cl2N2O	274.0640	275.0713	273.0567	ESI+	9
522 Nicarbazin	330-95-0	C19H18N6O6	426.1288	427.1361	425.1215	ESI-	3
523 Nicosulfuron	111991-09-4	C15H18N6O6S	410.1009	411.1082	409.0936	ESI+	6
524 Nicotine	54-11-5	C10H14N2	162.1157	163.1230	161.1084	ESI+	6
525 Nitenpyram	150824-47-8	C11H15ClN4O2	270.0884	271.0957	269.0811	ESI+	6

Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
526 Nitralin	4726-14-1	C13H19N3O6S	345.0995	346.1068	344.0922	ESI+	12
527 Nitrothal-isopropyl	10552-74-6	C14H17NO6	295.1056	296.1129	294.0983	ESI+	6
528 Norflurazon	27314-13-2	C12H9ClF3N3O	303.0386	304.0459	302.0313	ESI+	16
529 Norflurazon-desmethyl	23576-24-1	C11H7ClF3N3O	289.0230	290.0303	288.0157	ESI+	12
530 Novaluron	116714-46-6	C17H9ClF8N2O4	492.0123	493.0196	491.0050	ESI+	18
531 Noviflumuron	121451-02-3	C17H7Cl2F9N2O3	527.9690	528.9763	526.9617	ESI-	6
532 Nuarimol	63284-71-9	C17H12ClF8N2O	314.0622	315.0695	313.0549	ESI+	11
533 Ofurace	58810-48-3	C14H16ClNO3	281.0819	282.0892	280.0746	ESI+	17
534 Omethoate	1113-02-6	C5H12NO4PS	213.0225	214.0298	212.0152	ESI+	4
535 Orbencarb	34622-58-7	C12H16ClNOS	257.0641	258.0714	256.0568	ESI+	12
536 Orthosulfamuron	213464-77-8	C16H20N6O6S	424.1165	425.1238	423.1092	ESI+	6
537 Oryzalin	19044-88-3	C12H18N4O6S	346.0947	347.1020	345.0874	ESI+	5
538 Oxabetrinil	94593-79-0	C12H12N2O3	232.0848	233.0921	231.0775	ESI+	2
539 Oxadiargyl	39807-15-3	C15H14Cl2N2O3	340.0381	341.0454	339.0308	ESI+	14
540 Oxadiazon	19666-30-9	C15H18Cl2N2O3	344.0694	345.0767	343.0621	ESI+	6
541 Oxadixyl	77732-09-3	C14H18N2O4	278.1267	279.1340	277.1194	ESI+	12
542 Oxamyl	23135-22-0	C7H13N3O3S	219.0678	220.0751	218.0605	ESI+	3
543 Oxasulfuron	144651-06-9	C17H18N4O6S	406.0947	407.1020	405.0874	ESI+	6
544 Oxaziclomefone	153197-14-9	C20H19Cl2NO2	375.0793	376.0866	374.0720	ESI+	2
545 Oxyacboxin	5259-88-1	C12H13NO4S	267.0565	268.0638	266.0492	ESI+	3
546 Oxydemeton-methyl	301-12-2	C6H15O4PS2	246.0149	247.0222	245.0076	ESI+	6
547 Paclobutrazol	76738-62-0	C15H20ClN3O	293.1295	294.1368	292.1222	ESI+	8
548 Paraoxon-ethyl	311-45-5	C10H14NO6P	275.0559	276.0632	274.0486	ESI+	6
549 Paraoxon-methyl	950-35-6	C8H10NO6P	247.0246	248.0319	246.0173	ESI+	3
550 Paraquat	1910-42-5	C12H14Cl2N2	256.0534	257.0607	255.0461	ESI+	5
551 Parathion	56-38-2	C10H14NO5PS	291.0330	292.0403	290.0257	ESI+	3
552 Pebulate	1114-71-2	C10H21NOS	203.1344	204.1417	202.1271	ESI+	6
553 Penconazole	66246-88-6	C13H15Cl2N3	283.0643	284.0716	282.0570	ESI+	12
554 Pencycuron	66063-05-6	C19H21ClN2O	328.1342	329.1415	327.1269	ESI+	10
555 Pendimethalin	40487-42-1	C13H19N3O4	281.1376	282.1449	280.1303	ESI+	6
556 Penoxsulam	219714-96-2	C16H14F5N5O5S	483.0636	484.0709	482.0563	ESI+	6
557 Pentachlorophenol	87-86-5	C6HCl5O	263.8470	264.8543	262.8397	ESI-	3
558 Pentoxazone	110956-75-7	C17H17ClFNO4	353.0830	354.0903	352.0757	ESI+	2
559 Permethrin	52645-53-1	C21H20Cl2O3	390.0790	391.0863	389.0717	ESI+	12
560 Pethoxamid	106700-29-2	C16H22ClNO2	295.1339	296.1412	294.1266	ESI+	7
561 Phenmedipham	13684-63-4	C16H16N2O4	300.1110	301.1183	299.1037	ESI+	6
562 Phenothrin	26002-80-2	C23H26O3	350.1882	351.1955	349.1809	ESI+	9
563 Phenthoate	2597-03-7	C12H17O4PS2	320.0306	321.0379	319.0233	ESI+	12
564 Phorate	298-02-2	C7H17O2PS3	260.0128	261.0201	259.0055	ESI+	6
565 Phorate-oxon	2600-69-3	C7H17O3PS2	244.0357	245.0430	243.0284	ESI+	6
566 Phorate-sulfone	2588-04-7	C7H17O4PS3	292.0027	293.0100	290.9954	ESI+	6
567 Phorate-sulfoxide	2588-03-6	C7H17O3PS3	276.0077	277.0150	275.0004	ESI+	6
568 Phosalone	2310-17-0	C12H15ClNO4PS2	366.9869	367.9942	365.9796	ESI+	12
569 Phosfolan	947-02-4	C7H14NO3PS2	255.0153	256.0226	254.0080	ESI+	6
570 Phosmet	732-11-6	C11H12NO4PS2	316.9945	318.0018	315.9872	ESI+	12
571 Phosphamidon	13171-21-6	C10H19ClNO5P	299.0689	300.0762	298.0616	ESI+	12
572 Phoxim	14816-18-3	C12H15N2O3PS	298.0541	299.0614	297.0468	ESI+	6
573 Picloram	1918-02-1	C6H3Cl3N2O2	239.9260	240.9333	238.9187	ESI+	9
574 Picolinafen	137641-05-5	C19H12F4N2O2	376.0835	377.0908	375.0762	ESI+	6
575 Picoxystrobin	117428-22-5	C18H16F3NO4	367.1031	368.1104	366.0958	ESI+	6
576 Pinoxaden	243973-20-8	C23H32N2O4	400.2362	401.2435	399.2289	ESI+	6
577 Piperonyl-butoxide	51-03-6	C19H30O5	338.2093	339.2166	337.2020	ESI+	12
578 Piperophos	24151-93-7	C14H28NO3PS2	353.1248	354.1321	352.1175	ESI+	6
579 Pirimicarb	23103-98-2	C11H18N4O2	238.1430	239.1503	237.1357	ESI+	3
580 Pirimicarb-desmethyl	30614-22-3	C10H16N4O2	224.1273	225.1346	223.1200	ESI+	6
581 Pirimicarb-desmethyl-formamido	27218-04-8	C11H16N4O3	252.1222	253.1295	251.1149	ESI+	2
582 Pirimiphos-ethyl	23505-41-1	C13H24N3O3PS	333.1276	334.1349	332.1203	ESI+	6
583 Pirimiphos-methyl	29232-93-7	C11H20N3O3PS	305.0963	306.1036	304.0890	ESI+	6
584 Prallethrin	23031-36-9	C19H24O3	300.1725	301.1798	299.1652	ESI+	6
585 Pretilachlor	51218-49-6	C17H26ClNO2	311.1652	312.1725	310.1579	ESI+	6
586 Primisulfuron-methyl	86209-51-0	C15H12F4N4O7S	468.0363	469.0436	467.0290	ESI+	9
587 Probenazole	27605-76-1	C10H9NO3S	223.0303	224.0376	222.0230	ESI+	4
588 Prochloraz	67747-09-5	C15H16Cl3N3O2	375.0308	376.0381	374.0235	ESI+	15
589 Profenofos	41198-08-7	C11H15BrClO3PS	371.9351	372.9424	370.9278	ESI+	12
590 Profoxydim	139001-49-3	C24H32ClNO4S	465.1741	466.1814	464.1668	ESI+	24
591 Promecarb	2631-37-0	C12H17NO2	207.1259	208.1332	206.1186	ESI+	6
592 Prometon	1610-18-0	C10H19NO5O	225.1590	226.1663	224.1517	ESI+	6
593 Prometryn	7287-19-6	C10H19N5S	241.1361	242.1434	240.1288	ESI+	6
594 Propachlor	1918-16-7	C11H14ClNO	211.0764	212.0837	210.0691	ESI+	6
595 Propamocarb	24579-73-5	C9H20N2O2	188.1525	189.1598	187.1452	ESI+	6
596 Propanil	709-98-8	C9H9Cl2NO	217.0061	218.0134	215.9988	ESI+	9
597 Propaphos	7292-16-2	C13H21O4PS	304.0898	305.0971	303.0825	ESI+	10
598 Propaquizafop	111479-05-1	C22H22ClN3O5	443.1248	444.1321	442.1175	ESI+	12
599 Propargite	2312-35-8	C19H26O4S	350.1552	351.1625	349.1479	ESI+	6
600 Propazine	139-40-2	C9H16ClN5	229.1094	230.1167	228.1021	ESI+	6

Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
601 Propetamphos	31218-83-4	C10H20NO4PS	281.0851	282.0924	280.0778	ESI+	18
602 Propham	122-42-9	C10H13NO2	179.0946	180.1019	178.0873	ESI+	6
603 Propiconazole	60207-90-1	C15H17Cl2N3O2	341.0698	342.0771	340.0625	ESI+	12
604 Propisochlor	86763-47-5	C15H22ClNO2	283.1339	284.1412	282.1266	ESI+	12
605 Propoxur	114-26-1	C11H15NO3	209.1052	210.1125	208.0979	ESI+	6
606 Propoxycarbazone	181274-15-7	C15H18N4O7S	398.0896	399.0969	397.0823	ESI+	20
607 Propylene-thiourea	2122-19-2	C4H8N2S	116.0408	117.0481	115.0335	ESI+	6
608 Propyzamide	23950-58-5	C12H11Cl2NO	255.0218	256.0291	254.0145	ESI+	10
609 Proquinazid	189278-12-4	C14H17IN2O2	372.0335	373.0408	371.0262	ESI+	6
610 Prosulfocarb	52888-80-9	C14H21NOS	251.1344	252.1417	250.1271	ESI+	4
611 Prosulfuron	94125-34-5	C15H16F3N5O4S	419.0875	420.0948	418.0802	ESI+	7
612 Prothioconazole	178928-70-6	C14H15Cl2N3OS	343.0313	344.0386	342.0240	ESI+	10
613 Prothioconazole-desthio	120983-64-4	C14H15Cl2N3O	311.0592	312.0665	310.0519	ESI+	10
614 Prothiophos	34643-46-4	C11H15Cl2O2PS2	343.9628	344.9701	342.9555	ESI+	12
615 Prothoate	2275-18-5	C9H20NO3PS2	285.0622	286.0695	284.0549	ESI+	6
616 Pymetrozine	123312-89-0	C10H11N5O	217.0964	218.1037	216.0891	ESI+	4
617 Pyracarbolid	24691-76-7	C13H15NO2	217.1103	218.1176	216.1030	ESI+	3
618 Pyraclofos	89784-60-1	C14H18ClN2O3PS	360.0464	361.0537	359.0391	ESI+	12
619 Pyraclostrobin	175013-18-0	C19H18ClN3O4	387.0986	388.1059	386.0913	ESI+	11
620 Pyraflufen-ethyl	129630-19-9	C15H13Cl2F3N2O4	412.0204	413.0277	411.0131	ESI+	12
621 Pyrasulfotole	365400-11-9	C14H13F3N2O4S	362.0548	363.0621	361.0475	ESI+	9
622 Pyrazolynate	58011-68-0	C19H16Cl2N2O4S	438.0208	439.0281	437.0135	ESI+	2
623 Pyrazophos	13457-18-6	C14H20N3O5PS	373.0861	374.0934	372.0788	ESI+	12
624 Pyrazosulfuron-ethyl	93697-74-6	C14H18N6O7S	414.0958	415.1031	413.0885	ESI+	6
625 Pyrazoxyfen	71561-11-0	C20H16Cl2N2O3	402.0538	403.0611	401.0465	ESI+	12
626 Pyributicarb	88678-67-5	C18H22N2O2S	330.1402	331.1475	329.1329	ESI+	6
627 Pyridaben	96489-71-3	C19H25ClN2OS	364.1376	365.1449	363.1303	ESI+	12
628 Pyridalyl	179101-81-6	C18H14Cl4F3NO3	488.9680	489.9753	487.9607	ESI+	18
629 Pyridaphenthion	119-12-0	C14H17N2O4PS	340.0647	341.0720	339.0574	ESI+	6
630 Pyridate	55512-33-9	C19H23ClN2O2S	378.1169	379.1242	377.1096	ESI+	12
631 Pyrifenox	88283-41-4	C14H12Cl2N2O	294.0327	295.0400	293.0254	ESI+	8
632 Pyriftalid	135186-78-6	C15H14N2O4S	318.0674	319.0747	317.0601	ESI+	2
633 Pyrimethanil	53112-28-0	C12H13N3	199.1109	200.1182	198.1036	ESI+	6
634 Pyrimidifen	105779-78-0	C20H28ClN3O2	377.1870	378.1943	376.1797	ESI+	12
635 Pyriminobac-methyl (E)	136191-64-5	C17H19N3O6	361.1274	362.1347	360.1201	ESI+	6
636 Pyriproxyfen	95737-68-1	C20H19NO3	321.1365	322.1438	320.1292	ESI+	6
637 Pyroquilon	57369-32-1	C11H11NO	173.0841	174.0914	172.0768	ESI+	6
638 Pyroxsulam	422556-08-9	C14H13F3N6O5S	434.0620	435.0693	433.0547	ESI+	6
639 Quinalphos	13593-03-8	C12H15N2O3PS	298.0541	299.0614	297.0468	ESI+	6
640 Quinclorac	84087-01-4	C10H5Cl2NO2	240.9697	241.9770	239.9624	ESI+	11
641 Quinmerac	90717-03-6	C11H8ClNO2	221.0244	222.0317	220.0171	ESI+	12
642 Quinoclamine	2797-51-5	C10H6ClNO2	207.0087	208.0160	206.0014	ESI+	19
643 Quinoxifen	124495-18-7	C15H8Cl2FNO	306.9967	308.0040	305.9894	ESI+	12
644 Quizalofop (free acid)	76578-12-6	C17H13ClN2O4	344.0564	345.0637	343.0491	ESI+	24
645 Quizalofop-ethyl	76578-14-8	C19H17ClN2O4	372.0877	373.0950	371.0804	ESI+	12
646 Quizalofop-methyl	76578-13-7	C18H15ClN2O4	358.0720	359.0793	357.0647	ESI+	12
647 Quizalofop-P	94051-08-8	C17H13ClN2O4	344.0564	345.0637	343.0491	ESI+	9
648 Quizalofop-P-ethyl	100646-51-3	C19H17ClN2O4	372.0877	373.0950	371.0804	ESI+	12
649 Rabenzazole	40341-04-6	C12H12N4	212.1062	213.1135	211.0989	ESI+	12
650 Resmethrin	10453-86-8	C22H26O3	338.1882	339.1955	337.1809	ESI+	6
651 Rimsulfuron	122931-48-0	C14H17N5O7S2	431.0569	432.0642	430.0496	ESI+	9
652 Rotenone	83-79-4	C23H22O6	394.1416	395.1489	393.1343	ESI+	6
653 Saflufenacil	372137-35-4	C17H17ClF4N4O5S	500.0544	501.0617	499.0471	ESI+	8
654 Sebuthylazine	7286-69-3	C9H16ClN5	229.1094	230.1167	228.1021	ESI+	6
655 Sebuthylazine-desethyl	37019-18-4	C7H12ClN5	201.0781	202.0854	200.0708	ESI+	12
656 Secbumeton	26259-45-0	C10H19N5O	225.1590	226.1663	224.1517	ESI+	4
657 Sethoxydim	74051-80-2	C17H29NO3S	327.1868	328.1941	326.1795	ESI+	12
658 Siduron	1982-49-6	C14H20N2O	232.1576	233.1649	231.1503	ESI+	5
659 Silafluofen	105024-66-6	C25H29FO2Si	408.1921	409.1994	407.1848	ESI+	2
660 Silthiofam	175217-20-6	C13H21NOSSi	267.1113	268.1186	266.1040	ESI+	5
661 Simazine	122-34-9	C7H12ClN5	201.0781	202.0854	200.0708	ESI+	12
662 Simazine-2-hydroxy	2599-11-3	C7H13N5O	183.1120	184.1193	182.1047	ESI+	5
663 Simeconazole	149508-90-7	C14H20FN3OSi	293.1360	294.1433	292.1287	ESI+	6
664 Simetryn	1014-70-6	C8H15N5S	213.1048	214.1121	212.0975	ESI+	4
665 Spinetoram A	187166-40-1	C42H69NO10	747.4921	748.4994	746.4848	ESI+	2
666 Spinetoram B	187166-15-0	C43H69NO10	759.4921	760.4994	758.4848	ESI+	3
667 Spinosyn A	131929-60-7	C41H65NO10	731.4608	732.4681	730.4535	ESI+	4
668 Spinosyn D	131929-63-0	C42H67NO10	745.4765	746.4838	744.4692	ESI+	4
669 Spirodiclofen	148477-71-8	C21H24Cl2O4	410.1052	411.1125	409.0979	ESI+	11
670 Spiromesifen	283594-90-1	C23H30O4	370.2144	371.2217	369.2071	ESI+	4
671 Spirotetramat	203313-25-1	C21H27NO5	373.1889	374.1962	372.1816	ESI+	6
672 Spiroxamine	118134-30-8	C18H35NO2	297.2668	298.2741	296.2595	ESI+	6
673 Sulcotrione	99105-77-8	C14H13ClO5S	328.0172	329.0245	327.0099	ESI+	2
674 Sulfallate	95-06-7	C8H14ClNS2	223.0256	224.0329	222.0183	ESI+	9
675 Sulfaquinoxaline	59-40-5	C14H12N4O2S	300.0681	301.0754	299.0608	ESI+	6

Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
676 Sulfometuron-methyl	74222-97-2	C15H16N4O5S	364.0841	365.0914	363.0768	ESI+	6
677 Sulfosulfuron	141776-32-1	C16H18N6O7S2	470.0678	471.0751	469.0605	ESI+	6
678 Sulfotep	3689-24-5	C8H20O5P2S2	322.0227	323.0300	321.0154	ESI+	6
679 Sulprofos	35400-43-2	C12H19O2PS3	322.0285	323.0358	321.0212	ESI+	6
680 Tau-Fluvalinate	102851-06-9	C26H22ClF3N2O3	502.1271	503.1344	501.1198	ESI+	15
681 Tebuconazole	107534-96-3	C16H22ClN3O	307.1451	308.1524	306.1378	ESI+	10
682 Tebufenozide	112410-23-8	C22H28N2O2	352.2151	353.2224	351.2078	ESI+	6
683 Tebufenpyrad	119168-77-3	C18H24ClN3O	333.1608	334.1681	332.1535	ESI+	12
684 Tebupirimfos	96182-53-5	C13H23N2O3PS	318.1167	319.1240	317.1094	ESI+	6
685 Tebutam	35256-85-0	C15H23NO	233.1780	234.1853	232.1707	ESI+	6
686 Tebuthiuron	34014-18-1	C9H16N4O5	228.1045	229.1118	227.0972	ESI+	6
687 Teflubenzuron	83121-18-0	C14H6Cl2F4N2O2	379.9742	380.9815	378.9669	ESI-	12
688 Tembotrione	335104-84-2	C17H16ClF3O6S	440.0308	441.0381	439.0235	ESI+	12
689 Temphos	3383-96-8	C16H20O6P2S3	465.9897	466.9970	464.9824	ESI+	6
690 Teparloxydim	149979-41-9	C17H24ClNO4	341.1394	342.1467	340.1321	ESI+	8
691 Terbacil	5902-51-2	C9H13ClN2O2	216.0666	217.0739	215.0593	ESI-	10
692 Terbucarb	1918-11-2	C17H27NO2	277.2042	278.2115	276.1969	ESI+	12
693 Terbufos	13071-79-9	C9H21O2PS3	288.0441	289.0514	287.0368	ESI+	5
694 Terbufos-sulfone	56070-16-7	C9H21O4PS3	320.0340	321.0413	319.0267	ESI+	6
695 Terbufos-sulfoxide	10548-10-4	C9H21O3PS3	304.0390	305.0463	303.0317	ESI+	6
696 Terbumeton	33693-04-8	C10H19N5O	225.1590	226.1663	224.1517	ESI+	6
697 Terbumeton-desethyl	30125-64-5	C8H15N5O	197.1277	198.1350	196.1204	ESI+	3
698 Terbutylazine	5915-41-3	C9H16ClN5	229.1094	230.1167	228.1021	ESI+	10
699 Terbutylazine-2-hydroxy	66753-07-9	C9H17N5O	211.1433	212.1506	210.1360	ESI+	6
700 Terbutylazine-desethyl	30125-63-4	C7H12ClN5	201.0781	202.0854	200.0708	ESI+	12
701 Terbutryn	886-50-0	C10H19N5S	241.1361	242.1434	240.1288	ESI+	6
702 Tetrachlorvinphos (CVMP)	22248-79-9	C10H9Cl4O4P	363.8993	364.9066	362.8920	ESI+	16
703 Tetraconazole	112281-77-3	C13H11Cl2F4N3O	371.0215	372.0288	370.0142	ESI+	7
704 Tetraethylpyrophosphate	107-49-3	C8H20O7P2	290.0684	291.0757	289.0611	ESI+	6
705 Tetramethrin	7696-12-0	C19H25NO4	331.1784	332.1857	330.1711	ESI+	12
706 Thenylchlor	96491-05-3	C16H18ClNO2S	323.0747	324.0820	322.0674	ESI+	12
707 Thiabendazole	148-79-8	C10H7N3S	201.0361	202.0434	200.0288	ESI+	6
708 Thiocloprid	111988-49-9	C10H9ClN4S	252.0236	253.0309	251.0163	ESI+	6
709 Thiamethoxam	153719-23-4	C8H10ClN5O3S	291.0193	292.0266	290.0120	ESI+	12
710 Thiazafurion	25366-23-8	C6H7F3N4O5	240.0293	241.0366	239.0220	ESI+	6
711 Thiazopyr	117718-60-2	C16H17F5N2O2S	395.0931	397.1004	395.0858	ESI+	6
712 Thidiazuron	51707-55-2	C9H8N4O5	220.0419	221.0492	219.0346	ESI+	6
713 Thienencarbazone-methyl	317815-83-1	C12H14N4O7S2	390.0304	391.0377	389.0231	ESI+	3
714 Thifensulfuron-methyl	79277-27-3	C12H13N5O6S2	387.0307	388.0380	386.0234	ESI+	6
715 Thifluzamide	130000-40-7	C13H6Br2F6N2O2S	525.8421	526.8494	524.8348	ESI+	29
716 Thiobencarb	28249-77-6	C12H16ClNOS	257.0641	258.0714	256.0568	ESI+	11
717 Thiodicarb	59669-26-0	C10H18N4O4S3	354.0490	355.0563	353.0417	ESI+	6
718 Thiofanox	39196-18-4	C9H18N2O2S	218.1089	219.1162	217.1016	ESI+	2
719 Thiofanox-sulfone	39184-59-3	C9H18N2O4S	250.0987	251.1060	249.0914	ESI+	9
720 Thiofanox-sulfoxide	39184-27-5	C9H18N2O3S	234.1038	235.1111	233.0965	ESI+	12
721 Thiometon	640-15-3	C6H15O2PS3	245.9972	247.0045	244.9899	ESI+	2
722 Thionazin	297-97-2	C8H13N2O3PS	248.0384	249.0457	247.0311	ESI+	6
723 Thiophanate-ethyl	23564-06-9	C14H18N4O4S2	370.0769	371.0842	369.0696	ESI+	6
724 Thiophanate-methyl	23564-05-8	C12H14N4O4S2	342.0456	343.0529	341.0383	ESI+	6
725 Thiram	137-26-8	C6H12N2S4	239.9883	240.9956	238.9810	ESI+	6
726 Tolclofos-methyl	57018-04-9	C9H11Cl2O3PS	299.9544	300.9617	298.9471	ESI+	12
727 Tolyfluanid	731-27-1	C10H13Cl2FN2O2S2	345.9780	346.9853	344.9707	ESI+	22
728 Topramezone	210631-68-8	C16H17N3O5S	363.0889	364.0962	362.0816	ESI+	12
729 Tralkoxydim	87820-88-0	C20H27NO3	329.1991	330.2064	328.1918	ESI+	6
730 Tralomethrin	66841-25-6	C22H19Br4NO3	660.8098	661.8171	659.8025	ESI+	7
731 Triadimefon	43121-43-3	C14H16ClN3O2	293.0931	294.1004	292.0858	ESI+	12
732 Triadimenol	55219-65-3	C14H18ClN3O2	295.1088	296.1161	294.1015	ESI+	7
733 Tri-allate	2303-17-5	C10H16Cl3NO5	303.0018	304.0091	301.9945	ESI+	16
734 Triapenthenol	76608-88-3	C15H25N3O	263.1998	264.2071	262.1925	ESI+	12
735 Triasulfuron	82097-50-5	C14H16ClN5O5S	401.0561	402.0634	400.0488	ESI+	12
736 Triazamate	112143-82-5	C13H22N4O3S	314.1413	315.1486	313.1340	ESI+	4
737 Triazophos	24017-47-8	C12H16N3O3PS	313.0650	314.0723	312.0577	ESI+	6
738 Triazoxide	72459-58-6	C10H6ClN5O	247.0261	248.0334	246.0188	ESI+	11
739 Tribenuron-methyl	101200-48-0	C15H17N5O6S	395.0900	396.0973	394.0827	ESI+	5
740 Trichlorfon	52-68-6	C4H8Cl3O4P	255.9226	256.9299	254.9153	ESI+	10
741 Triclopyr	55335-06-3	C7H4Cl3NO3	254.9257	255.9330	253.9184	ESI-	2
742 Tricyclazole	41814-78-2	C9H7N3S	189.0361	190.0434	188.0288	ESI+	6
743 Tridemorph	81412-43-3	C19H39NO	297.3032	298.3105	296.2959	ESI+	6
744 Trietazine	1912-26-1	C9H16ClN5	229.1094	230.1167	228.1021	ESI+	6
745 Triethanolamine	102-71-6	C6H15NO3	149.1052	150.1125	148.0979	ESI+	6
746 Trifloxystrobin	141517-21-7	C20H19F3N2O4	408.1297	409.1370	407.1224	ESI+	6
747 Trifloxysulfuron	145099-21-4	C14H14F3N5O6S	437.0617	438.0690	436.0544	ESI+	9
748 Triflumizole	68694-11-1	C15H15ClF3N3O	345.0856	346.0929	344.0783	ESI+	9
749 Triflumizole Metabolite	131549-75-2	C12H14ClF3N2O	294.0747	295.0820	293.0674	ESI+	2
750 Triflumuron	64628-44-0	C15H10ClF3N2O3	358.0332	359.0405	357.0259	ESI+	8

	Compound	CAS	Formula	M	[M+H] ⁺	[M-H] ⁻	Ionisation Mode	MRM Transitions
751	Triflusulfuron-methyl	126535-15-7	C17H19F3N6O6S	492.1039	493.1112	491.0966	ESI+	8
752	Triforine	26644-46-2	C10H14Cl6N4O2	431.9248	432.9321	430.9175	ESI+	7
753	Trinexapac-ethyl	95266-40-3	C13H16O5	252.0998	253.1071	251.0925	ESI+	6
754	Triphenyl phosphate	115-86-6	C18H15O4P	326.0708	327.0781	325.0635	ESI+	6
755	Tris (2-chloro-1-(chloromethyl)ethyl) phosphate	13674-87-8	C9H15Cl6O4P	427.8839	428.8912	426.8766	ESI+	26
756	Triticonazole	131983-72-7	C17H20ClN3O	317.1295	318.1368	316.1222	ESI+	9
757	Tritosulfuron	142469-14-5	C13H9F6N5O4S	445.0279	446.0352	444.0206	ESI+	4
758	Valifenalate	283159-90-0	C19H27ClN2O5	398.1608	399.1681	397.1535	ESI+	16
759	Vamidothion	2275-23-2	C8H18NO4PS2	287.0415	288.0488	286.0342	ESI+	6
760	Vamidothion-sulfone	70898-34-9	C8H18NO6PS2	319.0313	320.0386	318.0240	ESI+	6
761	Vamidothion-sulfoxide	20300-00-9	C8H18NO5PS2	303.0364	304.0437	302.0291	ESI+	6
762	Vernolate	1929-77-7	C10H21NOS	203.1344	204.1417	202.1271	ESI+	5
763	Warfarin	81-81-2	C19H16O4	308.1049	309.1122	307.0976	ESI+	6
764	XMC (3, 5-xylyl methylcarbamate)	2655-14-3	C10H13NO2	179.0946	180.1019	178.0873	ESI+	12
765	Ziram	137-30-4	C6H12N2S4Zn	303.9175	304.9248	302.9102	ESI+	2
766	Zoxamide	156052-68-5	C14H16Cl3NO2	335.0247	336.0320	334.0174	ESI+	18

■ Further Information

Application News No.C136 describes the analysis of 646 pesticides in a single multi-residue method built using the Shimadzu Pesticide Library.

■ Scope and Legal Disclaimers

Whilst every effort has been made to ensure the accuracy of the Library, the method will need to be verified in a laboratory as conditions may differ marginally. The influence of sample matrices, extraction protocols, LC behaviour and technical experience may affect the performance of the LC/MS/MS analysis.

Shimadzu assumes no responsibility or contingent liability, including indirect or consequential damages, for any use to which the purchaser may put the referenced suppliers' products, or for any adverse circumstances arising therefrom.

First Edition: Jun. 2016



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2016

Shimadzu Corporation

www.shimadzu.com/an/

Application News

No. C136

Liquid Chromatography Mass Spectrometry

Expanding Capabilities in Multi-Residue Pesticide Analysis Using The LCMS-8060

David R. Baker¹, Laëtitia Fages², Eric Capodanno², Neil Loftus¹
¹Shimadzu Corporation, UK; ²Phytocontrol, France

■ Abstract

With an increasing global population, food security is increasingly under threat and there is a growing challenge for agriculture to produce more food, safely and more sustainably. The use of herbicides, insecticides, and fungicides reduce crop losses both before and after harvest, and increase crop yields. However, pesticide residues resulting from the use of plant protection products on crops may pose a risk to human health and require a legislative framework to monitor pesticide residues in food.

National programs for pesticide monitoring in the US, Europe and Japan have set Maximum Residue Levels (MRL's) or tolerance information (EPA) for pesticides in food products. A default value of 0.01 mg/kg is applied for MRL enforcement, which therefore requires highly sensitive and specific analytical technologies to monitor an increasing number of pesticides.

This application note describes the expanded capability of the LCMS-8060 to help accelerate method development workflows and support increased pesticide monitoring programs. Using the Shimadzu Pesticide MRM Library (the Library includes information on 766 certified reference materials) a single multi-residue LC/MS/MS method was developed for 646 pesticides (3 MRM transitions for over 99 % targeted pesticides resulting in 1,919 transitions in total, with a polarity switching time of 5 msec).

Keywords: Pesticides; food safety; LCMS-8060; Pesticide MRM Library, 776 compound library

■ Introduction

There are more than 1,000 pesticides used globally on soil and crops. With the ever increasing international trade of the food industry, regulatory bodies around the world have increased the number of regulated pesticides and the maximum residue levels (MRLs) allowed in food commodities. In the EU, regulation 396/2005/EC and its annexes set MRLs for over 500 pesticides in 370 food products.¹⁾ In the US, tolerances for more than 450 pesticides and other ingredients are established by the US EPA²⁾ and Japan's positive list system for agricultural chemical residues in foods contains MRLs for over 400 pesticides in various commodities.³⁾

National pesticide monitoring programs create new challenges for food safety laboratories as the number of pesticides required for analysis is increasing together with an expanded range of food products.

In this application paper we present the development of a LC-MS/MS method for screening and quantifying over 646 pesticides in a single method. The method

was quickly and efficiently set up using the Shimadzu Pesticide MRM Library. For each target pesticide analysis, up to 3 MRMs (Multiple Reaction Monitoring) transitions were imported from the library. 3 MRMs transitions provided additional data confidence in reporting results in comparison to the conventional 2 transitions used in most methods. As the LCMS-8060 has a high data acquisition speed 1,919 transitions were acquired using a polarity switching speed of 5 msec over a 10.5 minutes gradient elution.

To evaluate the method QuEChERS extracts of mint, tomato and apple were provided by a commercial laboratory as raw acetonitrile extracts and spiked with 646 pesticides (data is presented on the mint extract as it is the more complex sample matrix). The method was evaluated in matrix to ensure that the reporting limits were in agreement with recognised MRL's.

■ Experiment

Food extracts of mint, tomato and apple were supplied by Phytocontrol, France, following established QuEChERS protocols. Final extracts were prepared in acetonitrile without any dilution. Certified reference materials for the Shimadzu Pesticide MRM Library were obtained from ACSD, France as stock solutions. All solvents were of LCMS quality purchased from Sigma-Aldrich.

A six point calibration curve from 0.002 - 0.1 mg/kg (2 - 100 pg/μL) were generated using internal standard method. Two internal standards (Atrazine-d5 and Diuron-d6) were spiked in during the auto-sampler sequence for quantitation.

The robustness of the LCMS-8060 was assessed by peak area response for 646 pesticides spiked into mint, tomato and apple matrix extracts at 0.05 mg/kg.

■ LC/MS/MS method development

The Shimadzu Pesticide MRM Library has 766 pesticides in its database (Application News No. C135). For each pesticide several MRM's are included in the database and in this analysis the default value used was 3 MRM's. For this method, 1,919 transitions were selected in both positive and negative ionisation mode using a switching time of 5 msec (1,819 MRM transitions were in positive mode and 100 MRM transitions in negative mode).

To optimize ion source conditions (for example, DL temperature, interface temperature, heating block temperature, heating gas flow, drying gas flow and nebulizer gas flow) the interface setting software was used. This tool provides an optimized response for all compounds.

Table 1 LC and MS/MS Acquisition Parameters

Liquid chromatography		Mass spectrometry	
UHPLC	Nexera LC system	LC/MS/MS	LCMS-8060
Analytical column	Restek Raptor Biphenyl (2.1 mm I.D. × 100 mm L., 2.7 μm)	Ionisation mode	Heated electrospray
Column temperature	35 °C	Polarity switching time	5 msec
Flow rate	0.4 mL/min	Pause time	1 msec
Solvent A	2 mmol/L ammonium formate + 0.002 % formic acid - Water	Total MRM transitions	1,919 (1,819 positive; 100 negative)
Solvent B	2 mmol/L ammonium formate + 0.002 % formic acid - Methanol	MRM Dwell	4 msec (target ion); 1 msec (reference ion)
Binary Gradient B.Conc.	3 % (0 min) - 10 % (1.00 min) - 55 % (3.00 min) - 100 % (10.50 - 12.00 min) - 3 % (12.01 - 15.00 min)	Interface temperature	350 °C
Injection volume	2 μL sample (plus 40 μL water)	Heating block	300 °C
		Desolvation line	150 °C
		Heating gas	10 L/min
		Drying gas	10 L/min
		Nebulizer gas	3 L/min




Fig. 1 MRM chromatograms of 646 pesticides spiked into a mint extract at 0.01 mg/kg (Up to 3 MRMs per compound and 5 msec polarity switching time).




Fig. 2 MRM chromatograms for pesticides most commonly detected in plant products listed in the 2015 European Food Safety Journal. In this report, residues exceeding the legal limits were related to 58 different pesticides. Compounds such as boscalid, chlorpyrifos, cyprodinil, fenhexamid, fludioxonil, pyraclostrobin and tebuconazole (highlighted in the MRM chromatogram) are some of the most frequently detected compounds present in more than 4 % of the samples analyzed.

The MRM chromatograms show the response to each pesticide spiked into a food matrix at the default MRL of 0.01 mg/kg.

Results and Discussion**Shimadzu Pesticide MRM Library**

(Application News No. C135)

A flexible tool for expanding capabilities in pesticide monitoring programs

The Pesticide MRM Library has been created using 766 certified reference materials and is designed to help accelerate method development and compound management.

The library contains an average of 8 optimized MRM transitions for each compound (including positive and negative ion modes). In total, more than 6,000 MRM transitions are held within the 766 compound library. The library itself documents CAS#, formula, activity, mono-isotopic mass and adduct masses, rank of MRM transitions, synonyms, InChI, InChIKey, compound names translation (Japanese and Chinese) and links to websites offering further information (for example; alanwood.net, PAN pesticide database, Chemical Book, ChemSpider).

The library also serves as a powerful data repository for reporting and checking pesticide data sources.

Creating flexible pesticide monitoring methods**Building a new LC/MS/MS method**

To create new pesticide LC/MS/MS methods the user simply needs to select the target compounds from the library, identify the required number of MRMs for each compound and confirm the analytical column for the analysis. (The new method can be used to expand current capabilities or to create focused methods with a limited number of pesticides). The new method is simply imported into LabSolutions.

As the LCMS-8060 has a high data acquisition speed of 30,000 u/sec, high sensitivity and a polarity switching speed of 5 msec, the capabilities of the library can be expanded to meet the future needs of any laboratory.

Expanded capability of the LCMS-8060

The LCMS-8060 has a data acquisition speed of 30,000 u/sec which creates new opportunities for expanding compound lists.

As one example, between 6.45 and 6.60 minutes 25 pesticide compounds elute (Fig. 3). Even with high data density acquisitions the average variation in peak area response was less than 3 %RSD (varying between 1.1 - 5.9 %RSD).




Fig. 3 The LCMS-8060 can acquire MRM data at a high speeds and enables precise quantitation even with high data density. Between 6.45 and 6.60 minutes 25 compounds were monitored (Table 2).

Table 2 Peak area variation (%RSD; n=6) for 25 pesticides eluting over a nine-second time window (6.45 - 6.60 minutes) spiked into a mint matrix extract at the reporting limit of 0.01 mg/kg.

Compound Name	CAS number	Formula	M	Polarity	MRM Quantitation Ion	RT	Average Peak Area	%RSD (n=6)
Trinexapac-ethyl	95266-40-3	C ₁₃ H ₁₆ O ₅	252.0998	+	252.90 > 69.05	6.45	1,780,015	3.1
Iprovalicarb	140923-17-7	C ₁₈ H ₂₈ N ₂ O ₃	320.2100	+	321.20 > 119.15	6.46	1,442,486	2.8
Dodemorph	1593-77-7	C ₁₈ H ₃₅ N ₂ O	281.2719	+	282.30 > 116.15	6.47	658,920	4.2
Fluopyram	658066-35-4	C ₁₆ H ₁₁ ClF ₆ N ₂ O	396.0464	+	397.00 > 145.00	6.47	2,439,146	1.9
Flutolanil	66332-96-5	C ₁₇ H ₁₆ F ₃ N ₂ O ₂	323.1133	+	324.10 > 242.00	6.48	3,372,285	2.7
Trifloxysulfuron	145099-21-4	C ₁₄ H ₁₄ F ₃ N ₅ O ₆ S	437.0617	+	438.00 > 182.15	6.48	1,822,340	2.5
Azaconazole	60207-31-0	C ₁₂ H ₁₁ Cl ₂ N ₃ O ₂	299.0228	+	300.00 > 159.00	6.50	1,580,445	2.0
Terbutryn	886-50-0	C ₁₀ H ₁₉ N ₅ S	241.1361	+	242.10 > 157.95	6.50	755,446	3.4
Prometryn	7287-19-6	C ₁₀ H ₁₉ N ₅ S	241.1361	+	242.10 > 158.00	6.50	1,300,193	2.6
Azimsulfuron	120162-55-2	C ₁₃ H ₁₆ N ₁₀ O ₅ S	424.1026	+	425.10 > 182.10	6.50	2,498,050	1.8
Metominostrobin	133408-50-1	C ₁₆ H ₁₆ N ₂ O ₃	284.1161	+	285.10 > 193.95	6.51	2,929,500	1.7
Thifluzamide	130000-40-7	C ₁₃ H ₆ Br ₂ F ₆ N ₂ O ₂ S	525.8421	+	528.60 > 148.05	6.51	193,982	5.9
Nicarbazin	330-95-0	C ₁₃ H ₁₀ N ₄ O ₅	302.0651	-	301.10 > 137.15	6.52	973,101	2.6
Bromobutide	74712-19-9	C ₁₅ H ₂₂ BrNO	311.0885	+	312.10 > 194.10	6.53	1,829,781	2.1
Saflufenacil	372137-35-4	C ₁₇ H ₁₇ ClF ₄ N ₄ O ₅ S	500.0544	+	501.00 > 198.00	6.53	465,224	2.3
Cyproconazole	94361-06-5	C ₁₅ H ₁₈ CIN ₃ O	291.1138	+	292.10 > 70.05	6.54	1,174,967	1.7
Clomazone	81777-89-1	C ₁₂ H ₁₄ CINO ₂	239.0713	+	239.90 > 125.00	6.54	3,409,656	1.7
Fensulfothion	115-90-2	C ₁₁ H ₁₇ O ₄ PS ₂	308.0306	+	309.00 > 281.00	6.54	4,267,514	1.4
Oxasulfuron	144651-06-9	C ₁₇ H ₁₈ N ₄ O ₆ S	406.0947	+	407.10 > 150.15	6.54	2,911,533	1.1
Rimsulfuron	122931-48-0	C ₁₄ H ₁₇ N ₅ O ₇ S ₂	431.0569	+	432.00 > 182.00	6.55	4,722,065	1.8
Fenthion-oxon	6552-12-1	C ₁₀ H ₁₅ O ₄ PS	262.0429	+	263.10 > 231.00	6.55	3,075,195	1.4
Nitrothal-isopropyl	10552-74-6	C ₁₄ H ₁₆ NO ₆ Na	317.0875	+	295.10 > 230.95	6.56	2,199,581	3.0
Chlorantraniliprole	500008-45-7	C ₁₈ H ₁₄ BrCl ₂ N ₅ O ₂	480.9708	+	483.90 > 452.90	6.57	2,407,025	2.7
Fipronil-sulfone	120068-36-2	C ₁₂ H ₄ Cl ₂ F ₆ N ₄ O ₂ S	451.9336	-	451.00 > 414.90	6.57	2,843,708	2.0
Valifenalate	283159-90-0	C ₁₉ H ₂₇ CIN ₂ O ₅	398.1608	+	399.20 > 155.00	6.59	3,845,335	1.9

Final method performance for 646 pesticides

In order to test the performance of the developed method, linearity, repeatability and longer term robustness were assessed for all 646 pesticides.

Linearity

Linearity was assessed over a six point calibration curve from 0.002 - 0.1 mg/kg (2 - 100 pg/μL). All 646 pesticides achieved excellent R^2 values greater than 0.99 in both tomato and mint spiked extracts with typical values greater than 0.996. Calibration curves were generated using a linear curve fit type and 1/C weighting. Typical calibration curve data is presented below in Fig. 4.




Fig. 4 Calibration curves for selected pesticides spiked into a mint matrix extract in the range 0.002 - 0.1 mg/kg. The quantitation MRM chromatogram is shown in black (qualifier ion MRM chromatograms are shown in red and blue).

Repeatability

To assess the robustness of the system and the developed method during routine analysis, repeat injections of a mint matrix sample spiked with 646 pesticides at 0.05 mg/kg, were analyzed over a 24 hour period.
The results for selected compounds are displayed below in Fig. 5.

Compounds were selected throughout the run at equidistant points (closest elution points to 3, 4, 5, 6, 7, 8, 9 and 10 minutes), including positive and negative ion detection, (Table 3).
The peak area variance was less than 5.7 % for all pesticides measured.




Fig. 5 Peak area response for several pesticides following 100 repeat injections of a 0.05 mg/kg spiked into mint matrix extract.

Table 3 Peak area variance for selected following the repeated injection of a 0.05 mg/kg spiked into mint matrix extract (number of sample replicates was 100; the analysis sequence was 24 hours).

Compound Name	CAS Number	Formula	M	Polarity	MRM Quantitation Ion	RT (mins)	Average Peak Area	%RSD (n=100)
Butocarboxim-sulfoxide	34681-24-8	C7H14N2O3S	206.0725	+	207.10 > 75.10	3.042	1,220,391	2.6
Thiofanox-sulfone	39184-59-3	C9H18N2O4S	250.0987	+	268.10 > 57.00	4.001	442,724	5.7
Monolinuron	1746-81-2	C9H11ClN2O2	214.0509	+	215.10 > 99.10	4.985	2,904,116	3.7
Probenazole	27605-76-1	C10H9NO3S	223.0303	+	224.00 > 41.05	5.995	1,145,189	3.5
Dipropetryn	4147-51-7	C11H21N5S	255.1518	+	256.20 > 144.05	6.999	3,289,597	3.4
Pyraflufen-ethyl	129630-19-9	C15H13Cl2F3N2O4	412.0204	+	413.00 > 339.00	8.004	3,653,333	3.5
Emamectin B1a	138511-97-4	C56H81NO15	1007.5606	+	886.40 > 158.20	9.008	3,109,562	4.5
Pyridalyl	179101-81-6	C18H14Cl4F3NO3	488.9680	-	491.90 > 109.05	10.171	1,579,422	5.0

Response to differing matrices

One of the major challenges in the quantitative LC/MS/MS analysis for pesticides in food is that compound and matrix-dependent response suppression or enhancement may occur. Although matrix effects can affect the peak area response between different food types following a QuEChERS extraction protocol, the peak area variance should be minimized within a single matrix.

Food extracts of apple, mint and tomato following QuEChERS extraction were spiked with 646 pesticides at 0.05 mg/kg and were repeatedly injected on the LCMS-8060 (n=100 repeat injections for each matrix; 300 injections in the same batch sequence). Fig. 6 shows the response for 3 selected pesticides analyzed in a single batch sequence corresponding to a 72 hour analysis sequence. Within a matrix, variance was less than 5.9 %RSD for all compounds.

Although the absolute peak area changes with different food matrices, the response between injection 1 and injection 100 for 2 pesticides (probenazole and dipropetryn) within a single matrix has a variance less than 5.7 %RSD.




Fig. 6 Peak area response for three pesticides spiked into apple, mint and tomato matrix extracts at 0.05 mg/kg over 72 hours. As in Fig. 5, compounds were selected to reflect peak area response throughout the chromatographic run (Table 3).




Fig. 7 MRM chromatograms for probenazole (RT 5.995 minutes) and dipropetryn (RT 6.999 minutes) for injection 1 and injection 100 spiked into apple, mint and tomato matrix extracts. The extracts were spiked at 0.05 mg/kg and analyzed over 72 hours.

Reducing matrix effects by extensively diluting the sample

The need to test for more pesticides in a wider range of samples at high sensitivity is very challenging as matrix effects from the sample extraction will influence both ion suppression and enhancement. Ion suppression can lead to errors in the detection capability, accuracy and precision of the method.

To reduce the effect of interfering compounds in the quantitation of complex samples extensive sample dilution is now widely used in routine analysis. It is an approach which is simple to build into multi-residue extraction methods and is cost effective.

This approach leads to greater robustness as a consequence of a reduced sample injection in the LC/MS/MS, higher data quality and increased instrument uptime.

Fig. 8 shows the results of diluting a matrix sample spiked at 0.005 mg/kg with dilution factors of 1:5, 1:10, 1:20, 1:50 and 1:100.

As matrix effects can be both significant and variable for different compounds Table 4 shows recovery data for a series of pesticides diluted from 0 to a dilution factor of 1:100.

Matrix suppression was reduced for most compounds when the sample was diluted 1:10 with recoveries in the range of 70 - 120 % with an associated repeatability RSDr ≤ 20 %. Relative standard deviations in relation to the mean values were typically less than 10 %.

Diluting the sample by a factor of 20 or 50 resulted in acceptable signal suppression from the matrix.




Fig. 8 MRM chromatograms for 3 selected compounds spiked into a mint extract at 0.005 mg/kg and diluted 1:5, 1:10, 1:20, 1:50 and 1:100 with water.

Table 4 Diluting a sample matrix extract spiked with 0.005 mg/kg with water reduced matrix ion suppression.

				Dilution series					
Compound	CAS	Formula	M	0	1:5	1:10	1:20	1:50	1:100
				Recovery					
Bentazone	25057-89-0	C10H12N2O3S	240.0569	32.1	44.6	65.5	72.7	91.7	98.1
Demeton-S-methyl-sulfone	17040-19-6	C6H15O5PS2	262.0099	51.1	78.5	89.6	91.1	114.2	116.8
Dimethoate	60-51-5	C5H12NO3PS2	228.9996	36.2	65.3	88.5	92.2	92.4	94.2
Isocarbamid	30979-48-7	C8H15N3O2	185.1164	28.8	57.1	81.8	98.7	102.5	96.4
Vamidothion	2275-23-2	C8H18NO4PS2	287.0415	53.6	76.3	98.2	98.5	101.5	114.1
Thiazafluron	25366-23-8	C6H7F3N4OS	240.0293	32.8	62.9	80.5	84.2	87.1	97.4
Demeton-S-methyl	919-86-8	C6H15O3PS2	230.0200	57.8	82.1	93.1	87.6	108.5	102.4
Sebuthylazine	7286-69-3	C9H16ClN5	229.1094	28.7	53.3	69.8	79.8	88.5	95.8
Flutriafol	76674-21-0	C16H13F2N3O	301.1027	27.3	46.1	71.4	76.1	81.8	87.3
Furametpyr	123572-88-3	C17H20ClN3O2	333.1244	48.3	69.8	86.9	86.2	97.6	101.9
Fenobucarb	3766-81-2	C12H17NO2	207.1259	60.9	79.2	100.7	96.1	102.8	103.9
Benodanil	15310-01-7	C13H10INO	322.9807	50.9	69.8	86.3	96.5	102.4	94.8
Terbuthylazine	5915-41-3	C9H16ClN5	229.1094	50.4	66.6	83.2	87.2	89.8	91.0
Dimethachlor	50563-36-5	C13H18ClNO2	255.1026	75.1	86.1	106.0	107.1	106.2	108.0
Dimethenamid	87674-68-8	C12H18ClNO2S	275.0747	72.6	84.9	102.9	100.0	103.6	97.3
Furalaxyl	57646-30-7	C17H19NO4	301.1314	82.2	89.1	106.6	108.6	106.2	102.4
Bixafen	581809-46-3	C18H12Cl2F3N3O	413.0310	66.8	79.3	99.0	95.6	103.7	97.1
Triflumuron	64628-44-0	C15H10ClF3N2O3	358.0332	54.2	71.8	95.5	84.9	95.3	101.7
Epoxiconazole	133855-98-8	C17H13ClFN3O	329.0731	61.6	77.2	98.8	95.3	90.0	101.2
Teflubenzuron	83121-18-0	C14H6Cl2F4N2O2	379.9742	41.8	50.9	80.1	86.8	100.0	97.7

■ Conclusion

A fast, selective and highly sensitive method has been developed for the quantitation of 646 pesticides using a single method with 1,919 transitions (corresponding to up to 3 MRM transitions per compound) and a LC gradient time of only 10.5 minutes.

As the LCMS-8060 has a rapid polarity switching time of 5 msec, the single multi-residue LC/MS/MS method supported the analysis of 34 pesticides in negative ion mode and 612 compounds in positive ion mode.

The enhanced performance and higher sensitivity of the LCMS-8060 has created new opportunities in sample dilution to reduce ion signal suppression and matrix effects. For most compounds a dilution factor of 1:20 or 1:50 was sufficient to provide recoveries in the range 70 - 120 %.



■ References

1. Commission Regulation (EC). 2005. No 396/2005 of the European Parliament and of the Council, maximum residue levels of pesticides in or on food and feed of plant and animal origin. Official Journal of the European Union, L 70: 1-16
2. US Environmental Protection Agency, Electronic code of federal regulation: Title 40: Part 180 - tolerances and exemptions for pesticide chemical residues in food. http://www.ecfr.gov/cgi-bin/text-idx?c=ecfr&tpl=/ecfrbrowse/Title40/40cfr180_main_02.tpl
3. Japanese Ministry of Health, Labour and Welfare, Department of Food Safety. 2006. Director Notice about Analytical Methods for Residual Compositional Substances of Agricultural Chemicals, Feed Additives, and Veterinary Drugs in Food (Syoku-An No. 0124001 January 24, 2005; amendments May 26, 2006).

First Edition: Jun. 2016



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2016

Application News

No. C154

Liquid Chromatography Mass Spectrometry

Applying 'MRM Spectrum Mode' and Library Searching for Enhanced Reporting Confidence in Routine Pesticide Residue Analysis

David R. Baker¹, Chris Titman¹, Jonathan Horner², Neil Loftus¹

¹Shimadzu Corporation, UK; ²Scientific Analysis Laboratories, UK

Abstract

To help reduce the incidence of false positive and false negative reporting in pesticide residue monitoring routine multiple-reaction monitoring (MRM) methods have been enhanced to monitor a higher number of fragment ion transitions to increase specificity and reporting confidence. In this workflow, typically 6-10 fragment ion transitions were monitored for each target pesticide as opposed to a conventional approach using 2-3 fragment ions. By acquiring a high number of fragment ion transitions, each target pesticide had a corresponding fragmentation spectra which could be used in routine library searching and compound verification using reference library match scores. This 'MRM Spectrum Mode' was applied to quantify and identify 193 pesticides using 1,291 MRM transitions without compromising limits of detection, linearity or repeatability.

Introduction

Multiple Reaction Monitoring (MRM) based LC-MS/MS techniques are widely used on triple quadrupole platforms for targeted quantitation as a result of high selectivity, sensitivity and robustness. In a regulated environment such as food safety there is a growing need to enhance the capability in routine monitoring programs by increasing the number of pesticides measured in a single analysis and at the same time delivering the highest confidence in compound identification to reduce false detect reporting. For pesticide analysis in the EU, identification criteria in SANTE/11945/2015 requires the retention time and the ion ratio from at least 2 MRM transitions to be within acceptable tolerance limits.¹ However, even applying this criteria it is well reported that false positives can occur in certain pesticide/commodity combinations.²⁻⁴

To reduce false negative and false positive reporting a higher number of MRM transitions were used for each target pesticide to increase the level of confidence in assay specificity. The number of fragment ion transitions monitored for each target pesticide was dependent upon the chemical structure with typically between 6-10 fragment ions for each compound. MRM Spectrum mode combines conventional MRM quantitation with the generation of a high quality MRM product ion spectrum which can be used in routine library searching and compound verification and identification.

In this application paper we present the development of a method for 193 pesticides, with 1,291 MRM transitions, and a 15 minute cycle time. In order to acquire this number of MRM transitions using a short run time a 3 msec dwell time was applied to each MRM transition and a 5 msec polarity switch was used. On average 7 MRM transitions were applied to each compound. The method was quickly set up using the Shimadzu Pesticide Method Package, a data base with more than 750 pesticides and over 6,000 MRM transitions designed to accelerate method set-up and help compound verification. MRM Spectrum mode was also compared to a conventional pesticide monitoring method with 2 MRMs per compound (386 MRMs in total) in order to assess the effect on data quality when adding additional MRM transitions to the method. Several different food commodities were analysed with varying complexity (turmeric, plum, peppermint, parsnip, cherry, lime, pumpkin, tomato, potato). Data was processed using LabSolutions Insight software which provides automated library searching of target MRM spectrum.




Fig. 1 Using a higher number of fragment ions in MRM data acquisition increases the specificity of detection and reduces false negative and false positive reporting. In the case of linuron, 9 precursor-fragment ion transitions were used to increase confidence in assay specificity. There is no compromise in data quality between methods despite a higher number of fragment ions monitored. Signal intensity, linearity, reproducibility are in good agreement between both methods.

Experimental

Pesticide spiked samples, extracted using established QuEChERS based methods, were provided by Scientific Analysis Laboratories, UK. In order to test the performance of the MRM Spectrum Mode database and library searching a number of matrices were tested including turmeric, plum, peppermint, parsnip, cherry, lime, pumpkin, tomato and potato. Final extracts were prepared in acetonitrile without any dilution and directly injected into the LC-MS/MS. A water co-injection method, performed automatically in the auto-sampler, was used to improve early eluting peak shapes in addition to a sub 2 micron particle size column to improve peak capacity (Table 1).

Calibration curves were prepared in the range 0.01 to 0.2 mg/kg. Repeatability of the method was tested using avocado matrix at 0.1 mg/kg. In the final method samples were analysed in ESI +/- using a polarity switching time of 5 msec.

On average 7 MRM transitions were applied to each compound, with more than 10 MRM transitions applied to 34 compounds. All MRM transitions were acquired throughout the MRM window without the need for triggering thresholds. The method includes a total of 1,291 MRM transitions for 193 pesticides in a run time of only 15 minutes. A dwell time of 3 msec was applied to every MRM transition. In order to evaluate the data quality from the MRM Spectrum Mode method, the same method was set up with 2 MRMs applied to each compound (386 MRMs in total) using the same acquisition method (Table 2).

LabSolutions software was used to automatically optimize the fragmentation for all pesticides and generate a MRM Spectrum mode method. The MRM Spectrum Mode method for library searching and compound verification could be simply and quickly set up using the Shimadzu pesticide database. This database contains more than 6,000 MRM transitions for over 750 pesticides.

LabSolutions Insight v3.0 software was used to review quantitative data and MRM Spectrum mode library searching with advanced filtering tools to review by exception and to reduce false detect reporting.

Table 1 LC acquisition parameters

Liquid chromatography		
UHPLC	Nexera LC system	
Analytical column	HSS T3 (100 × 2.1, 1.7 µm)	
Column temperature	40 °C	
Flow rate	0.4 mL/minute	
Solvent A	5 mmol/L ammonium formate and 0.004 % formic acid	
Solvent B	5 mmol/L ammonium formate and 0.004 % formic acid in methanol	
Binary Gradient	Time (mins)	%B
	1.50	35
	11.50	100
	13.00	100
	13.01	3
	15.00	Stop
Injection volume	0.1 µL (plus 30 µL water)	

Table 2 MS/MS methods used to acquire data in MRM Spectrum Mode and a conventional MRM method with 2 MRM transitions per compound. As part of the comparative study, the same LC conditions were used for both methods.

LC-MS/MS Mass spectrometry	MRM Spectrum Mode: generating library searchable spectra	2 MRM method
Target number of compounds	193	193
Total number of MRM transitions	1,291 transitions (1,229 in ESI+ and 62 in ESI-)	386 (374 in ESI+ and 12 in ESI-)
Pause time/dwell time	1 msec./3 msec.	1 msec./3 msec.
Ionisation mode	ESI +/-	ESI +/-
Polarity switching time	5 msec	5 msec
Interface temperature	350 °C	350 °C
Heat block temperature	300 °C	300 °C
Desolvation line temperature	150 °C	150 °C
Nebulising gas	3 L/min	3 L/min
Heating gas	10 L/min	10 L/min
Drying gas	10 L/min	10 L/min

Results and Discussion

In developing monitoring programs for chemical contamination methods are designed to determine a list of known analytes with a focus on delivering a rapid, cost-effective analysis that generates no false-negative or false-positive results. Guidelines for compound identification have been published by the EU in directive SANTE/11945/2015. This identification criteria requires at least two MRM transitions with an ion ratio and retention time within defined tolerance limits.

To help reduce false detect reporting in pesticide monitoring programs, a MRM method was developed with a higher number of MRM transitions for each target pesticide to increase the level of confidence in assay specificity. By combining multiple MRM transitions for a compound into a product ion spectrum, pesticide identification can be verified and confirmed against a MS/MS reference spectral library. Using MRM Spectrum mode can help markedly reduce false detect reporting without affecting the data quality for optimized quantitation or identification.

Fig. 2, shows the MRM chromatogram for all 193 pesticides spiked at 0.010 mg/kg measured with MRM Spectrum mode. Using this mode 1,291 MRM transitions were measured for 193 pesticides. Despite the high data density acquired with MRM Spectrum Mode (for example, 151 MRM transitions were registered in the same time window during the analysis, see Fig. 3) sensitivity was not affected by the high data acquisition rate.

Method performance




Fig. 2 Histogram showing the number of MRM transitions monitored at each time point and chromatogram showing all 193 target compounds. The highest number of overlapping MRM's acquired was 151. Even at such a high data sampling rate the response was in agreement with a conventional 2 MRM method with peak area variation less than 5.2% (n=5). This data is displayed below in more detail, Fig. 3.

Table 3 Between 8.80 mins and 9.30 mins 151 MRM transitions in both positive and negative ion were monitored. Peak area repeatability for the 22 compounds eluting in this time period is shown below.

	Ret. Time	# MRMs	Polarity	Peak Area %RSD (n=5)
Dichlofluanid	8.80	6	ESI+	2.2
Dichlofluanid 2	8.80	6	ESI+	3.4
Dichlofluanid 1	8.80	5	ESI+	2.6
Fluoxastrobin	8.82	12	ESI+	2.0
Fenhexamid	8.83	11	ESI+	2.2
Iprovalicarb	8.88	6	ESI+	2.3
Spirotetramat	8.89	6	ESI+	2.6
Azinphos-ethyl	8.90	5	ESI+	3.1
Chromafenozide	8.91	5	ESI+	3.2
Triticonazole	8.93	5	ESI+	2.1
Cyazofamid	9.01	5	ESI+	2.1
Prothioconazole	9.07	10	ESI+	1.9
desthio				
Diflubenzuron	9.09	4	ESI+	2.0
Pyriphenox	9.11	8	ESI+	2.0
Dodemorph	9.17	6	ESI+	2.1
Fenoxycarb	9.17	6	ESI+	2.0
Rotenone	9.17	6	ESI+	2.4
Fipronil	9.20	10	ESI-	5.2
Bixafen	9.25	8	ESI-	2.8
Tebufenozide	9.27	6	ESI+	3.9
Bensulide	9.27	6	ESI+	2.6
Neburon	9.30	9	ESI+	1.7
		Total MRM's 151	Average 2.6 %RSD	




Fig. 3 Between 8.80 mins and 9.30 mins 151 MRM transitions in both positive and negative ion were monitored. During this time period 22 target pesticides eluted with a peak area variation less than 5.2 % RSD. Data was acquired in an avocado sample matrix at a concentration of 0.1 mg/kg.

Method performance




Fig. 4 MRM chromatograms for ethirimol (positive ion) and lufenuron (negative ion) acquired using a conventional 2 fragment ion MRM method and compared to a method with a higher number of precursor-fragment ions to increase confidence in assay specificity and reporting.

Despite acquiring a higher number of MRM transitions the library searchable MRM approach (acquiring 1,291 transitions in a single method) results in the same signal intensity compared to a conventional 2 fragment ion MRM method (acquiring 386 MRM transitions in a single method). The repeatability for each MRM method was evaluated by repeatedly injecting (n=5) an avocado extract corresponding to a concentration of 0.1 mg/kg. In each MRM method the %RSD was less than 3.5% for both compounds.

To minimize the possibility of false positive and false negative reporting LC-MS/MS methods were developed with a high number of MRM transitions for each pesticide. The performance of this approach was compared with a conventional MRM method monitoring 2 transitions for each pesticide.

In Fig. 4, the MRM chromatograms for 2 compounds, ethirimol and lufenuron, are shown for the same sample extract acquired using different MRM methods (the sample is avocado spiked at 0.1 mg/kg). The MRM chromatograms show un-smoothed data and are scaled to the same signal intensity for each compound. Ethirimol and lufenuron elute at 7.02 and 10.75 mins corresponding to time windows of high data density with more than one hundred MRM transitions monitored in the same time segment. However, regardless of the high number of fragment ions monitored, the absolute signal intensity for both approach's is near identical in positive and negative ion mode.

Fig. 5 shows the correlation between the peak areas for all pesticides measured using 2 different MRM methods. The linear regression curve shows a good agreement between the peak areas measured for all pesticides spiked into sample matrix with a slope value near unity and an intercept near zero.




Fig. 5 Absolute peak area response for all 193 pesticides acquired using a conventional MRM method with 386 transitions compared to a MRM method with 1,291 transitions designed for library searchable verification. Both approaches result in near identical peak areas regardless of the number of fragment ions used to verify and identify each pesticide.

■ Spectrum based identification

In this study, the number of qualifier fragment ion transitions was increased for each pesticide and the combined transitions were used to create a MRM product ion spectrum. This product ion spectrum derived from MRM acquisitions was used in conventional library matching routines comparing against a reference spectrum to generate a similarity score.

In Fig. 6, demeton-S-methyl sulphone was to highlight library matching in different matrices including cumin, potato, mucuna pruriens powder, tomato, black pepper, peppermint tea and turmeric. Even in the presence of complex spice matrices the library matching approach identified demeton-S-methyl sulphone with a high similarity score and a high degree of confidence for data reporting.




Fig. 6 MRM spectrum identification in different matrices for demeton-S-methyl sulphone

■ Spectrum based identification

To increase the confidence in reporting results the number of qualifier transitions was increased for each pesticide and the combined MRM transitions were used to create a product ion spectrum. This MRM product ion spectrum can then be automatically compared against a reference spectrum to generate a product ion spectrum match score using conventional library matching.

Fig. 7 highlights the advantage of using a library searchable fragment ion spectrum in identifying and quantifying desmedipham and phenmedipham. Both desmedipham and phenmedipham share several common fragment ions and have similar retention times. Using MRM Spectrum Mode and comparing to a library searchable spectra, both desmedipham and phenmedipham are positively identified (fragment ions at m/z 154 and 182 are absent in product ion spectrum for phenmedipham).




Fig. 7 MRM chromatogram for desmedipham and phenmedipham spiked into a cumin extract at 0.1 mg/kg. As phenmedipham shares common transitions and elutes at a similar retention time as desmedipham the MRM spectrum can be used to distinguish between both pesticides to avoid false positive reporting.

Quantitation

As one example, carbendazim was spiked into a matrix at three different concentration levels. In Fig. 8, all MRM transitions were detected even at the reporting level of 0.010mg/kg with a signal to noise for all fragment ion transitions greater than 9. The response was linear for all transitions throughout the calibration range (0.010-0.200mg/kg) as shown Fig. 9.




Fig. 8 By applying a range of collision energies to carbendazim 12 precursor-fragment ions are generated. MRM 192.10>159.95 was used in generating sensitive and robust quantitation whilst the product ion spectrum using all 12 fragment ions was used in confirming peak identification.




Fig. 9 Calibration curve for carbendazim using the optimized quantitation ion transition (MRM 192.10>159.95). The response was linear for all calibration and QC samples. All 12 fragment ions were above a signal to noise ratio of 10 even at the reporting level of 0.010mg/kg.

The limit on the number of MRM transitions used to generate a product ion spectrum is dependent on the chemical structure of the pesticide molecule. In the case of carbendazim, several bonds could be broken using collision energies between 10-60V resulting in a product ion spectrum of 12 fragment ions. The product ion spectrum can then be used for library search and analyte confirmation as shown in Fig. 10. For each calibration level ranging from 0.010-0.200mg/kg the library similarity score was greater than 99 confidently confirming the target analyte. The advantage of this technique is that library searchable product ion spectrum data is used in target compound identification without compromising sensitivity, accuracy and robustness in quantitative data reporting.




Fig. 10 MRM Product ion spectrum data for carbendazim in 3 calibration levels (0.010-0.200mg/kg) spiked into a food matrix was compared with an authentic library spectrum of carbendazim. In all library searches the similarity score was greater than 99 indicating a very high confidence in compound verification and reporting.

Data Reporting

Automated reference library matching and quantitation results can be simply viewed using LabSolutions Insight software (Fig 11).

LabSolutions Insight software helps to review by exception and to reduce false positive reporting by verifying compound identification using library matching scores and retention time variation from a calibration standard.




Fig. 11 LabSolutions Insight software helps to review quantitative and reference library matching results quickly and easily. Flexible filtering and sorting tools can be used to help reduce reporting false detects, especially in high throughput laboratories by filtering results based upon a similarity score with a reference library product ion spectrum.

Conclusions

False positive results are a major issue for all pesticide residue monitoring laboratories. EU regulations require that retention time and the ion ratio between 2 MRM transitions are within a set threshold. However, even applying this criteria false positives may occur for certain pesticide/commodity combinations.

In this application paper, we have applied MRM Spectrum Mode to identify and quantify 193 target pesticides in a number of different sample matrices. The library score is used as an additional identification criterion in order to improve identification confidence.

Acquisition of the MRM Spectrum mode method (1,291 MRM transitions) did not compromise data quality when compared to a conventional 2 MRM per compound method (386 MRM transitions) with consistent signal response and repeatability in both methods. The MRM product ion spectrums were demonstrated to be consistent across the linear range and between different matrices. The method acquired data in both positive and negative ion modes with a polarity switching time of 5 msec enabling fast cycle times and a high data collection rate.

All 1,291 MRM transitions were acquired throughout the MRM window. No 'triggering' of MRM transitions was necessary due to the short dwell times that were applied using the LCMS-8060. Therefore, MRM transitions can be swapped between qualifier and qualifier if needed and the peak shape of the additional MRM transitions can be assessed.

References

- 1 European Commission SANTE/11945/2015. Guidance document on analytical quality control and method validation procedures for pesticides residues analysis in food and feed.
- 2 Schürmann A., Dvorak V., Crüzer C., Butcher P., Kaufmann A., False-positive liquid chromatography/tandem mass spectrometric confirmation of sebutylazine residues using the identification points system according to EU directive 2002/657/EC due to a biogenic insecticide in tarragon. Rapid Communications Mass Spectrometry, Volume 23, Issue 8, April 2009, Pages 1196-1200.
- 3 Kaufmann A., Butcher P., Maden K., Widmer M., Giles K., Uriá D.. Are liquid chromatography/electrospray tandem quadrupole fragmentation ratios unequivocal confirmation criteria? Rapid Communications, Mass Spectrometry, Volume 23, Issue 7, April 2009, Pages 985-998.
- 4 Pozo Ó., Sancho J., Ibáñez M., Hernández F., Niessen W., Confirmation of organic micropollutants detected in environmental samples by liquid chromatography tandem mass spectrometry: Achievements and pitfalls, TrAC Trends in Analytical Chemistry, Volume 25, Issue 10, November 2006, Pages 1030-1042.

First Edition: Mar. 2017



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2017

Shimadzu Corporation

www.shimadzu.com/an/

For LabSolutions Version 5.82 and Later

LC/MS/MS Method Package for Residual Pesticides Ver. 2

- World's largest compound panel (646 pesticides in a single method)
- Ultra-high-speed method, with detection in 10.5 minutes
- Ultra-high-speed detection delivers robust and reproducible data quality for extended pesticide programs
- Pretreatment program to improve the shape of peaks for polar pesticides
- Excellent data stability in combination with the LCMS-8000 series




Pesticide Screening Analysis for 646 Pesticides in 10.5 Minutes

Using this method package in combination with the LCMS-8050/8060 enables performing a pesticide analysis of 646 pesticides in 10.5 minutes.

Since the LCMS-8050/8060 have a rapid polarity switching time of 5 msec, a single multi-residue LC/MS/MS method supports the analysis of 612 pesticides in positive ion mode and 34 compounds in negative ion mode (Three MRM transitions for over 99% targeted pesticides, resulting in 1,919 transitions in total).


MRM Chromatograms for 646 Pesticides Spiked into a Mint Extract at 0.01 mg/kg



Maintains High Data Quality Even with High-Speed Analysis


The LCMS-8050/8060 delivers exceptional sensitivity and enhanced capability in pesticide screening. By developing unique technologies for fast data acquisition, the LCMS-8050/8060 systems create new opportunities for high data density analysis and expanded pesticide screening programs. For example, the peak area variation for 25 pesticides eluting over a nine-second time window (6.45-6.60 minutes) spiked into a mint matrix extract at the reporting limit of 0.01mg/kg varied between 1.1-5.9% (%RSD; n=6).

MRM Chromatograms for 25 Components Detected with Retention Times of 6.45 to 6.60 Minutes
(Three MRM Events Per Compound, One Polarity Switching) and
Ratio of Dwell Time (green), Pause Time (red) and Polarity Switching Time (blue) in 460 msec of
Loop Time When 25 Pesticides are Simultaneously Monitored (Top right)




LC/MS/MS Method Package for Residual Pesticides Ver. 2

Co-Injection Improves the Peak Shapes for Polar Compounds




Samples prepared by the QuEChERS method are typically extracted with 100% acetonitrile. While this approach is effective for most pesticides, it results in poor peak shape and response for highly polar pesticides. To enhance the peak shape and detection of highly polar pesticides, the acetonitrile extract was co-injected with a larger volume of water. This technique achieves band compression and delivers a marked improvement in peak shape and response for early eluting compounds, particularly for highly polar pesticides. A pretreatment program for co-injection is incorporated in the method for 646 pesticides.

Note: Co-Injection mode is only supported by the SIL-30AC.




Maintains Stability with Different Matrices



To assess the robustness of the system, three matrices (apple, mint and tomato) were spiked with 646 pesticides at 0.05 mg/kg and repeatedly injected (each matrix sample was analyzed with 100 repeat injections; the total batch run was 300 injections). Thiofanox-sulfone, dipropetryn, and probenazole were selected as representative response markers with retention times close to four, six, and seven minutes, respectively. Within a matrix, variance was less than 5.9%RSD for all compounds.

Although the absolute peak area changes with different food matrices, the response between injection 1 and injection 100 for two pesticides (probenazole and dipropetryn) within a single matrix has a variance less than 5.7%RSD.



The new LC/MS/MS Method Package for Residual Pesticides Ver.2 extends the capability of the LCMS-8000 series in pesticide analysis and creates new opportunities for expanded pesticide screening programs supported by a pre-packaged single LC/MS/MS method for 646 pesticides.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.

LC/MS/MS Screening Software

LabSolutions Insight Library Screening



Software Helps with Screening Target Components

Screening Using the MRM-Automatic Product Ion Scan Method


By using a method with three product ion scan events specified at different collision energies (CE) for MRM events, three MS/MS spectra can be obtained while simultaneously detecting compounds with high sensitivity. Consequently, by searching the library using a merged spectrum of the three spectra, highly accurate identification results can be obtained.

Type	Event#	+/-	Compound Name	m/z	Time (1.623 min - 8.982 min)
MRM	21	+	Bromovalerylurea	223.00>17	
- Product Ion Scan	22	+	> CE-15.0, 20.00:1000.00		
(Product Ion Scan)	23	+	> CE-30.0, 20.00:1000.00		
(Product Ion Scan)	24	+	> CE-45.0, 20.00:1000.00		
MRM	25	+	Chlorpheniramine	275.00>23	
- Product Ion Scan	26	+	> CE-15.0, 20.00:1000.00		
(Product Ion Scan)	27	+	> CE-30.0, 20.00:1000.00		
(Product Ion Scan)	28	+	> CE-45.0, 20.00:1000.00		
MRM	29	+	7-aminonimetazepam	266.10	
- Product Ion Scan	30	+	> CE-15.0, 20.00:1000.00		

MRM	Acq. Time: 3.986 - 5.986 min	Compound Name: Bromovalerylurea
-----	------------------------------	---------------------------------

Ch	Precursor m/z	Product m/z	Pause Time	Dwell Time	Q1 Pre Bias	CE	Q3 Pre Bias	Use for Survey
Ch1	223.00	179.95	1.0	1.0	-11.0	-12.0	-20.0	<input checked="" type="checkbox"/>
Ch2	223.00	57.10	1.0	1.0	-11.0	-20.0	-23.0	<input checked="" type="checkbox"/>
Ch3								
Ch4								
Ch5								

The Insight Library Screening software can display library search results for multiple analytes at the same time. It can display the MRM chromatogram, MS/MS spectrum, library search results for compounds contained in the sample, and structural formulas of identified compounds.



Screening Results for Toxicological Substances Added to Whole Blood

The sedative hypnotic drug bromovalerylurea was identified from the MRM chromatogram and library search results (similarity score 87).

Screening Using Multi-MRM Method


MS/MS spectra can be obtained by specifying multiple channels for MRM events and plotting the intensity of each channel as MS/MS spectrum peaks. By searching the library using the MS/MS spectrum obtained, samples can be screened with high accuracy.

Type	Event#	+/-	Compound Name	m/z	Time (0.800 min - 13.247 min)
MRM	194	+	Decoquinat	418.25>372.20	
MRM	195	+	Virginiamycin S1	824.35>205	
MRM	196	+	Malachite Green	329.20>313.15	
MRM	197	+	Norgestimate	370.25>124.10	
MRM	198	+	Coumaphos	363.00>227.00	

MRM	Acq. Time:	11.562	-	12.062	min	Compound Name:	Malachite Green
-----	------------	--------	---	--------	-----	----------------	-----------------

Ch	Precursor m/z	Product m/z	Pause Time	Dwell Time	Q1 Pre Bias(V)	CE	Q3 Pre Bias(V)
Ch1	329.20	313.15	1.0	3.0	-13.0	-38.0	-22.0
Ch2	329.20	208.15	1.0	3.0	-12.0	-36.0	-21.0
Ch3	329.20	165.10	1.0	3.0	-16.0	-65.0	-29.0
Ch4	329.20	241.10	1.0	3.0	-10.0	-55.0	-25.0
Ch5	329.20	239.10	1.0	3.0	-16.0	-65.0	-23.0
Ch6	329.20	192.10	1.0	3.0	-16.0	-60.0	-19.0
Ch7	329.20	284.10	1.0	3.0	-10.0	-45.0	-13.0
Ch8	329.20	235.15	1.0	3.0	-16.0	-47.0	-24.0
Ch9	329.20	285.15	1.0	3.0	-12.0	-36.0	-13.0
Ch10	329.20	297.10	1.0	3.0	-12.0	-53.0	-30.0
Ch11	329.20	268.15	1.0	3.0	-10.0	-62.0	-27.0
Ch12							
Ch13							
Ch14							

The Insight Library Screening software can display library search results for multiple analytes at the same time. It can display the MRM chromatogram, multi-channel MRM-based MS/MS spectrum, library search results, structural formulas, and MRM-based quantitative calibration curves for compounds contained in the sample.




Results from Screening for Veterinary Drugs Added to Foods (Salmon Extract)


Malachite green disinfectant identified from the MRM chromatogram and library search results (similarity score 99).

Printing Sophisticated Reports

Using the Library Screening software, separate reports can be printed for each data file. Summary reports or detailed results for each compound can also be printed.



Detailed results



Summary report



For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.

Shimadzu Corporation
www.shimadzu.com/an/

Application News

No.L497

Supercritical Fluid Extraction / Chromatography

Using the Nexera UC Online SFE-SFC-MS System to Analyze Residual Pesticides in Agricultural Products

The Nexera UC online SFE-SFC-MS system combines supercritical fluid extraction (SFE) and supercritical fluid chromatography (SFC) in one online system, so that the entire process from extraction of target components to acquisition of data can be performed completely automatically. Furthermore, the system can add polar organic solvents (modifiers) to the supercritical carbon dioxide fluid during SFE and SFC, so that the system can be used to extract and analyze components with a wide range of polarities.

Meanwhile, ever since the positive list system was enacted in 2006 in Japan for residual pesticides in foods, which applies to more than 800 types of pesticides, there has been increasing demand for a system able to simultaneously analyze multiple pesticides with a wide range of properties, including pretreating samples.

This article describes an example of using the Nexera UC online SFE-SFC-MS system to analyze residual pesticides in agricultural products.

■ Online SFE-SFC-MS System

The operating principle of the Nexera UC online SFE-SFC-MS system is shown in Fig. 1. The extraction vessel filled with the sample is placed in the SFE unit and heated to an internal temperature of 40 °C (Fig. 1A). Then supercritical carbon dioxide fluid is pumped into the extraction vessel. After filling the vessel, the flow is stopped to allow static extraction of target components (Fig. 1B). After static extraction, the fluid is pumped through the extraction vessel for dynamic extraction (Fig. 1C). During dynamic extraction, extracted substances flow from the extraction vessel and into the analytical column. However, due to the high level of contaminant components in agricultural products, passing all the extract substances through the analytical column or mass spectrometer could damage the column or contaminate the mass spectrometer. Therefore, the Nexera UC online SFE-SFC-MS system splits the flow to send only a portion of the substances extracted from dynamic extraction through the analytical column. After dynamic extraction, fluid is only sent through the analytical flow line, where the analytical column is used for gradient separation and the mass spectrometer for detecting the target components (Fig. 1D).

■ Sample Preparation

The QuEChERS is a well-known method that prioritizes simplicity and speed and is commonly used to pretreat agricultural products for residual pesticide analysis. However, the method involves many steps, such as adding reagents, solvent extraction, purification by dispersive solid phase extraction, and centrifugal separation. In contrast, the online SFE-SFC-MS system requires only mixing 1 g of agricultural product crushed with a mixer with 1 g of a dehydrating agent* and placing the mixture in the extraction vessel, as shown in Fig. 2. Consequently, the system improves analytical productivity, reduces the environmental impact, and also avoids human errors involved in the pretreatment steps. Using a dedicated rack changer, the system can continuously extract and analyze up to 48 samples at a time.

* "Miyazaki Hydro-Protect" Patent No. 3645552




Fig. 1 Analysis Flow by Online SFE-SFC-MS





Fig. 2 Sample Preparation

Table 1 Analytical Conditions

[SFE]		[SFC]	
Solvent	: A) Super critical fluid of CO ₂ B) 0.1 % Ammonium formate in methanol	Column	: Shim-pack UC-RP (250 mm L. × 4.6 mm I.D., 5 μm)
Flowrate	: 5 mL/min	Mobile Phase	: A) Super critical fluid of CO ₂ B) 0.1 % Ammonium formate in methanol
Extraction	: 0-3 min. Static mode (B.Conc. 5 %) 3-6 min. Dynamic mode (B.Conc. 5 %)	Time Program	: B.Conc. 0 % (0 min.) → 10 % (11 min.) → 30 % (14 min.) → 40 % (14.01-17 min.)
Extraction		Flowrate	: 3 mL/min
Vessel Temp.	: 40 °C	Make-up	: 0.1 % Ammonium formate in methanol (0.1 mL/min.)
BPR Pressure	: A) 14.8 MPa, B) 15 MPa (split rate: 3 %)	Column Temp.	: 40 °C
Make-up	: 0.1 % Ammonium formate in methanol (0.4 mL/min.)	BPR Pressure	: A) 15 MPa, B) 40 MPa
		Detector	: LCMS-8050 MRM mode

■ Analysis of Standard Mixture of Pesticides


The standard mixture sample of 510 pesticide components were mixed with a dehydrating agent and analyzed using the analytical conditions indicated in Table 1. Fig. 3 shows the results. Using the system, we were able to accomplish the entire process, from extraction to data acquisition, in about 45 minutes per analysis. For 327 components, we obtained good repeatability for the concentration range from 1 to 100 ng/g (less than 30 %RSD for relative standard deviation for peak area at respective concentrations) and good linearity (contribution ratio of at least $R^2 = 0.99$). Table 2 also shows how pesticides with a wide range of polarities were analyzed with good repeatability and linearity.

**Fig. 3 Mass Chromatogram of Standard Pesticide Mixture Solution****Table 2 Repeatability and Linearity for Representative Pesticides**

Compounds	LogPow	Repeatability (%RSD, n=5)	Range (ng/g)	R ²
Ethofenprox	6.9	6.1	1-100	0.9991
Hexaflumuron	5.68	6.8	1-100	0.9992
Benzofenap	4.69	1.4	2-200	0.9990
Mepronil	3.66	4.6	1-100	0.9993
Prometryn	3.34	2.7	1-100	0.9994
Fenamidone	2.8	3.0	2-200	0.9991
Ethylchlozate	2.5	3.0	1-100	0.9996
Imazosulfuron	1.6	6.2	1-100	0.9998
Bensulfuron methyl	0.79	8.1	1-100	0.9996
Primisulfuron methyl	0.2	5.5	1-100	0.9994
Halosulfuron methyl	-0.02	5.5	1-100	0.9996
Azimsulfuron	-1.4	4.2	1-100	0.9998


<Acknowledgments>

This Application News bulletin includes results obtained in cooperation with Osaka University, Kobe University, and the Miyazaki Agricultural Research Institute from the program for the "Development of Systems and Technology for Advanced Measurement and Analysis," sponsored by the Japan Science and Technology Agency (JST). We are deeply grateful to all those involved.

**Fig. 4 MRM Chromatograms of Representative Pesticides**

■ Analysis of a Tomato

Analysis of 10 ng/g of 510 pesticide components added to a tomato resulted in good repeatability (less than 20 %RSD for the relative standard deviation of the peak area) and a good recovery rate (70 to 120 %) for 248 components. Plots of LogPow and recovery rate results are shown in Fig. 5. It shows that pesticides with a wide range of polarities were analyzed with good recovery.

**Fig. 5 LogPow vs. Recovery Rate for Tomato Analysis**

Application News

No. G294

Gas Chromatograph

Analysis of Organophosphorus Pesticides Using Nexis GC-2030

Cases have been reported of health problems due to foods contaminated with pesticides, and there is currently heightened interest in food safety countermeasures. Using a detector with high selectivity for specific components, or a mass spectrometer highly capable of qualitative analysis are effective when analyzing trace components in foods and other samples in which there are many impurities.

The FPD-2030 flame photometric detector, which is installed in Nexis GC-2030 gas chromatograph, has the world's highest level of sensitivity* thanks to the optimized nozzle shape and the advanced dual focus system.

In the analysis of pesticides in foods, this detector provides high sensitivity and high stability.

In this Application News, we introduce an analysis of organophosphorus pesticides using Nexis GC-2030 gas chromatograph, which is equipped with the FPD-2030.

*As of February 2017

E. Kobayashi, T. Murata

Analysis Results

A mixture standard solution of 54 organophosphorus pesticides* (20 mg/L) was introduced via split injection, and the elution positions of each pesticide were confirmed.

Table 1 Analytical Conditions

Model	: Nexis GC-2030
Detector	: FPD-2030 (P-mode)
Column	: SH-Rtx-1701 (0.25 mm I.D. × 30 m, d.f. = 0.25 μm)
Column Temperature	: 60 °C (2 min) - 25 °C/min - 150 °C (0 min) - 5 °C/min - 200 °C (12 min) - 5 °C/min - 280 °C (7 min) Total 50.6 min
Injection Mode	: Split 1 : 20
Carrier Gas Controller	: Constant Linear Velocity (He)
Linear Velocity	: 30 cm/sec
Injection Temperature	: 250 °C
Detector Temperature	: 275 °C
Injection Volume	: 1 μL





Fig. 1 Chromatogram of 20 mg/L Organophosphorus Pesticides

Trace Level Analysis

Table 2 and Fig. 2 show the analysis conditions and the chromatogram respectively for a trace level analysis of 5 µg/L organophosphorus pesticides via high-pressure splitless injection.

Table 2 Analysis Conditions for Low-Concentration Organophosphorus Pesticides

Model	: Nexis GC-2030
Detector	: FPD-2030 (P-mode)
Column	: SH-Rtx-1701 (0.25 mm I.D. × 30 m, d.f. = 0.25 µm)
Column Temperature	: 60 °C (1 min) - 20 °C/min - 180 °C (0 min) - 5 °C/min - 200 °C (10 min) - 7 °C/min - 280 °C (5 min) Total 37.4 min
Injection Mode	: High Pressure Splitless (300 kPa, 1 min)
Carrier Gas Controller	: Constan Linear Velocity (He)
Linear Velocity	: 46.8 cm/sec
Injection Temperature	: 260 °C
Detector Temperature	: 300 °C
Injection Volume	: 2 µL



	S/N		S/N
1: Ethoprophos	: 42	12: MEP (Fenitrothion)	: 17
2: Phorate	: 39	13: Isofenphos	: 15
3: Thiometon	: 42	14: PAP (Phenthoate)	: 16
4: Terbufos	: 30	15: Prothiofos	: 18
5: Etrimfos	: 33	16: DMTP (Mathidathion)	: 21
6: ECP (Dichlofenthion)	: 24	17: Butamifos	: 19
7: Dimethoate	: 26	18: Sulprofos	: 25
8: Tolclofos-methyl	: 23	19: Fensulfothion	: 25
9: Chlorpyrifos	: 16	20: EPN	: 25
10: Formothion	: 5	21: PMP (Phosmet)	: 22
11: MPP (Fenthion)	: 18	22: Pyraclofos	: 15

Fig. 2 Chromatogram of Low-Concentration (5 µg/L) Organophosphorus Pesticides

First Edition: Jul. 2017



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2017

Shimadzu Corporation

www.shimadzu.com/an/

Application News

No. C99

Liquid Chromatography Mass Spectrometry

Quantitative Analysis of Veterinary Drugs Using the Shimadzu LCMS-8050 Triple Quadrupole Mass Spectrometer

Foods in which chemical residues, like pesticides, feed additives, and veterinary drugs found in excess of maximum residue levels have been banned from sale in many countries around the world. Compounds that are subject to residue standards vary widely and the list is expected to grow. Because of this, there is a need for a

highly sensitive and rapid analytical technique to analyze as many of these compounds as possible in a single run. This Application News introduces an example of the high-sensitivity analysis of 89 veterinary drugs in a crude extract of livestock and fishery products.

Sample Preparation

The typical samples used in the analysis of veterinary drugs contain large amounts of lipids because they are commonly meat and fish samples. Sample preparation is extremely important to ensure excellent sensitivity and repeatability. To avoid the typical time-consuming and laborious solid phase extraction sample preparation procedure, the QuEChERS method, which is typically used for the preparation of vegetables, was selected to simplify sample preparation.

The QuEChERS method normally consists of two steps, the first is an acetonitrile extraction and the second a cleanup step, but this time only the acetonitrile extraction step was used.

* QuEChERS Extraction Salts kit: Restek Q-sep™ AOAC2007.01




Fig. 1 Sample Preparation Procedure

Improved Peak Shape Using Sample / Water Co-Injection

When conducting reversed phase chromatography, the peaks of polar compounds may split or collapse depending on the relationship between the sample solvent and mobile phase. In cases where the sample solvent is rich in organic solvent, the elution strength must be lowered (by substitution or dilution) with the addition of water. As the pretreated sample solvent in this analysis consists of 100 % acetonitrile, injection in that state into the LC/MS will result in split peaks for some of the substances (Fig. 2 left).

To eliminate as much of the time and effort typically associated with sample preparation, the pretreatment features of the autosampler (SIL-30A) were utilized to conduct co-injection of sample and water, which resulted in improved peak shapes.





Fig. 2 Comparison of Peak Shape

MRM Analysis of Matrix Standards

Fig. 3 shows the MRM chromatogram of the matrix standard solution consisting of the sample solution with added standard solution (data obtained using pork extract solution). Table 1 shows the lower limits of quantitation for the standard solution without added matrix and with added matrix, respectively. In a crude extract obtained by acetonitrile extraction alone, sensitivity was comparable to that obtained for most of

the compounds using only standard solution. Although there were several compounds for which the lower limit of quantitation was different in the standard solution than the matrix-added solution, rather than attributing this to matrix effects, it is thought to be caused by elevated background due to ions derived from contaminating components (Refer to Fig. 5).




Fig. 3 MRM Chromatograms of 89 Veterinary Drugs (10 µg/L pork extract solution with added standard solution)

Table 1 LOQs of Veterinary Drugs in Neat Standards and Matrix Standards and Calibration
Range of Veterinary Drugs in Matrix Standards

	Std. Solution		Matrix-Added Std. Solution	
	Min. Conc.		Min. Conc.	Max. Conc.
Gentamicin	0.5		1	50
Sulfanilamide	1		1	50
Levamisole	0.05		0.05	50
Lincomycin	0.01		0.01	10
5-Propylsulfonyl-1-benzimidazole-2-amine	0.05		0.05	10
Diaveridine	0.01		0.01	10
Trimethoprim	0.02		0.02	20
Marbofloxacin	0.01		0.01	50
Sulfisomidine	0.02		0.02	20
Norfloxacin	0.5		0.5	50
Ormetoprim	0.02		0.02	10
Thiabendazole	0.01		0.01	10
Ciprofloxacin	0.05		0.5	10
Neospiramycin I	0.01		0.05	10
Danofloxacin	0.1		0.1	10
Enrofloxacin	0.05		0.1	50
Oxytetracycline	0.01		0.1	50
Xylazine	0.01		0.01	10
Orbifloxacin	0.05		0.05	50
Sulfacetamide	1		1	50
Clenbuterol	0.01		0.01	10
Tetracycline	0.05		0.01	50
Spiramycin I	0.01		0.01	50
Sarafloxacin	0.5		0.5	50
Difloxacin	0.05		0.1	50
Sulfadiazine	0.02		0.1	20
Sulfathiazole	0.02		0.1	20
Sulfapyridine	0.02		0.1	20
Carbadox	0.05		0.05	10
Pyrimethamine	0.02		0.02	20
Sulfamerazine	0.02		0.02	20
Chlortetracycline	0.1		0.1	50
Tilmicosin	0.1		0.1	50
Thiamphenicol	1		1	50
Sulfadimidine	0.02		0.02	20
Sulfamethoxydiazine	0.01		0.02	10
Sulfamethoxypyridazine	0.02		0.02	20
Sulfisozole	0.01		0.01	50
Trichlorfon (DEP)	0.05		0.05	50
Sulfamonomethoxine	0.02		0.02	20
Furazolidone	1		1	50
Difurazone	0.05		0.05	50
Erythromycin A	0.01		0.01	50
Cefazolin	0.5		0.5	50
Sulfachloropyridazine	0.02		0.02	20
Sulfadimethoxine	0.02		0.02	10
Tylosin	0.05		0.05	50
Sulfamethoxazole	0.02		0.1	10
Sulfaethoxypyridazine	0.02		0.02	10
Tiamulin	0.01		0.01	50
Florfenicol	0.5		10	50
2-Acetylamino 5-nitrothiazole	0.05		0.05	50
Sulfatrazoxazole	0.01		0.01	5
Leucomycin	0.01		0.01	50
Sulfisoxazole	0.01		0.05	50
Oxolinic acid	0.01		0.1	50
Chloramphenicol	0.5		1	50
Clorsulon	0.5		1	50
Sulfabenzamide	0.01		0.01	10
Ethopabate	0.01		0.01	10
Sulfadoxine	0.02		0.02	20
Sulfaquinoxaline	0.02		0.02	10
Prednisolone	0.1		0.05	20
Ofloxacin	0.5		0.5	50
Flubendazole	0.01		0.01	50
Methylprednisolone	0.5		0.5	50
Nalidixic acid	0.01		0.01	50
Dexamethasone	0.5		0.5	50
Flumequine	0.01		0.01	50
Benzylpenicillin	0.5		0.5	50
Sulfantran	0.2		0.2	50
Sulfabromomethazine	0.01		0.01	50
beta-Trenbolone	0.02		0.1	50
Emamectin B1a	0.01		0.01	50
alpha-Trenbolone	0.02		0.1	50
Piromidic acid	0.01		0.05	50
Zeranol	1		0.1	50
Ketoprofen	0.01		0.05	50
Testosterone	0.01		0.05	10
Famphur	0.05		0.05	50
Fenobucarb (BPMC)	0.01		0.01	50
Clostebol	0.05		0.05	50
Dichlofenac	0.01		0.01	50
Melengestrol Acetate	0.05		0.05	50
Temephos (Abate)	0.01		0.5	50
Allethrin	0.1		1	50
Closantel	0.01		0.01	10
Monensin	0.01		0.01	10

(Unit: µg/L)




Fig. 4 MRM Chromatograms in the Vicinity of the LOQ and Calibration Curves of Typical Compounds

Recoveries of Veterinary Drugs in Crude Extracts from Livestock and Fishery Products (Matrix Effect Verification)

We examined whether or not the matrix affected measurement of actual samples. This time, four types of food product samples were used, including shrimp, chicken meat, pork, and salmon. Standard solution was added to the acetonitrile extraction solution of each of these to obtain a final concentration of 10 $\mu\text{g/L}$, after

which the rates of recovery were determined. The results indicated that 90 % of the compounds were recovered at rates of 70 to 120 % and measurement was accomplished without any adverse matrix effects even though the crude extract solution was subjected only to acetonitrile extraction.




Fig. 5 Recoveries of Veterinary Drugs in Each of the Matrices

■ **Acetonitrile Extraction Efficiency Using QuEChERS Method**

To check the efficiency of acetonitrile extraction by the QuEChERS method, standard solution was added at stage (2) of Fig. 1 to obtain a concentration of 10 µg/L, and the recoveries were determined. Good recoveries of approximately 80 % were obtained in cases both

with and without the addition of matrix. However, relatively poor recoveries were seen for highly polar compounds such as tetracycline and quinolone. For these compounds, it is necessary to examine the use of a separate extraction solvent and extraction reagent.

Table 2 Recoveries (Pre-Spike)

Recovery	Without Matrix	With Matrix (Pork)	Compounds with Poor Recovery
< 50 %	17 (19 %)	13 (15 %)	Tetracyclines Quinolones
50 % - 70 %	1 (1 %)	8 (9 %)	
> 70 %	71 (80 %)	68 (76 %)	

■ **Robustness**

We checked the long-term stability of the instrument using a solution of pork crude extract (spiked with 10 µg/L standard solution). Even after continuous

measurement of an extremely complex matrix over a period of 3 days, we were able to obtain stable data.




Fig. 6 Area Plot and %RSD of Typical Compounds with Continuous Analysis

Table 3 Analytical Conditions

Column	: Shim-pack XR-ODS II (75 mm × 2.0 mm I.D., 2.2 µm)
Mobile Phase A	: 0.1 % Formic Acid - Water
Mobile Phase B	: Acetonitrile
Time Program	: 1 %B (0 min) → 15 %B (1 min) → 40 %B (6 min) → 100 %B (10-13 min) → 1 %B (13.01-16 min)
Flowrate	: 0.2 mL/min.
Injection Volume	: 2 µL (2 µL sample solution + 10 µL water)
Oven Temperature	: 40 °C
Ionization Mode	: ESI (Positive / Negative)
Probe Voltage	: +2.0 kV / -1.0 kV
Neburizing Gas Flow	: 3.0 L/min.
Drying Gas Flow	: 10.0 L/min.
Heating Gas Flow	: 10.0 L/min.
Interface Temperature	: 400 °C
DL Temperature	: 200 °C
Block Heater Temperature	: 400 °C

Application News

No. C161

Liquid Chromatograph Mass Spectrometry

Multi-Residue Veterinary Drug Analysis of >200 Compounds using MRM Spectrum Mode by LC-MS/MS

Veterinary drugs are used for therapeutic, metaphylactic, prophylactic and growth promotion purposes. To provide an assurance that food from animals is safe with regards to residues of veterinary medicines, regulatory authorities have established Maximum Residue Limits (MRL's) for certain drugs in target tissues and animal species. Some pharmacologically active compounds identified by regulatory authorities have been prohibited and their hazardousness at all levels are being considered (EU regulation EC 37/2010; Commission Decision 2003/181/EC; 21CFR Part 556 Tolerances for Residues of New Animal Drugs in Food). In this article, we describe how a triple quadrupole mass spectrometer, which is both highly sensitive and selective, contributes to reducing false positive and false negative reporting when using a measurement mode called MRM Spectrum mode. MRM Spectrum mode acquires a high number of fragment ion transitions for each target compound and generates fragmentation spectra that can be used in routine library searching and compound verification using reference library match scores.

David Baker ^{*1}, Laetitia Fages ^{*2}, Eric Capodanno ^{*2}, Neil Loftus ^{*1}
^{*1} : Shimadzu, Manchester, UK
^{*2} : Phytocontrol, Nimes, France

■ Samples and Analysis Conditions

Samples of beef, egg, honey, milk and salmon were extracted and spiked with veterinary drugs in the calibration range of 0.001 to 0.1 mg/kg. Repeatability was assessed at low and high concentrations. Samples were measured using Shimadzu's Nexera X2 UHPLC and LCMS-8060 triple quadrupole mass spectrometer (Table 1 and 2). Over 200 veterinary drugs were targeted and over 2,000 MRM transitions in both ESI +/- were monitored during a gradient elution time of 12 minutes.

Table 1 UHPLC Conditions

Liquid chromatography			
UHPLC	Nexera LC system		
Analytical column	Restek Biphenyl (100 × 2.1, 2.7 μm)		
Column temperature	40 °C		
Flow rate	0.4 mL/minute		
Solvent A	0.1 % formic acid 0.5 mM ammonium formate solution		
Solvent B	0.1 % formic acid in methanol		
Binary Gradient	Time (mins)	%B	Time (mins) %B
	0.00	2	14.60 2
	12.50	100	17.50 Stop
	14.50	100	

Table 2 MS/MS Acquisition Parameters


Mass spectrometry	
Mass spectrometer	Shimadzu LCMS-8060
Pause time/dwell time	1 msec/3 msec
Polarity switching time	Pos/neg switching time set to 5 msec
Scope	218 drugs in positive ion mode (including internal standards) 11 drugs in negative ion mode Structure Analytics (in house development tool)
Source temperatures (interface; heat block; DL)	350 °C; 300 °C; 150 °C
Gas flows (nebulising; heating; drying)	3 L/min; 10 L/min; 10 L/min

■ Advantages of MRM Spectrum Mode


The measurement method can be easily set using the MRM optimization tool and measurement window (MRM Synchronization) settings of LabSolutions LCMS. The method achieves high data densities and a high data sampling rate across each elution peak. This approach generates a consistent loop time and sampling rate producing reliable quantitation and peak integration. It also provides great operator-friendliness in routine simultaneous analysis of veterinary drugs by enhancing flexibility in qualifier and quantifier ion selection. The number of fragment ion transitions generated from a single precursor ion is limited only by the chemical structure of the veterinary drug.

■ Results

MRM Spectrum mode was used to acquire a high number of fragment ion transitions for each veterinary drug target. For chlortetracycline, 11 precursor-fragment ion transitions were acquired using optimized collision energies (Fig. 1). Acquiring a high number of fragment ion transitions enables generation of fragmentation spectra which can be used in library searching and compound verification for each veterinary drug. (Chlortetracycline is a tetracycline class of antimicrobials. According to the Sixth ESVAC report published in 2016, of the overall sales of antimicrobials in the 29 EU countries in 2014, the largest amount, expressed as a proportion of mg/PCU, was accounted for by tetracyclines (33.4 %). This is followed by penicillins (25.5 %) and sulfonamides (11.0 %). Chlortetracycline was selected as a representative target).



Compound name Chlortetracycline
Accurate mass 479.1216 [M+H]⁺
Formula C₂₂H₂₃ClN₂O₈
CAS 57-62-5



MRM Spectrum Mode

11 MRM's acquired for chlortetracycline at 10pg/uL in egg.

1:479.10>444.00 (+)	CE: -23V	7:479.10>300.80(+)	CE: -45V
2:479.10>461.95 (+)	CE: -35V	8:479.10>287.90(+)	CE: -53V
3:479.10>154.00 (+)	CE: -34V	9:479.10>274.95(+)	CE: -44V
4:479.10>98.05(+)	CE: -45V	10:479.10>370.95(+)	CE: -31V
5:479.10>260.05(+)	CE: -60V	11:479.10>285.85(+)	CE: -56V
6:479.10>303.05(+)	CE: -37V		

MRM Spectrum mode

Higher specificity

Higher reporting confidence

Library searchable fragment data.

The number of precursor-fragment ion transitions monitored is limited only by the structural chemistry of the molecule. Typically more than 10 precursor-fragment ion transitions were monitored for each veterinary drug.

Fig. 1 Utilization of MRM Spectrum Mode (Chlortetracycline)

Fig. 2 shows the MRM reference spectrum for chlortetracycline with assigned fragment structures. The MRM Spectrum mode is a measurement mode which combines MRM with the generation of a product ion spectrum. The product ion spectrum can be used for compound identification by searching a library.

As the collision energy was optimized for each fragment ion to generate a product ion spectrum, the library spectrum is highly specific and selective.




Fig. 2 MRM Reference Spectrum with Assigned Fragment Structures (Chlortetracycline)

Library Identification using MRM Spectrum Mode





Fig. 3 Library Searchable MRM Spectra in Different Matrices Spiked at 10 pg/μL (Chlortetracycline)

Fig. 4 shows the MRM spectra and the n=10 measurement results of four compounds for salmon extract spiked with virginiamycin S1 at a concentration of 10 pg/μL. The library match score was above 99 in all injections (MRM spectra of injections 1, 5 and 10 are

indicated). Also, the %RSD for oxytetracycline, sulfadimethoxine, ormetoprim, and virginiamycin spiked into salmon extract (n=10; 10 pg/uL) acquired using a conventional 2-MRM method was compared with that of the MRM spectrum method.




Compound name	Oxytetracycline		Sulfadimethoxine		Ormetoprim		Virginiamycin	
Number of MRM's	2MRM's	8MRMs	2MRM's	11MRMs	2MRM's	11MRMs	2MRM's	11MRMs
Mean peak area								
Quantitation ion	1890170	1729171	7809989	7227748	8291171	8160952	2232967	1956045
%RSD	3.74	3.04	1.49	1.46	1.54	1.18	0.91	1.65

Fig. 4 MRM Spectra and n=10 Results of Salmon Extract Spiked with Virginiamycin S1 at 10 pg/μL


Quantitation Results using MRM Spectrum Mode

To assess the robustness of the MRM Spectrum mode, the same sample was repeatedly injected. The method used complies with the identification criteria set out in the EU guidelines SANTE/11945/2015 that require the retention time and the ion ratio from at least 2 MRM ion ratios to be within acceptable tolerance limits. The absolute response and signal variability were

compared to those of the MRM Spectrum mode (Fig. 4). Both methods resulted in a variance of less than 4 %RSD (n=10 for each method; 10 pg/uL spiked into salmon matrix). Fig. 5 indicates MRM spectra and the calibration curve obtained for sulfamerazine as an example of quantitation results.



Compound name Sulfamerazine
Precursor 265.10 [M+H]⁺
Formula C₁₁H₁₂N₄O₂S
CAS 127-79-7



MRM Spectrum Mode

11 MRM's acquired for sulfamerazine in beef.

1:265.10>156.05 (+)	CE: -17V
2:265.10>92.05 (+)	CE: -30V
3:265.10>172.05 (+)	CE: -15V
4:265.10>108.05 (+)	CE: -26V
5:265.10>65.00 (+)	CE: -11V
6:265.10>110.05 (+)	CE: -22V
7:265.10>80.05 (+)	CE: -15V
8:265.10>38.95 (+)	CE: -14V
9:265.10>199.10 (+)	CE: -21V
10:265.10>93.10 (+)	CE: -30V
11:265.10>66.05 (+)	CE: -24V

Fig. 5 MRM Spectra and Calibration Curve of Sulfamerazine (1 pg/uL to 100 pg/uL)

Conclusion

The level of confidence in compound identification and verification was increased by using a higher number of MRM transitions for each veterinary drug target and thereby reducing false negative and false positive reporting. Although the number of transitions for each target is dependent upon the chemical structure of the target, typically more than 10 transitions can be monitored for each compound. MRM Spectrum mode combines conventional quantitation with the

generation of a high quality product ion spectrum which can be used to achieve highly reliable compound identification and verification by library searching. In this research, use of the MRM Spectrum mode was examined by quantifying and identifying 212 veterinary drugs (the method included 2,009 MRM transitions). Limits of detection, linearity or repeatability were not compromised compared to a conventional 2-MRM method.

First Edition: Aug. 2017



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2017

Application News

No. C103

Liquid Chromatography Mass Spectrometry

Analysis of Nivalenol, Deoxynivalenol, 3-Acetyldeoxynivalenol and 15-Acetyldeoxynivalenol Using Triple Quadrupole LC/MS/MS (LCMS-8050)

Nivalenol and deoxynivalenol are mycotoxins which are produced by the fusarium fungi. A provisional reference value of 1.1 ppm was established in Japan for deoxynivalenol (Notification No. 0521001 issued by the Pharmaceutical and Food Safety Bureau, Japanese Ministry of Health, Labour and Welfare on May 21, 2002). The test methods specified for deoxynivalenol are HPLC for both qualitative and quantitative analysis, and LC/MS for verification testing (Notification No.

0717001 issued by the Dept. of Food Safety, Pharmaceutical and Food Safety Bureau, Japanese Ministry of Health, Labour and Welfare on July 17, 2003).

This paper describes an LC-MS/MS method for high-sensitivity simultaneous analysis of the four compounds, nivalenol, deoxynivalenol and the deoxynivalenol metabolites, 3-acetyl-deoxynivalenol and 15-acetyl-deoxynivalenol.

■ Analysis of a Standard Mixture

Fig. 1 shows the chromatograms obtained using a 2 μ L injection of the four-component standard mixture (each 10 ppb), and Table 1 shows repeatability of retention time and peak areas for the four substances, respectively, using six repeat measurements.

Nivalenols are detected using the heated electrospray ionization (hESI) method in negative mode. Although water and acetonitrile alone can be used as the LC eluent for LC/MS analysis, higher sensitivity was obtained for each compound by adding low-concentration ammonium acetate (in this case, 0.5 mmol/L) to eluent A. Fig. 1 shows the mass chromatograms for the highest sensitivity MRM transitions for each compound. The analytical conditions are shown in Table 2.

Next, six repeat analyses of a 10 ppb standard solution were conducted, corresponding to approximately 1/100 the concentration of the provisional reference value. The relative standard deviations (%RSD) for the measured retention times and peak areas are shown in Table 1. Good repeatability was obtained for both retention time and peak area.

Table 1 Repeatability (10 ppb, n=6)

	R.T. %RSD	Area %RSD
Nivalenol	0.04	2.57
Deoxynivalenol	0.04	6.52
15-Acetyldeoxynivalenol	0.06	4.09
3-Acetyldeoxynivalenol	0.05	2.58




Fig. 1 MRM Chromatograms of a Standard Mixture (10 ppb each)

■ Linearity of Calibration Curves

Fig. 2 shows the calibration curves generated using the analytical conditions of Table 2. Excellent linearity with a coefficient of determination greater than $R^2 = 0.999$

was obtained for calibration curves using a concentration range from 1 to 250 ppb for each component.




Fig. 2 Linearity of Calibration Curves: a) Nivalenol b) Deoxynivalenol c) 15-Acetyldeoxynivalenol d) 3-Acetyldeoxynivalenol

■ Analysis of Wheat

Fig. 3 describes the sample pretreatment procedure for wheat. The wheat extract solution was purified using either the MultiSep #227 multi-function column (Romer Labs) or the Autoprep MF-T column (Showa Denko K.K.). The chromatograms generated using the samples prepared using the MultiSep #227 (unspiked samples) and the standard-spiked samples, respectively, are shown in Fig. 4. The standard mixture was added to obtain a final concentration of 25 ppb for the four components (about 1/40 of the provisional reference

value), respectively. No large contaminant peaks were detected in the chromatograms of the pretreated samples. Furthermore, although deoxynivalenol was detected, it was at a level below that of the provisional reference value. The spike-and-recovery rates for the four components were excellent, from 101 to 107 %, without any particular matrix effects. Even in samples pretreated using Autoprep MF-T, comparable spike-and-recovery test results were obtained.




Fig. 3 Pretreatment




Fig. 4 Chromatograms of Wheat
(Dotted line: Unspiked Sample, Solid line: Spiked Sample, Spiked at 25 ppb each)
a) Nivalenol b) Deoxynivalenol
c) 15-Acetyldeoxynivalenol d) 3-Acetyldeoxynivalenol

Table 2 Analytical Conditions

Column	: Shim-pack XR-ODS III (150 mm L. × 2.0 mm I.D., 2.2 μm)		
Mobile Phases	: A 0.5 mmol/L Ammonium Acetate - Water : B Acetonitrile		
Time Program	: 5 % B (0 min) → 45 % B (5.0 min) → 95 % B (5.01-7.0 min) → 5 % B (7.01 min) → STOP (12 min)		
Flowrate	: 0.3 mL/min		
Column Temperature	: 40 °C		
Injection Volume	: 2 μL		
Probe Voltage	: -3.0 kV (ESI-negative mode)		
DL Temperature	: 100 °C		
Block Heater Temperature	: 200 °C		
Interface Temperature	: 200 °C		
Nebulizing Gas Flow	: 2 L/min		
Drying Gas Flow	: 10 L/min		
Heating Gas Flow	: 10 L/min		
MRM Transition	: Nivalenol	371.05 > 281.20	CE: 16.0 V
	: Deoxynivalenol	355.10 > 59.00	CE: 22.0 V
	: 15-Acetyldeoxynivalenol	397.10 > 59.00	CE: 22.0 V
	: 3-Acetyldeoxynivalenol	397.10 > 59.00	CE: 26.0 V

Application News

No. C138

Liquid Chromatography Mass Spectrometry

Multi-Residue Analysis of 18 Regulated Mycotoxins by LC/MS/MS

D. Baker¹, C. Titman¹, J. Horner², N. Loftus¹:

¹ Shimadzu UK, ² Scientific Analysis Laboratories

Mycotoxins are one of the most important contaminants in food and feed due to their widespread distribution in the environment and toxic effects on humans and animals.¹⁾ Structurally, mycotoxins are a very diverse group with a wide range of physicochemical properties and low molecular weights.²⁾ They are produced by fungi (mould) frequently found on agricultural produce, and are often not visible to the naked eye.³⁾ Some of the most commonly contaminated food stuffs include wheat, oats, rye, corn, barley, rice, nuts and milk.⁴⁾

Due to the risks posed by mycotoxins in food they are regulated globally, including, the EU, US, China, Singapore and Brazil.⁵⁾ In the EU, reporting limits are harmonised in Regulation (EC) No 1886/2006 (amended by (EC) No 1126/2007) and sampling and analysis in Regulation (EC) No 401/2006.

LC/MS/MS is the technique most commonly employed for mycotoxin quantitation in order to achieve the necessary low reporting limits in complex food and feed matrices.

Experimental

Solvent extracts were provided by Scientific Analysis Laboratories (SAL, UK) following validated extraction protocols. Samples were analysed using the Nexera UHPLC and the LCMS-8060 triple quadrupole detector (Table 1). Calibration was performed using ¹³C internal standards spiked during sample extraction. All MRM transitions and associated internal standards for each compound are listed in Table 2. All solvents used during analysis were LCMS quality from Sigma-Aldrich.

Due to the wide range of physical and chemical properties of mycotoxins, different LC/MS/MS methods are typically developed for small groups of compounds with similar properties.

In this application paper a single LC/MS/MS method has been developed for the determination of 18 mycotoxins in food safety. Limits of quantification were at or below the maximum levels set in the EC/1886/2006 document. The scope of the method included Aflatoxins (B1, B2, G1, G2), Fumonisin (B1, B2, B3), Ochratoxin A (OTA) and Trichothecenes (3-acetyldeoxynivalenol (3AcDON), 15-acetyldeoxynivalenol (15AcDON), Deoxynivalenol (DON), Diacetoxyscripanol (DAS), Fusarenon-X (FUS X), HT-2, Neosolaninol (NEO), Nivalenol (NIV), T2, Zearalenone (ZON)) with an analysis cycle time of 12.5 minutes.

Table 1 Analytical Conditions

UHPLC	: Nexera LC System
Mobile Phase	: A; Water with additives B; Methanol with additives
Column	: Reversed phase column (100 mm L. x 2.1 mm I.D.)
Column Temperature	: 40 °C
Flowrate	: 0.4 mL/minute
Gradient	: B. Conc 15 % (0 min) → 25 % (1 min) → 40 % (2 min) → 41 % (4.5 min) → 100 % (7.5 - 10.0 min) → 15 % (10.10 min) → Stop (12.5 min)
LC-MS/MS	: LCMS-8060
Dwell Time	: 10 to 40 msec.
Pause Time	: 1 msec.
Ionisation Mode	: ESI +/-
Polarity Switching	: 5 msec.





Fig. 1 MRM Chromatograms of 18 Mycotoxins

AFB1 (aflatoxin B1; 1 µg/kg), AFB2 (aflatoxin B2; 1 µg/kg), AFG1 (aflatoxin G1; 1 µg/kg), AFG2 (aflatoxin G2; 1 µg/kg), OTA (ochratoxin A; 4 µg/kg), FB1 (fumonisin B1; 100 µg/kg), FB2 (fumonisin B2; 100 µg/kg), FB3 (fumonisin B3; 100 µg/kg), 15-AcDON (15-acetyldeoxynivalenol; 100 µg/kg), 3-AcDON (3-acetyldeoxynivalenol; 100 µg/kg), DON (deoxynivalenol; 100 µg/kg), DAS (diacetoxyscripanol; 100 µg/kg), FUS-X (fusarenon-X; 100 µg/kg), HT-2 (100 µg/kg), T2 (100 µg/kg), NEO (neosolaninol; 100 µg/kg), NIV (nivalenol; 100 µg/kg), ZON (zearalenone; 100 µg/kg).
For clarity only 2 MRM transitions are displayed per compound and the following MRM chromatograms were changed; neosolaninol (x0.3), T2 (x0.3), aflatoxins (x3), fumonisins (x2).

Table 2 All MRM's Measured in the Mycotoxin Method and Corresponding Calibration Range and R² Result

	Compound name	Parent ion	Ret. Time (mins)	MRM 1	MRM 2	MRM 3	ISTD	Calibration range µg/kg	R ²
1	Aflatoxin B1	[M+H] ⁺	6.773	313 > 241	313 > 285	313 > 269	¹³ C Aflatoxin B1	0.1 - 10	0.9988
2	Aflatoxin B2	[M+H] ⁺	6.621	315 > 259	315 > 287	315 > 243	¹³ C Aflatoxin B2	0.1 - 10	0.9995
3	Aflatoxin G1	[M+H] ⁺	6.453	329 > 243	329 > 200		¹³ C Aflatoxin G1	0.1 - 10	0.9998
4	Aflatoxin G2	[M+H] ⁺	6.219	331 > 245	331 > 285		¹³ C Aflatoxin G2	0.1 - 10	0.9965
5	Ochratoxin A	[M+H] ⁺	7.509	404 > 239	404 > 221	404 > 358	¹³ C Ochratoxin A	0.4 - 40	0.9969
6	Fumonisin B1	[M+H] ⁺	6.811	722 > 352	722 > 334	722 > 704	¹³ C Aflatoxin B2	10 - 1000	0.9937
7	Fumonisin B2	[M+H] ⁺	7.260	706 > 318	706 > 354	706 > 688	¹³ C Aflatoxin B2	10 - 1000	0.9998
8	Fumonisin B3	[M+H] ⁺	7.073	706 > 318	706 > 354	706 > 688	¹³ C Aflatoxin B2	10 - 1000	0.9991
9	Deoxynivalenol	[M+H] ⁺	2.372	297 > 279	297 > 249		¹³ C Deoxynivalenol	10 - 1000	0.9992
10	Diacetoxyscirpenol	[M+NH ₄] ⁺	6.349	384 > 229	384 > 307	384 > 247	¹³ C T2 Toxin	10 - 1000	0.9994
11	T2	[M+NH ₄] ⁺	7.206	484 > 185	484 > 215	484 > 245	¹³ C T2 Toxin	10 - 1000	0.9989
12	HT-2	[M+Na] ⁺	6.822	447 > 345	447 > 285		¹³ C T2 Toxin	10 - 1000	1.0000
13	Nivalenol	[M-CH ₃ COO] ⁻	1.684	371 > 281	371 > 311		¹³ C HT-2	10 - 1000	0.9991
14	Neosolaniol	[M+NH ₄] ⁺	3.227	400 > 215	400 > 305	400 > 185	¹³ C Deoxynivalenol	10 - 1000	0.9995
15	Fusarenon X	[M+H] ⁺	2.986	355 > 247	355 > 277		¹³ C Deoxynivalenol	10 - 1000	0.9987
16	Zearalenone	[M-H] ⁻	7.711	317 > 175	317 > 131	317 > 273	¹³ C T2 Toxin	10 - 1000	0.9985
17	15-Acetyldeoxynivalenol	[M+H] ⁺	4.406	339 > 261	339 > 297		¹³ C Deoxynivalenol	10 - 1000	1.0000
18	3-Acetyldeoxynivalenol	[M+H] ⁺	4.618	339 > 261	339 > 297		¹³ C Deoxynivalenol	10 - 1000	0.9986
19	¹³ C HT-2	[M+NH ₄] ⁺	6.844	464 > 278					
20	¹³ C T2	[M+NH ₄] ⁺	7.228	508 > 322					
21	¹³ C Aflatoxin B1	[M+H] ⁺	6.754	330 > 301					
22	¹³ C Aflatoxin B2	[M+H] ⁺	6.614	332 > 303					
23	¹³ C Aflatoxin G1	[M+H] ⁺	6.435	346 > 212					
24	¹³ C Aflatoxin G2	[M+H] ⁺	6.219	348 > 259					
25	¹³ C Ochratoxin A	[M+H] ⁺	7.516	424 > 250					

**Fig. 2 Calibration Curves for Selected Compounds**

Calibration Curves for Aflatoxin (0.1 – 10 µg/kg), Deoxynivalenol (10 – 1000 µg/kg), and Zearalenone (10 – 1000 µg/kg).

Conclusions

In this study a single method has been developed for the analysis of 18 regulated mycotoxins with an injection to injection cycle time of 12.5 minutes. This method achieves the required EU reporting limits (between 0.1 -10 µg/kg) with linear regression

coefficients R² typically greater than 0.998 (Fig. 2 and Table 1). The LC mobile phase, column and gradient were all optimised and provided chromatographic resolution of 15-acetyldeoxynivalenol and 3-acetyldeoxynivalenol.

References

- 1) A. Rahmani, S. Jinap, and F. Soleimany. 2009. Qualitative and Quantitative Analysis of Mycotoxins. Comprehensive reviews in food science and food safety. 8: 202-251.
- 2) M. E. Zain. 2011. Impact of mycotoxins on humans and animals. Journal of Saudi Chemical Society. 15 (2): 129-144.
- 3) M. Sameni, A. Dübecke and J. F. Weber. 2014. Simultaneous Multi- Residue Determination of Mycotoxins in Foods Using LC-MS/MS. Journal of Environmental & Analytical Toxicology. 5 (2): 1000259
- 4) J. W. Bennett and M. Klich. 2003. Mycotoxins. Clinical Microbiology Reviews. 16 (3): 497-516
- 5) Mycotoxins regulations for Food. http://www.mycotoxins.info/myco_info/consum_regu.html Accessed 6th September 2016

First Edition: Oct. 2016

**For Research Use Only. Not for use in diagnostic procedure.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2016

Shimadzu Corporation

www.shimadzu.com/an/

Application News

No.L506

High Performance Liquid Chromatography

Assay of Aflatoxin M₁ in Milk Based on Notification Test Methodology, Using Prominence-i and the RF-20A_{XS} Fluorescence Detector

Aflatoxin M₁ (AFM₁) is a mycotoxin suspected of carcinogenicity in humans that is detected in the milk of mammals that eat food contaminated with aflatoxin B₁. The notification "Handling of Aflatoxin M₁ in Milk" issued on July 23, 2015 (Notification No. 0723-[1] of the Department of Food Safety, PFSB, MHLW)¹⁾ sets a regulatory level for AFM₁ in milk of 0.5 µg/kg, and came into force on January 23, 2016.

The assay methodology for AFM₁ in milk was included in "Test Methodology for Aflatoxin M₁ in Milk" (Notification No. 0723-[5] of the Department of Food Safety, PFSB, MHLW)²⁾, which was announced on the same day and describes two test methodologies.

(1) Test method consisting of quantitation by HPLC with attached fluorescence detector and confirmation by LC/MS or LC/MS/MS.

(2) Screening method using an assay kit.

We describe an analysis of commercially available milk that is compliant with test method (1). We analyzed for AFM₁ in bovine milk using the Prominence-i integrated HPLC and the RF-20A_{XS} fluorescence detector. Under these conditions we were able to measure AFM₁ at a concentration of 1/10th Japan's regulatory level for AFM₁ in milk.

■ Analysis of Standard Aflatoxin M₁ Solutions

Chromatograms obtained after analysis of standard AFM₁ solutions (0.1 µg/L, equivalent to 1/100th the regulatory concentration) are shown in Fig. 1, and the analytical conditions used are shown in Table 1. The relative standard deviation (%RSD) of peak areas after repeating analysis six times was 3.4 %. Fig. 2 shows the calibration curve for 0.1 to 20 µg/L. Good linearity was achieved with a contribution ratio R² of ≥ 0.9999 within the concentration range. These results show the RF-20A_{XS} fluorescence detector can be used to analyze trace quantities of AFM₁ with high sensitivity and high precision.

When the standard AFM₁ solution of 0.1 µg/L is processed according to the pretreatment procedure shown in Fig. 3, which follows the notification methodology, it produces a sample equivalent to 1/100th the regulatory level for AFM₁ in milk (0.005 µg/kg).




Fig. 1 Chromatograms for Standard AFM₁ Solution Equivalent to 1/100th the Regulatory Concentration (0.1 µg/L, Test Repeated 6 Times)




Fig. 2 AFM₁ Calibration Curve (0.1-20 µg/L)

Table 1 HPLC Analytical Conditions

System	:Prominence-i
Column	:Shim-pack VP-ODS (150 mm L. × 4.6 mm I.D., 5 µm)
Mobile Phase	:Water/Acetonitrile = 3/1 (v/v)
Flowrate	:1.0 mL/min
Column Temp.	:40 °C
Detection	:RF-20A _{XS} , Ex. at 365 nm, Em. at 435 nm
RF Cell.	:Conventional Cell
Cell Temp.	:25 °C
Injection Volume	:100 µL

■ Analysis of Aflatoxin M₁ in Milk

We analyzed commercially available milk and milk with added AFM₁. AFM₁ was added to make up a concentration of 0.05 µg/kg in milk (1/10th the regulatory level), and pretreatment was performed according to the notification methodology.²⁾ The pretreatment procedure is shown in Fig. 3. Refer to the notification methodology²⁾ for further details.

An AflaStarTM R* immunoaffinity column from Romer Labs was used to remove contaminant constituents. The chromatograms obtained after analysis of these samples are shown in Fig. 4. (A) is the chromatogram for milk with added AFM₁, and (B) is the chromatogram for milk with no added AFM₁.

The analytical conditions were the same as those used in Fig. 1, which are shown in Table 1.

* "AflaStar" is a registered trademark of Romer Labs.
The AflaStarTM R can be purchased from Shimadzu GLC Ltd.




Fig. 3 Milk Pretreatment Procedure

- *1 A pretreatment centrifugation step is needed for raw milk and other milks that are not homogenized. Refer to reference²⁾ for details.
- *2 Refer to the annotations in reference²⁾ for detailed information on use of the immunoaffinity column.
- *3 AFM₁ can adhere to the container during drying, so it is recommended that silane-treated containers be washed with 20 % to 30 % aqueous acetonitrile then dried before use.
- *4 AFM₁ can adhere to glass containers used to hold samples for HPLC even when these containers have been treated with silane, so it is recommended that plastic containers be used.

The percentage recovery calculated according to Eqn. 1 shown below was 98 %. We found that using the RF-20Axs fluorescence detector allows for analysis at concentrations 1/10th the regulatory level with high sensitivity and good precision.

A small peak was observed at the AFM₁ elution position when milk with no added AFM₁ was analyzed. Using LC/MS/MS to analyze the milk with no added AFM₁ suggested this peak was derived from AFM₁, and the concentration of the substance present was below 1/100th Japan's regulatory level.

$$\text{Recovery rate (\%)} = \frac{(\text{Peak area of milk with added standard AFM}_1) - (\text{Peak area of milk with no added standard AFM}_1)}{\text{Peak area of standard AFM}_1 \text{ sample}} \times 100$$

Eqn. 1 Percentage Recovery Equation




Fig. 4 HPLC Chromatograms for Commercially Available Milk (A) With added standard AFM₁, (B) with no added standard AFM₁

References

- 1) "Handling of Aflatoxin M₁ in Milk" (July 23, 2015, Notification No. 0723-[1] of the Department of Food Safety, PFSB, MHLW)
- 2) "Test Methodology for Aflatoxin M₁ in Milk" (Notification No. 0723-[5] of the Department of Food Safety, PFSB, MHLW)

First Edition: Jun. 2016



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

No.L512

High Performance Liquid Chromatography

Analysis of Mycotoxins in Grain Using Mycotoxin Screening System

Mycotoxins are chemical products produced by organisms in the fungus kingdom and are toxic to humans, animals, and crops. As an example, aflatoxins are a type of mycotoxin that are some of the most carcinogenic naturally occurring substances in the world. They are classified as Group 1 carcinogens (carcinogenic to humans) by the WHO International Agency for Research on Cancer (IARC), and subject to strict regulations in many countries and regions of the world.

This Application News describes the screening analysis for mycotoxins in grain products (soft wheat flour and rice flour) using the i-Series Solution Package mycotoxin screening system.

■ i-Series Solution Package Mycotoxin Screening System

The screening system comprises a compact and easy to use integrated i-Series HPLC system together with analysis methods including sample pretreatment methods. The system comes ready to use and capable of data acquisition and analysis, including columns and method files designed for mycotoxin analysis, an instruction manual with analysis methods, and report templates. For screening applications, the system can determine whether mycotoxin levels in food are in excess of reference levels.




Fig. 1 Mycotoxin Screening System

Currently, HPLC and LC/MS are the most common techniques used to identify aflatoxins in food. With HPLC, fluorescent derivatization is often performed to improve sensitivity, though disadvantages of derivatization procedures are the time required and their complexity. Meanwhile, though LC/MS is more selective in terms of sensitivity, major financial investment into system is required.

The i-Series Solution Package comes with a built-in PDA detector, and can be further enhanced with an RF-20Axs fluorescence detector that offers world-class sensitivity. The package can also detect aflatoxins directly without derivatization.

■ Analysis of a Standard Solution

Mycotoxin targets of the screening system are shown in Table 1, chemical structures of some of these targets are shown in Fig. 2, and analytical conditions are shown in Table 2. The package includes analysis method files that eliminate the need to configure analytical conditions. An RF-20Axs was used to perform analysis with on-time excitation wavelength/emission wavelength switching.

Table 1 Screening Target Compounds

	Mycotoxin	Abbreviation	Matrix
1	Aflatoxin M ₁	AFM ₁	Milk
2	Aflatoxin G ₂	AFG ₂	Grain
3	Aflatoxin G ₁	AFG ₁	
4	Aflatoxin B ₂	AFB ₂	
5	Aflatoxin B ₁	AFB ₁	
6	Zearalenone	ZON	
7	Ochratoxin A	OTA	
8	Nivalenol	NIV	
9	Deoxynivalenol	DON	Apple
10	Patulin	PAT	




Fig. 2 Target Mycotoxin Structures

Table 2 Analytical Conditions

System	: Nexera-i 3D, RF-20Axs
Column	: Shim-pack GIST C18 (75 mm L. x 3.0 mm I.D., 2 μm)
Mobile Phase	: A) 20 mmol/L (Sodium) phosphate buffer (pH 2.5) B) Acetonitrile C) Methanol (Gradient elution)
Flowrate	: 1.0 mL/min
Column Temp.	: 55 °C
Injection Vol.	: 10 μL
Detection (RF-20Axs)	: AFB ₁ , AFB ₂ , AFG ₁ , AFG ₂ , AFM ₁ : Ex 365 nm, Em 450 nm : OTA, ZON : Ex 320 nm, Em 465 nm
Detection (Nexera-i 3D)	: NIV, DON : 220 nm (ch 1) : PAT : 276 nm (ch 2)

Although regulatory limits for mycotoxin levels in food can vary by country and region, the screening system is compatible with the strictest regulatory limits that are found in the EU (excluding regulatory limits in baby food). Chromatograms of a standard mixture with mycotoxin levels equivalent to EU reference levels¹⁾ is shown in Fig. 3.




Fig. 3 Chromatograms of a Standard Mixture

■ Analysis of Mycotoxins in Grain

This section describes an analysis of milled grains after pretreatment. Fig. 4 shows an overview of the sample pretreatment method. Further details can be found in the mycotoxin screening system instruction manual. Chromatograms of pretreated samples of soft wheat flour and rice flour and of pretreated samples of soft wheat flour and rice flour spiked with a standard mixture of mycotoxins that are produced in grains (shown among the screening target compounds listed in Table 1) are shown in Fig. 5 and Fig. 6. Comparing the area of each peak in the standard mixture that contains mycotoxins at EU reference levels and each peak in the flour samples allows identification of whether the mycotoxins present in flour samples are in excess and violation of reference levels. The system makes this determination without the need for complex analysis of results by the user, allowing for easy screening of target compounds.




Fig. 4 Sample Pretreatment




Fig. 5 Chromatograms of Soft Wheat Flour




Fig. 6 Chromatograms of Rice Flour

Footnotes
1) Converted concentrations in the standard mixture were obtained according to a pretreatment method described in the i-Series Solution Package Mycotoxin Screening System instruction manual.
2) MycoSpin™ 400 and MultiSep® 227 are registered trademarks of Romer Labs.



Mycotoxin Screening System

i-Series Solution Package




EU Criteria Concentrations of 10 Mycotoxin Components Detected with High Sensitivity in Only 14 Minutes

Mycotoxin (mold toxin) is a generic term for metabolites produced by molds on food products that are hazardous to human and animal health. To improve food safety, food processors have been inspecting ingredients for the presence of such mycotoxins. In only 14 minutes, this mycotoxin screening system is able to detect the presence of 10 mycotoxin components with high sensitivity at concentration levels specified by EU standards, which are the strictest in the world. Furthermore, because the system does not involve a sample derivatization process, samples can be measured much more efficiently.

Three Key Features

- Detects mycotoxins with high sensitivity at criteria concentrations specified by EU standards, which are the strictest standards in the world.
- Rapid screening detects the 10 components in only 14 minutes.
- Screening results and reports are available immediately after each analysis is finished.



Criteria Values for Respective Standards*1

Mycotoxin	High-Risk Foods	EU	Codex	Japan
Aflatoxins B1, B2, G1, and G2	Grains (wheat, etc.)	Total 4 to 15 µg/kg AFB1 2 to 12 µg/kg	Total 10 to 15 µg/kg	Total 10 µg/kg
Aflatoxin M1	Milk	0.05 µg/kg	0.5 µg/kg	0.5 µg/kg
Ochratoxin A	Wheat, etc.	2 to 10 µg/kg	5 µg/kg	Not specified
Patulin	Apples	25 to 50 µg/kg	50 µg/kg	50 µg/kg
Deoxynivalenol	Wheat	500 to 1750 µg/kg	1000 µg/kg	1100 µg/kg (tentative criteria value)
Nivalenol	Wheat	Not specified		
Zearalenone	Grains	20 to 400 µg/kg (2 to 3 mg/kg in feed)	Not specified	Not specified (1 mg/kg in feed)

*1: Excluding foods intended for infants
For more details, see the precautions indicated below.

Ready to Use - Measure Samples Immediately After Switching the System ON

This system is designed for rapid screening for 10 mycotoxin components in grain products, such as wheat and rice flours, apples, and milk. Using the provided kit, which includes a column and CD-ROM, containing optimized pretreatment methods that minimize effects from contaminant components and analytical parameters for regulated components, sample measurements can be started immediately.

In addition to the above screening kit and Nexera-i system, standard samples for the regulated components, pretreatment cartridges, and mobile phases are also required.



Pretreatment

Analysis and data acquisition


Detect EU Criteria Concentrations of Mycotoxins with High Sensitivity

Aflatoxins Detectable Directly without Fluorescent Derivatization


By using the i-Series built-in PDA detector with an RF-20Ax fluorescence detector, which offers the highest sensitivity levels in the world, the system can detect aflatoxins at concentration levels specified as criteria values*² in EU directives, without using fluorescent derivatization. The instruction manual provided with the system includes pretreatment methods optimized for target samples, which are grains (soft wheat and rice flours), milk, and apples. The troubleshooting section includes key considerations for each process step, from extraction to acquisition, which helps ensure that reliable data can be acquired even when analyzing samples for the first time.

*2: Excluding foods intended for infants


Grain (Soft Wheat Flour)




Milk



Grain (Soft Wheat Flour)



Apples



All added mycotoxin concentrations are converted for foods.


XR Control Chart

Mycotoxin Screening System

i-Series Solution Package

Preventing the Loss of Data Integrity Due to Alteration or Replacement of Data

The LabSolutions report set function resolves important issues faced by analytical laboratories, such as preventing data alteration due to intentional or unintentional operating errors, and ensures data integrity. In addition to the extensive security functions provided in previous versions, LabSolutions report sets in DB and CS versions of LabSolutions achieve the visibility of software operations. The report set function also serves as a useful tool for addressing concerns that were difficult to resolve using previous reports that were printed out (such as alternation or replacement of data within reports).



Precautions

1. The mycotoxin screening kit provides unmodified information and other content obtained by Shimadzu for applications involving screening for mycotoxins in grains (soft wheat and rice flours), milk, and apples. Therefore, it is not recommended for other applications.
2. Customers should use results obtained using the kit based on their own judgment.
3. Indicated criteria values were obtained from the following laws and regulations.
 - Commission Regulation (EC) No 1881/2006 of 19 December 2006
 - Commission Regulation (EU) No 165/2010 of 26 February 2010
 - Codex Standard 193-1995
 - Distribution of the report of the ninth session of the codex committee on contaminants in foods (REP15/CF)
 - Handling foods that contain aflatoxins (Notification No. 0331-5 of Pharmaceutical and Food Safety Bureau, Ministry of Health, Labour and Welfare of Japan, March 31, 2011)
 - Method for testing total aflatoxins (Notification No. 0816-1 of Pharmaceutical and Food Safety Bureau, Ministry of Health, Labour and Welfare of Japan, August 16, 2011)
 - Handling milk that contains aflatoxin M1 (Notification No. 0723-1 of Pharmaceutical and Food Safety Bureau, Ministry of Health, Labour and Welfare of Japan, July 23, 2015)
 - Method for testing aflatoxin M1 in milk (Notification No. 0723-6 of Pharmaceutical and Food Safety Bureau, Ministry of Health, Labour and Welfare of Japan, July 23, 2015)
 - Setting tentative criteria values for deoxynivalenol in wheat (Notification No. 0521001 of Pharmaceutical and Food Safety Bureau, Ministry of Health, Labour and Welfare of Japan, May 21, 2002)
 - Partial revision of ministerial ordinance on milk and milk products concerning compositional standards, etc. and standards and regulations of foods, food additives, etc. (Notification No. 1126001 of Pharmaceutical and Food Safety Bureau, Ministry of Health, Labour and Welfare of Japan, November 26, 2003)



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.

© Shimadzu Corporation, 2016

First Edition: August 2016, Printed in Japan

Application News

No. C104

Liquid Chromatography Mass Spectrometry

Analysis of Diarrhetic Shellfish Toxin Using Triple Quadrupole LC/MS/MS (LCMS-8050)

The Japanese Ministry of Health, Labour and Welfare (JMHLW) specified in July, 1980 that the mouse bioassay (MBA) be used as the official method for diarrhetic shellfish toxin, and that the permissible exposure limit be 0.05 MU per gram of edible shellfish*. Shellfish in which the toxin exceeds this limit are prohibited from being sold at market according to the Japanese Food Sanitation Law Article 6, Item 2.

Due to significant technological advances since 1980, the sensitivity and accuracy obtained using the MBA method are significantly inferior compared to the high-precision, high-sensitivity possible using liquid chromatography mass spectrometry analytical instrumentation, which is currently used for this application. A complete transition to instrumental analysis for lipophilic marine biotoxins is scheduled to be implemented by January 2015 throughout the EU.

Based on this international trend, the JMHLW is currently considering migration to an instrumental analysis assay and setting new reference values to be used with instrumental analysis, in addition to the introduction of the Codex standard for okadaic acids (OA, Reference 1).

Table 1 CODEX Standard 292-2008

	Reference Value
OA Acids (OA and DTX group)	Permissible ingestion limit of 0.16 mg OA per kg of edible shellfish

Fig. 1 shows examples of LC/MS/MS high-sensitivity analysis of okadaic acid (OA), dinophysistoxin 1 (DTX1) and pectenotoxins (PTX1, 2, 6) and yessotoxin 1 (YTX1). Thus, it is possible to conduct high-sensitivity, high-separation analysis of each component.

Fig. 2 and Fig. 3 show MRM chromatograms of standard samples of OA and DTX1, respectively.

* The amount of toxin resulting in the death of two out of three mice following intraperitoneal administration of the equivalent of 20 g per edible shellfish.




Fig. 1 MRM Chromatograms of Diarrhetic Shellfish Toxin (1 ng/mL)




Fig. 2 MRM Chromatograms of Okadaic Acid (OA)




Fig. 3 MRM Chromatograms of Dinophysistoxin 1 (DTX1)

In addition, the calibration curves of OA and DTX1 are shown in Fig. 4. In both cases, the coefficient of determination R^2 was greater than 0.9999, indicating excellent linearity. Comparable linearity was also obtained for the other four substances. Thus, instrumental analysis of shellfish by LC/MS/MS offers high sensitivity and accuracy, making it a highly effective analytical method. For this reason it is attracting attention as an alternative to the traditional MBA method.




Fig. 4 Calibration Curves of OA and DTX1

Table 2 Analytical Conditions

Column	: InertSustain C8 (50 mm L. × 2.1 mm I.D., 3 μm)
Mobile Phases	: A 2 mmol/L Ammonium Formate – Water (pH adjusted to 8.5 with ammonia water) : B Methanol
Time Program	: 20 %B (0 min) – 100 %B (10 min) – 20 %B (10.01 min) – STOP (15 min)
Flowrate	: 0.2 mL/min
Column Temperature	: 40 °C
Injection Volume	: 10 μL
Probe Voltage	: +4.0 kV/-3.0 kV (ESI-positive / negative mode)
DL Temperature	: 200 °C
Block Heater Temperature	: 400 °C
Interface Temperature	: 350 °C
Nebulizing Gas Flow	: 3 L/min
Drying Gas Flow	: 10 L/min
Heating Gas Flow	: 10 L/min
MRM Transition	: (+) PTX6 906.50 > 835.40, PTX1 892.60 > 821.40, PTX2 876.50 > 805.40 : (-) OA 803.50 > 255.20, YTX1 1141.50 > 1061.30, DTX1 817.50 > 255.20

The diarrhetic shellfish toxin standards were provided courtesy of Dr. Toshiyuki Suzuki of the Japanese National Research Institute of Fisheries Science.

Reference 1: July, 2014, Food Safety Commission of Japan "Natural Poison Evaluation Report – Okadaic Acid Group Among Bivalves"
<http://www.fsc.go.jp/fscis/evaluationDocument/list?itemCategory=009>



Application News

No.L509

High Performance Liquid Chromatography

Analysis of Residual Antimicrobials in Meat with Antimicrobial Screening System (Part 1)

In May 2006, the positive list system took effect in Japan that, in principle, prohibited the sale of food products with residual levels of pesticides, animal feed additives, and veterinary drugs (collectively referred to as agricultural chemicals, etc.) above the level determined by the Minister of Health, Labour and Welfare.¹⁾

Antimicrobials are a type of veterinary drug and animal feed additive, and used for the treatment and prevention of disease in livestock and marine products. Quinolones and sulfonamides are two common groups of synthetic antimicrobials.

Shimadzu's quick and simple antimicrobial screening system is capable of screening 24 antimicrobials compounds. An example screening analysis targeting 12 widely used quinolones (old quinolones, new quinolones) is described here. Application News No.L510 also describes an example screening analysis targeting 12 antimicrobials including sulfonamides (also including antifolates).

Antimicrobial Screening System

Shimadzu's antimicrobial screening system is able to determine whether levels of antimicrobials subject to regulation in various countries are above a maximum residue limit (MRL). Table 1 shows MRLs for the target quinolones.

The system uses an i-Series integrated HPLC instrument and RF-20Axs high-sensitivity fluorescence detector, and comes with a sample pretreatment method, analytical column, analytical method files, and a UV spectral library that allow for immediate operation after installation. When the analysis method capable of simultaneous component analysis is used, the system can be used for simultaneous screening of multiple components. The determination of whether MRL have been exceeded can be viewed immediately after the system completes analysis. The photodiode array (PDA) detector built into the i-Series instrument supports highly accurate screening with compound identification based on retention times as well as UV spectra.

Table 1 Maximum Residue Limits and Sample Solution Concentration of Screening Target Compounds

Compound	MRL (mg/kg)	Sample Solution Concentration (mg/L)
1 Marbofloxacin	0.01	0.025
2 Ofloxacin	0.01	0.025
3 Ciprofloxacin	0.01	0.025
4 Danofloxacin	0.01	0.025
5 Enrofloxacin	0.01	0.025
6 Orbifloxacin	0.01	0.025
7 Sarafloxacin	0.01	0.025
8 Difloxacin	0.01	0.025
9 Oxolinic acid	0.01	0.025
10 Nalidixic acid	0.01	0.025
11 Flumequine	0.01	0.025
12 Piromidic acid	0.01	0.025

Sample Pretreatment

Sample pretreatment was performed based on Simultaneous Analysis Method I for Veterinary Drugs by HPLC (Livestock and Marine Products).^{2),3)} After acetonitrile extraction and removing fat by acetonitrile/hexane partitioning, sample solution was prepared by evaporation then redissolution. Fig. 1 shows the sample pretreatment protocol, and Table 1 shows sample solution concentrations after pretreatment. Refer to the instruction manual of the system for the details of the sample pretreatment procedure.




Fig. 1 Sample Pretreatment Protocol

Table 2 Analytical Conditions

System	: LC-2040C 3D, RF-20Axs
Column	: Shim-pack FC-ODS (150 mm L. × 4.6 mm I.D., 3 μm)
Mobile Phase	: A) 20 mM (Sodium) Phosphate Buffer Containing 0.1 M Sodium Perchlorate B) Acetonitrile/Methanol=90/10
Time Program	: Gradient Elution
Flowrate	: 1.0 mL/min
Column Temp.	: 40 °C
Injection Volume	: 5 μL
Detection	: <LC-2040C 3D> 280 nm <RF-20Axs> Ex at 290 nm, Em at 495 nm Ex at 325 nm, Em at 365 nm
Cell Temp.	: 40 °C (PDA), 30 °C (RF)

■ Analysis of Quinolones in Meat

Chicken and pork were used as samples. Chromatograms of the pretreated matrix solutions (blue line), matrix solutions spiked with standard solution to create matrix standard solutions (red line), and neat standard solution (black line) are shown in Fig. 2. Standard solution was added to matrix solutions to create matrix standard solutions with quinolone concentrations of 0.01 mg/kg. The analytical conditions are shown in Table 2. Analysis was performed with the fluorescence detector in dual wavelength mode. New quinolones (compounds 1 to 8 in Table 1) were detected at an excitation wavelength of 290 nm and fluorescence wavelength of 495 nm, and old quinolones (compounds 9 to 11 in Table 1) were detected at an excitation wavelength of 325 nm and fluorescence wavelength of 365 nm. Piromidic acid (compound 12 in Table 1) differs from other quinolones in exhibiting no fluorescence characteristics, and was detected using the PDA detector. Employing the analytical conditions shown, all 12 compounds were separated and eluted in approximately 22 minutes.

■ Similarity Calculation Using UV Spectral Library

The PDA-detected compound (piromidic acid) can be analyzed qualitatively based on UV spectra as well as retention times. Its spectrum can be checked for similarity against the library spectra. Fig. 3 shows a UV spectrum of piromidic acid in pork matrix spiked with a standard solution of piromidic acid at threshold concentration. The degree of similarity with the library spectrum was 0.998.




Fig. 3 Spectra of Piromidic Acid




Fig. 2 Chromatograms of Chicken and Pork:
Matrix Standard Solution (Red Line), Matrix Solution (Blue Line), Neat Standard Solution (Black Line)

<References>

- 1) The Japanese Positive List System for Agricultural Chemical Residues in Foods, Japan's Ministry of Health, Labour and Welfare
- 2) Multiresidue Method I for Veterinary Drugs, Etc. by HPLC (Animal and Fishery products)
Director Notice about Analytical Methods for Residual Compositional Substances of Agricultural Chemicals, Feed Additives, and Veterinary Drugs in Food (Syoku-An No. 0124001, January 24, 2005. Final amendments were made on May 26, 2006.), Japan's Ministry of Health, Labour and Welfare
- 3) "Standard methods of analysis in food safety regulation (for veterinary drugs and animal feed additives)" p.26-43, Japan Food Hygiene Association (2003), edited under the supervision of the Japan's Ministry of Health, Labour and Welfare

First Edition: Oct. 2016



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2016

Shimadzu Corporation

www.shimadzu.com/an/

Application News

No.L510

High Performance Liquid Chromatography

Analysis of Residual Antimicrobials in Meat with Antimicrobial Screening System (Part 2)

Antimicrobials are a type of veterinary drug and animal feed additive, and are used for the treatment and prevention of disease in livestock and marine products. Residual antimicrobials are often found in livestock and marine products, so threshold levels for antimicrobials are set by regulation to ensure the safety of the consumer based on amounts that do not harm human health.

Due to ongoing reports of recent cases of regulatory violations in various countries and the large number of compounds targeted for testing, there is a demand for quick and simple antimicrobial screening.

While Application News No.L509 described an example of using the antimicrobial screening system for screening 12 quinolone compounds, this Application News describes an example screening analysis of 12 antimicrobial target compounds including sulfanomides.

■ Sample Pretreatment

Sample pretreatment for analysis of residual antimicrobials in meat usually employs liquid-liquid extraction (and sometimes solid phase extraction), but this process takes time and effort. In this article, we employed a QuEChERS method designed to be more efficient and reduce pretreatment times. The QuEChERS method is used to pretreat vegetables and fruits for residual pesticide analysis.

After using the QuEChERS method to perform extraction and fat removal, sample solutions were prepared by evaporation and redissolution steps. Table 1 shows the maximum residue limits (MRLs) of target compounds and sample solution concentrations after sample pretreatment, and Fig. 1 shows the sample pretreatment protocol. Refer to the instruction manual of the system for the details of the sample pretreatment procedure.

Table 1 Maximum Residue Limits and Sample Solution Concentration of Screening Target Compounds

Compound	MRL (mg/kg)	Sample Solution Concentration (mg/L)
1 Sulfadiazine	0.01	0.025
2 Sulfamerazine	0.01	0.025
3 Sulfadimidine	0.01	0.025
4 Sulfamonomethoxine	0.01	0.025
5 Trimethoprim	0.01	0.025
6 Sulfamethoxazole	0.01	0.025
7 Ormetoprim	0.01	0.025
8 Sulfadimethoxine	0.01	0.025
9 Sulfaquinoxaline	0.01	0.025
10 Pyrimethamine	0.01	0.025
11 Difurazon	0.01	0.025
12 Nicarbazin*	0.01	0.025

*1: Concentration of N, N'-Bis(4-nitrophenyl)urea, the main constituent of nicarbazin.




Fig. 1 Sample Pretreatment Protocol

■ Analysis of Antimicrobials Including Sulfonamides in Meat

Chicken and beef were used as samples. The analytical conditions are shown in Table 2. Chromatograms of the pretreated matrix solutions (blue line), matrix solutions spiked with standard solution to create matrix standard solutions (red line), and neat standard solution (black line) are shown in Fig. 2.

Standard solution was added to matrix solutions to make up antimicrobial concentrations, including sulfonamide concentrations, of 0.01 mg/kg in matrix standard solutions. Standard solutions were prepared to the sample solution concentrations listed in Table 1.

The photodiode array (PDA) detector (six-wavelength) built in the i-Series instrument was used for detecting all target compounds. Employing the analytical conditions shown, all 12 compounds were separated and eluted in approximately 25 minutes.

Table 2 Analytical Conditions

System	: LC-2040C 3D
Column	: Shim-pack FC-ODS (150 mm L. × 4.6 mm I.D., 3 μm)
Mobile Phase	: A) 20 mM (Sodium) Phosphate Buffer Containing 0.1 M Sodium Perchlorate B) Acetonitrile/Methanol=80/20
Time Program	: Gradient Elution
Flowrate	: 1.0 mL/min
Column Temp.	: 50 °C
Injection Volume	: 20 μL
Detection	: 240 nm 270 nm 280 nm 285 nm 350 nm 380 nm
Cell Temp.	: 40 °C

■ Similarity Calculation Using UV Spectral Library

All target compounds in this Application News can be analyzed qualitatively based on UV spectra as well as retention times. Sample spectra can be also checked for similarity against library spectra. Fig. 3 shows a UV spectrum of sulfaquinoxaline in a beef matrix spiked with a standard solution of sulfaquinoxaline at threshold concentration. Degree of similarity with the library spectrum was 0.997.




Fig. 3 Spectra of Sulfaquinoxaline




Fig. 2 Chromatograms of Chicken and Beef:
Matrix Standard Solution (Red Line), Matrix Solution (Blue Line), Neat Standard Solution (Black Line)

First Edition: Oct. 2016



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2016

Antimicrobial Screening System

i-Series Solution Package



Screening of 24 Synthetic Antimicrobial Compounds That Remain in Meat

Synthetic antimicrobial agents are used as animal drugs and feed additives and they are known to remain in the body of farm animals. Residue limits prescribing the safe ingestible amounts of such agents have been set so that people can safely eat meat and processed meat products. Thanks to a pretreatment method*1 that minimizes the influence of impurities in meat (beef, pork, and chicken muscle), LC system enables screening of 24 synthetic antimicrobial compounds, which are regulated in Japan, Europe, and other regions, remaining in meat.

*1: Consumables for pretreatment are not included in the package. They need to be provided by the customer.

Three Features

- Optimized pretreatment method and two detectors support screening of 24 synthetic antimicrobial compounds.
- Detects synthetic antimicrobial agents of the standard residual concentration with high sensitivity.
- Displays screening results immediately after analysis is finished.

Target Compounds of the Antimicrobial Screening System

#	Target Compound
1	Marbofloxacin
2	Ofloxacin
3	Ciprofloxacin
4	Danofloxacin
5	Enrofloxacin
6	Orbifloxacin
7	Sarafloxacin
8	Difloxacin

#	Target Compound
9	Oxolinic acid
10	Nalidixic acid
11	Flumequine
12	Piromidic Acid
13	Sulfadiazine
14	Sulfamerazine
15	Sulfadimidine
16	Sulfamonomethoxine

#	Target Compound
17	Trimethoprim
18	Sulfamethoxazole
19	Ormethoprim
20	Sulfadimethoxine
21	Sulfaquinoxaline
22	Pyrimethamine
23	Difurazon
24	Nicarbazin

Antimicrobial Screening System


i-Series Solution Package

Detects Synthetic Antimicrobial Agents of the Standard Residual Concentration with High Sensitivity


Comprehensively Supports Processes from Pretreatment to Checking Results

This screening system consists of an i-Series LC system and a special kit that contains analytical methods and analytical columns in order to provide ready-to-use screening of synthetic antimicrobial agents. It can reduce the time it takes for users to learn pretreatment procedures and examine conditions. Compound names, retention times, UV spectra and other data of the regulated compounds are preregistered in the analytical methods for higher data reliability.

Quinolone Agents



Sulfa Drugs



(A standard solution was added to chicken meat by an amount equivalent to the standard concentration.)

Displays Results Immediately After Analysis Is Finished

Measurement results for each sample can be automatically created as quantitative reports and summary reports and output in PDF. Use of the multi-data report creation function makes it easier to view large amounts of data.

Precautions

1. The synthetic antimicrobial agents screening kit includes information and data that have been obtained by Shimadzu Corporation for the purpose of screening for synthetic antimicrobial agents remaining in meat. Avoid using the kit for purposes other than its intended purpose.
2. Customers are responsible for how any results from the use of this kit are used.
3. Standard residual concentrations are determined according to the laws and regulations below.
 - Ministerial Notification No. 297 of the Japan's Ministry of Health, Labour and Welfare
 - Development of applicable laws and regulations for the enforcement of Item 3, Article 11 of Japan's Food Sanitation Act revised based on the partial revision of the Food Sanitation Act Issued by Dept. of Food Safety, Pharmaceutical and Food Safety Bureau, Ministry of Health, Labour and Welfare (No. 1129001) on November 29, 2005 (final revision: May 31, 2007)

Data File	Judgement	Marbofloxacin	Ofloxacin
Beef_001.lcd	Pass	○ N.D.	○ N.D.
Beef_002.lcd	Fail	○ N.D.	× 10.70



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.

© Shimadzu Corporation, 2016
First Edition: August 2016, Printed in Japan

Application News

No.C105

Liquid Chromatography Mass Spectrometry

Analysis of Sulfamic Acid in Fertilizers Using LC/MS (LCMS-2020)

Sulfamic acid, due to its plant growth inhibiting effects, is subject to maximum limits in fertilizers as specified in the official standard¹⁾ for ordinary fertilizers according to the Japanese Fertilizers Regulation Act. According to the Testing Methods for Fertilizers²⁾ supervised by Japan's Food and Agricultural Materials Inspection Center (FAMIC), the ion chromatography (IC) method is specified as the test method for sulfamic acid in ammonium sulfate. It has been reported, however, that when applying this IC method with byproduct compound fertilizer (fertilizer produced by concentrating and drying liquid byproducts obtained from fermentation plants involved in amino acid production, etc.) samples that contain large amounts of organic matter, it is difficult to separate the sulfamic acid peaks from contaminant peaks generated from sample matrix.³⁾

In this application, we investigated the analytical conditions for LC/MS that would permit acquisition of mass information and provide high selectivity in order to eliminate the effects of contaminating components. The LCMS-2020 single quadrupole mass spectrometer was used for the analysis.

Good quantitative results were obtained, confirming the applicability of this method using byproduct compound fertilizer as the actual sample.

■ Analysis of Standard Solution

Table 1 shows the analytical conditions, and Fig. 1 shows chromatogram obtained using a standard solution (0.1 mg/L aqueous solution) of sulfamic acid.

As retention of a zwitterionic compound such as sulfamic acid is difficult using reversed phase conditions, we adopted conditions using a HILIC column. Isocratic analysis was conducted using a mobile phase consisting of acetonitrile / ammonium formate + formic acid (pH 3.2).

Applying the LC/MS method (ESI-Negative), we conducted selected ion monitoring (SIM) analysis using the deprotonated molecule at m/z 95.9. Fig. 2 shows the calibration curve. Excellent linearity was obtained over the entire concentration range of 0.001 to 0.1 mg/L, with a correlation coefficient greater than 0.999.




Fig. 1 Mass Chromatogram (SIM) of Sulfamic Acid (0.1 mg/L)




Fig. 2 Calibration Curve (0.001 – 0.1 mg/L)

Table 1 Analytical Conditions

Column	: Phenomenex Luna HILIC 20A (100 mm L. × 2.0 mm I.D., 5 μm)
Mobile Phases	: Acetonitrile/100 mmol/L Ammonium Formate+ Formic Acid (pH 3.2) = 90:10, v/v
Flowrate	: 0.2 mL/min
Column Temperature	: 40 °C
Injection Volume	: 1 μL
Probe Voltage	: -3.5 kV (ESI-negative mode)
DL Temperature	: 250 °C
Block Heater Temperature	: 400 °C
Nebulizing Gas Flow	: 1.5 L/min
Drying Gas Flow	: 15 L/min
Monitoring Ion (SIM)	: m/z 95.9

■ Analysis of Sulfamic Acid in Fertilizers

We verified the applicability of the LC/MS method using a byproduct compound fertilizer as an actual sample. The permissible content level of sulfamic acid is set based on the total amount of the principal component in each type of fertilizer. Here, taking the lower limit of quantitation of sulfamic acid in fertilizer as 1/5 the value of the minimum concentration permissible (sulfamic acid concentration 0.005 % per principal component 1 %), we conducted spike and recovery testing using a spike quantity equivalent to the lower limit of quantitation.

Fig. 3 shows the sample pretreatment procedure. The extraction method conforms to the Testing Methods for Fertilizers (2013) supervised by FAMIC. After weighing out 1 g of byproduct compound fertilizer, extraction was conducted using 100 mL of water, and after further diluting this 100 to 1 with water, the mixture was filtered to complete preparation of the fertilizer measurement solution.

As the total quantity of the principal component represented 5 % of the fertilizer content, the concentration of sulfamic acid corresponding to the lower limit of quantitation is calculated as 50 mg/kg of fertilizer. In the spike and recovery test, 0.5 mL of 100 mg/L standard sample was added to the fertilizer, and after letting the mixture stand for 30 minutes, a measurement solution was prepared using the same procedure. The concentration of sulfamic acid in the measurement solution is therefore 0.005 mg/L.

Representative chromatograms are shown in Fig. 4 including chromatograms of the standard sample (0.005 mg/L), the sample spiked with sulfamic acid, and the byproduct compound fertilizer measurement solution. Table 2 shows the analytical results. Sulfamic acid was not detected in the byproduct compound fertilizer, nor were there any noticeable peaks associated with contaminant components.

In the spike and recovery test, excellent results were obtained in continuous analysis (n=5), with an average recovery rate of 101 %. The LC/MS method investigated here in the analysis of highly contaminated byproduct compound fertilizer was demonstrated to permit quantitation by simply adding a dilution step following extraction, as opposed to the IC method which requires tedious processing to address the issue of high-contaminant content.




Fig. 3 Preparation Flow




Fig. 4 SIM Chromatograms of STD and Fertilizer Sample

Table 2 Repeatability of Peak Area and Retention Time in Spike and Recovery Test

	R.t (min)	Peak Area	Recovery (%)
1st	4.217	1564	103
2nd	4.252	1561	102
3rd	4.229	1508	99
4th	4.224	1511	99
5th	4.219	1534	100
Ave	4.228	1535	101
%RSD	0.336	1.735	

[References]

- 1) Notification Regarding Determination of the Official Standard for Ordinary Fertilizer Based on the Fertilizers Regulation Act, February 22, 1986, the Japan's Ministry of Agriculture, Forestry and Fisheries Notification No. 284, Final Revision
December 5, 2013 the Ministry of Agriculture, Forestry and Fisheries Notification No. 2939 (2013)
- 2) Food and Agricultural Materials Inspection Center (FAMIC): Testing Methods for Fertilizers
< <http://www.famic.go.jp/ffis/fert/sub9.html> >
- 3) Hiroi T., Shirai Y.: Simultaneous Determination of Sulfamic Acid and Ammonium Thiocyanate in Ammonium Sulfate by Nonsuppressed Ion Chromatography, Research Report of Fertilizer, 5, 1 – 12 (2012)

Application News

No.L513

High Performance Liquid Chromatography

Analysis of Nitrous Acid and Ammonium Thiocyanate in Fertilizers

After it is spread onto agricultural land, nitrogen fertilizer converts to nitrate-nitrogen or nitrite-nitrogen and leaches into subsoil and river water polluting the groundwater. Drinking water with a high nitrate-nitrogen or nitrite-nitrogen content is a potential public health risk. It causes methemoglobinemia that inhibits the oxygen-carrying capacity of blood and has caused deaths in infants outside Japan.

Fertilizers with high concentrations of nitrous acid and ammonium thiocyanate have a negative effect on plant growth, therefore maximum content levels (permitted content levels) for toxic substances are prescribed in official specifications for commercial fertilizers according to the Fertilizer Control Law¹⁾.

An example of simultaneous analysis of the nitrous acid and ammonium thiocyanate content of fertilizer by HPLC is described. Analysis was performed in conformance with the test method that appears in Testing Methods for Fertilizers (5.8.b and 5.9.a, 2016)²⁾, published by the Food and Agricultural Materials Inspection Center (FAMIC).

■ Analysis of a Standard Mixture

Analytical conditions are shown in Table 1. A chromatogram of a standard mixture of nitrous acid and ammonium thiocyanate (20 mg/L each) is also shown in Fig. 1. Please refer to the test method²⁾ for details on mobile phase preparation. Standard solution was prepared by dissolution and dilution with water. An NH2 column was used for analysis.

Table 1 Analytical Conditions

System	: Prominence
Column	: Shodex Asahipak NH2P-50 4E (250 mm L. × 4.6 mm I.D., 5 μm)
Guard Column	: Shodex Asahipak NH2P-50G 4A (10 mm L. × 4.0 mm I.D., 5 μm)
Mobile Phase	: (Sodium) phosphate buffer containing sodium perchlorate
Flowrate	: 1.0 mL/min
Column Temp.	: 40 °C
Detection	: UV-VIS detector (SPD-20AV) at 210 nm
Injection Vol.	: 10 μL




Fig. 1 Chromatogram of Standard Mixture (20 mg/L each)

■ Linearity of Calibration Curves

Fig. 2 shows calibration curves for nitrous acid and ammonium thiocyanate analyzed under the conditions shown in Table 1. The range used for calibration curves was 1 to 20 mg/L. Good linearity was obtained for both compounds with contribution rates (R^2) of 0.9999 or higher.




Fig. 2 Linearity of Calibration Curves
I) Nitrous acid (1 - 20 mg/L),
II) Ammonium thiocyanate (1 - 20 mg/L)

■ Repeatability

The relative standard deviation (%RSD) of retention times and peak areas obtained from an analysis of each compound at 0.1 mg/L repeated six times is shown in Table 2 and 3. The concentration analyzed (0.1 mg/L) was equivalent to 1/10 the lowest concentration on the calibration curve of either compound. Good repeatability was obtained for the retention times and peak areas of both compounds.

Table 2 Repeatability of Retention Time and Peak Area for Nitrous Acid Analysis

	R.T. (min)	Area
1st	6.452	6,752
2nd	6.452	6,801
3rd	6.450	6,722
4th	6.452	6,794
5th	6.452	6,727
6th	6.451	6,823
Ave.	6.452	6,770
%RSD	0.012	0.62

Table 3 Repeatability of Retention Time and Peak Area for Ammonium Thiocyanate Analysis

	R.T. (min)	Area
1st	16.868	2,551
2nd	16.870	2,534
3rd	16.882	2,524
4th	16.881	2,519
5th	16.885	2,546
6th	16.868	2,553
Ave.	16.876	2,538
%RSD	0.046	0.57

■ Analysis of Nitrous Acid and Ammonium Thiocyanate in Fertilizers

The method used to prepare fertilizer samples is shown in Fig. 3. The method of pretreatment differed depending on whether the sample was a powder or liquid, and on extraction liquid pH. Standard additions of nitrous acid and ammonium thiocyanate were made to the samples for analysis (fertilizer) before being further prepared according to the procedure shown in Fig. 3^{3), 4)}. The results of analysis are shown in Fig. 4.

Testing Methods for Fertilizers (2016) describes a method that uses ion chromatography for the analysis of ammonium thiocyanate (5.8.a). However, because this method produces a complex eluent, it takes some time for the baseline to stabilize. Another problem with the method is that it has resulted in multiple system peaks and peaks close to the elution position of ammonium thiocyanate. Nevertheless, on this occasion, the results show that good separation was achieved, including for contaminating constituents in the fertilizer.




Fig. 3 Pretreatment




Fig. 4 Chromatograms of Fertilizer Sample Solutions
(I) Powdered fertilizer A, (II) Powdered fertilizer B,
(III) Liquid fertilizer A, (IV) Liquid fertilizer B

References

- 1) Ministry of Agriculture, Forestry and Fisheries notification: Establishment of official specifications for commercial fertilizers according to the Fertilizer Control Law, dated February 22, 1986. Ministry of Agriculture, Forestry and Fisheries notification no. 284: Final revision, dated March 30, 2016. Ministry of Agriculture, Forestry and Fisheries notification no. 884: Enactment, dated April 1, 2016. [In Japanese]
- 2) Testing Methods for Fertilizers (2016), Food and Agricultural Materials Inspection Center (FAMIC): http://www.famic.go.jp/ffis/fert/obj/shikhenho_2016.pdf [In Japanese]
- 3) Masanori Hasegawa, Yasuharu Kimura: Determination of Nitrous Acid and Ammonium Thiocyanate in Fertilizer by High Performance Liquid Chromatography (HPLC): A Collaborative Study, Research Report of Fertilizer, vol. 8, 70-78 (2015) [in Japanese]
- 4) Kohei Ito, Yasuharu Kimura, Masanori Hasegawa, Yuji Shirai: Simultaneous determination of nitrous acid and thiocyanate in fertilizer by HPLC, Japanese Journal of Soil Science and Plant Nutrition, 87 (3), 120-124 (2016) [In Japanese]

First Edition: Jan. 2017



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2017

Shimadzu Corporation

www.shimadzu.com/an/

Application News

No.L514

High Performance Liquid Chromatography

Analysis of Melamine and Its Related Substances in Fertilizers

Calcium cyanamide is effective as a fertilizer, pesticide, soil amendment, and for many other uses, and it is essential compound for producing high quality vegetables. Recently, high levels of melamine were discovered as a byproduct in some calcium cyanamide hydrate products, pelletized by adding water to calcium cyanamide. Due to the risk of agricultural products absorbing the melamine from the soil, it has been identified as a potential public health risk. For example, if both melamine and its related substance cyanuric acid are ingested at the same time, they can form crystals that can impede kidney function.¹⁾

As a result, the Food Safety and Consumer Affairs Bureau in the Ministry of Agriculture, Forestry and Fisheries in Japan issued a notice specifying a 0.4 % provisional maximum allowable concentration of melamine in calcium cyanamide.²⁾

This article describes an example of pretreating and analyzing melamine and its related substances, namely ammeline, ammelide, and cyanuric acid, in fertilizer, in accordance with the testing methods supervised by the Food and Agricultural Materials Inspection Center (FAMIC) in Japan for fertilizers and other substances (2016, 8.1.c).^{3), 4), 5)}

■ Analysis of Standard Solution

The structure of melamine and its related substances is shown in Fig. 1. The analytical conditions are indicated in Table 1. The chromatogram of the standard mixture solution of melamine and its related substances (1 mg/L each) is shown in Fig. 2. For more details regarding the procedures used to prepare the standard solution and mobile phase, refer to the applicable test methods. Calibration curves for melamine and its related substances are shown in Fig. 3. Calibration curves were prepared for a concentration range of 0.05 to 5 mg/L. The results indicated good linearity, with a contribution rate (R^2) over 0.9999.




Fig. 1 Chemical Structure of Melamine and Its Related Substances

Table 1 Analytical Conditions

System	: Prominence
Column	: TOSOH, TSKgel Amide-80 (250 mm L. × 4.6 mm I.D., 5 μm)
Guard Column	: TOSOH, TSKgel guardgel Amide-80 (15 mm L. × 3.2 mm I.D.)
Mobile Phase	: (Sodium) phosphate buffer pH 6.7±0.2 / Acetonitrile = 1/4 (v/v)
Flowrate	: 1.0 mL/min
Column Temp.	: 40 °C
Injection Vol.	: 10 μL
Detection	: UV-VIS detector (SPD-20A) at 214 nm




Fig. 2 Chromatogram of Standard Mixture (1 mg/L each)




Fig. 3 Linearity (0.05 to 5 mg/L)

■ Repeatability

The relative standard deviation (%RSD) results for peak area after analyzing the standard solution (0.1 mg/L) six consecutive times were very good, with 0.41 % for cyanuric acid, 0.42 % for ammelide, 0.52 % for melamine, and 0.56 % for ammeline. When pretreated as indicated in Fig. 4, the 0.1 mg/L concentration of the standard solution is equivalent to a 0.02 % concentration of melamine and other related substances in fertilizer.

■ Analysis of Melamine and its Related Substances in Fertilizer

The analytical sample (fertilizer) was pretreated in accordance with the test method by adding a standard quantity of melamine and its related substances. The pretreatment procedure is indicated in Fig. 4 and the analytical results are shown in Fig. 5. In this example, five types of samples were tested, including two types of commercially available nitrolime, a synthetic fertilizer that contains calcium cyanamide, a synthetic fertilizer, and ammonium sulfate. The quantities of the substances added to the analytical samples, as a percentage of mass, were equivalent to about 0.035 to 2.8 % melamine, about 0.035 to 1.6 % ammelide, about 0.035 to 1.1 % ammelide, and about 0.037 to 1.2 % cyanuric acid. These results demonstrate that the Prominence system provides more than adequate performance for measuring the provisional 0.4 % melamine limit issued by the Food Safety and Consumer Affairs Bureau, Ministry of Agriculture, Forestry and Fisheries for calcium cyanamide and fertilizers that contain calcium cyanamide as an ingredient.




Fig. 4 Pretreatment Procedure




Fig. 5 Chromatograms for (I) Nitrolime A, (II) Nitrolime B, (III) Synthetic Fertilizer Containing Calcium Cyanamide, (IV) Synthetic Fertilizer, and (V) Ammonium Sulfate

References

- 1) Health effects of melamine, etc.: Food Safety Commission of Japan, October 9, 2008, updated April 30, 2009
- 2) Setting Allowable Limit Values for Melamine in Calcium Cyanamide, Food Safety and Consumer Affairs Bureau, Ministry of Agriculture, Forestry and Fisheries in Japan, Notice No. 6116, 2012, issued March 25, 2013 and partially revised March 30, 2013
- 3) Testing Methods for Fertilizers (2016), Food and Agricultural Materials Inspection Center (FAMIC): http://www.famic.go.jp/ffis/fert/obj/shikhenho_2016.pdf (in Japanese)
- 4) Etsuko Bando and Yuji Shirai: Validation of High Performance Liquid Chromatography (HPLC) for Determination of Melamine and Its Related Substances in Fertilizer, Research Report of Fertilizer Vol. 6, pp. 27 to 35 (2013)
- 5) Etsuko Bando and Shigehiro Kai: Determination of Melamine and Its Related Substances in Fertilizer by High Performance Liquid Chromatography (HPLC): A Collaborative Study, Research Report of Fertilizer Vol. 7 10-21 (2014)

First Edition: Jan. 2017



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2017

Application News

No.C140

Liquid Chromatography Mass Spectrometry

Ultra-Sensitive and Rapid Assay of Neonicotinoids, Fipronil and Some Metabolites in Honey by UHPLC-MS/MS [LCMS-8060]

Neonicotinoids are a class of insecticides widely used to protect fields as well as fruits and vegetables.

Recently the use of these compounds became very controversial as they were pointed as one cause of the honeybees colony collapse disorder. Since pollination is essential for agriculture, extensive studies have been conducted to evaluate the impact of neonicotinoids on bee health. Following this the European Food Security Authority (EFSA) limited the use of thiamethoxam, clothianidin and imidacloprid. Fipronil, a pesticide from a different chemical class, has been also banned by EFSA for maize seed treatment due to its high risk for honeybee health.

In order to better understand the effect of these compounds on bees and their contamination in pollen and honey, a highly sensitive assay method was necessary. A method was set up using Nexera X2 with LCMS-8060.

Sample Preparation

Thiamethoxam-d3, imidacloprid-d4 and clothianidin-d3 were used as internal standards.

Compound extraction was performed using a QuEChERS (Quick, Easy, Cheap, Effective, Rugged and Safe) method with an additional dispersive Solid Phase Extraction (dSPE) step.

5 g of honey (± 1 %) were weighted in a 50 mL polypropylene tube. 5 μ L of internal standard solution at 5 μ g/mL of each compound in acetonitrile was added on honey and let dry for 10 minutes. 10 mL of ultra pure water were added and the samples were homogenized by vortex mixing for 1 minute. 10 mL of acetonitrile were then added followed by vortex mixing for 1 minute.

After incubation at room temperature for one hour with gentle shaking, a commercially available salt mix from Biotage (4 g MgSO₄, 1 g Sodium Citrate, 0.5 g Sodium Citrate sesquihydrate, 1g NaCl) was added. After manual shaking, samples were centrifuged at 3000 g for 5 minutes at 10 °C.

Supernatant (6 mL) was transferred into a 15 mL tube containing 1200 mg of MgSO₄, 400 mg PSA and 400 mg C18 from Biotage. After centrifugation at 3000 g and 10 °C for 5 minutes the supernatant was transferred into a LCMS certified inert glass vial for analysis (Shimadzu LabTotal 227-34001-01).

Recovery

An "all-flowers" honey from the local supermarket was extracted with or without spike at 50 ppt. A blank extract (no honey) was prepared to evaluate losses or non specific interactions. Results are presented in Table 1.

Calculated recoveries are within acceptance values 70-120 % from EU SANTE/11945/2015.

Table 1 Measured Recoveries in Honey

Compound	Recovery	Compound	Recovery
Acetamiprid	78.8 %	Fipronil sulfone	74.2 %
Acetamiprid-N-desmethyl	93.4 %	Imidacloprid	83.2 %
Chlothianidin	70.6 %	Nitenpyram	87.0 %
Dinotefuran	76.5 %	Thiacloprid	82.2 %
Fipronil	78.1 %	Thiamethoxam	75.6 %




Fig. 1 Chromatogram of the Target Compounds at Their Lower Limit of Quantification

Table 2 Analytical Conditions

System	: Nexera X2	System	: LCMS-8060
Column	: ACE SuperC18 (100 mm L. × 2.1 mm I.D., 2 µm)	Ionization	: Heated ESI
Column Temperature	: 30 °C	Probe Voltage	: +1 kV (positive ionization) / -1.5 kV (negative ionization)
Mobile Phases	: A: Water = 0.05 % ammonia B: Methanol + 0.05 % ammonia	Temperature	: Interface: 400 °C Desolvation Line: 200 °C Heater Block: 400 °C
Flowrate	: 600 µL/min	Gas Flow	: Nebulizing Gas: 3 L/min Heating Gas: 10 L/min Drying Gas: 5 L/min
Gradient	: 5 %B to 100 %B in 3 min 100 %B to 5 %B in 0.1 min		
Total Run Time	: 4 min		
Injection Volume	: 2 µL (POISe mode with 10 µL of water)		

Table 3 MS/MS Acquisition Parameters

MRM Transitions	Name	Polarity	MRM Quan	MRM Qual	ISTD
	Acetamidrid	+	223.1 > 126.0	223.1 > 56.1	2
	Acetamidrid-N-desmethyl	+	209.1 > 126.0	211.1 > 128.0	2
	Clothianidin	+	250.1 > 169.1	250.1 > 132.0	3
	Dinotefuran	+	203.0 > 114.0	203.0 > 87.0	1
	Fipronil	-	435.0 > 330.0	435.0 > 250.0	3
	Fipronil sulfone	-	451.0 > 415.0	451.0 > 282.0	3
	Imidacloprid	+	256.1 > 175.1	258.1 > 211.1	2
	Nitenpyram	+	271.0 > 126.0	271.0 > 225.0	3
	Thiacloprid	+	253.1 > 126	253.1 > 90.1	1
	Thiamethoxam	+	292.1 > 211.1	292.1 > 181.1	1
	Thiamethoxam-D3	+	295.1 > 214.05	---	1
	Imidacloprid-D4	+	260.1 > 179.1	---	2
	Clothianidin-D3	+	253.1 > 132.05	---	3
Dwell Time	: 3 to 34 msec depending upon the number of concomitant transitions to ensure to have at least 30 points per peak (max total loop time 140 msec).				
Pause Time	: 1 msec				
Quadrupole Resolution	: Q1: Unit Q3: Unit				

Calibration


Calibration curves were prepared in acetonitrile to obtain final concentrations ranging from 0.5 pg/mL (1 fg on column) to 5 ng/mL. These concentrations corresponds to 1 ng/kg and 10 µg/kg in honey, respectively.

For each compound, the lower limit of quantification was selected to give an accuracy between 80-120 % (see table 4).

A typical calibration curve is shown in Fig. 2.

Table 4 Limits of Quantification in Honey

Compound	LOQ (µg/kg)	Compound	LOQ (µg/kg)
Acetamidrid	0.005	Fipronil sulfone	0.001
Acetamidrid-N-desmethyl	0.005	Imidacloprid	0.020
Chlothianidin	0.020	Nitenpyram	0.020
Dinotefuran	0.010	Thiacloprid	0.005
Fipronil	0.001	Thiamethoxam	0.005

**Fig. 2 Calibration Curve of Acetamidrid**

■ **Real Samples Analysis**

Nine honey samples purchased at the local supermarket or used as raw materials in cosmetics (orange tree honey) were assayed as unknowns.

All tested honeys showed concentrations far below the authorized maximum residue limit. But thanks to the very high sensitivity reached, even low concentrations of neonicotinoids were quantified. Results are presented in table 5. A representative chromatogram of a sample honey is shown in Fig. 3.

Table 5 Honey Samples Results (concentrations in µg/kg)

Honey	Acetamidrid	Clothianidin	Imidacloprid	Thiacloprid	Thiamethoxam
1. Provence creamy	---	---	0.20	---	0.010
2. Italy creamy	0.15	---	0.17	---	---
3. Pyrenees liquid	0.38	---	0.043	0.020	---
4. French-Spanish creamy	0.27	---	0.047	0.020	---
5. Thyme liquid	---	---	---	---	---
6. Lemon tree creamy	1.7	---	0.15	0.033	---
7. Orange tree liquid	1.2	---	0.62	---	---
8. Flowers creamy	0.14	---	0.055	0.39	---
9. Flowers liquid	0.34	---	0.11	0.010	---

Honey	Dinotefuran	Nitenpyram	Acetamidrid-N-desmethyl	Fipronil	Fipronil sulfone
1. Provence creamy	---	0.052	0.005	---	---
2. Italy creamy	---	0.040	---	---	---
3. Pyrenees liquid	---	---	0.015	0.004	---
4. French-Spanish creamy	---	0.032	---	---	---
5. Thyme liquid	---	---	---	---	---
6. Lemon tree creamy	---	---	0.020	---	---
7. Orange tree liquid	---	0.024	0.018	---	---
8. Flowers creamy	---	---	0.016	---	---
9. Flowers liquid	---	---	0.006	---	---




Fig. 3 Chromatogram of a Sample Honey (Pyrenees)

■ Stability

The thyme honey sample with no detectable target compound was spiked at 50 ng/kg with all compounds prior to extraction. The extract obtained was then consecutively injected 150 times in the system.

The results presented in Fig. 4 show excellent stability of the signal even at these low concentrations. This demonstrates that the excellent sensitivity can be maintained over long series of real sample analysis thanks to the ion source ruggedness.




Fig. 4 Stability of Peak Areas in Real Honey Samples

■ Conclusion

A method for ultra sensitive assay of neonicotinoids in honey was set up. The sample preparation was simple but provided excellent recoveries. The injection mode used prevented the use of tedious evaporation/reconstitution or dilution steps.

Thanks to the high sensitivity obtained enabled assay in real samples at very low levels far under the regulated residue levels. Furthermore, even at low measured concentrations, the system demonstrated its stability after long analytical series of real samples.

This method can be a very efficient support tool to better understand the impact of neonicotinoids on honey bee colonies and could be easily transposed to pollen or bee samples.

First Edition: Dec. 2016



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2016

Application News

No. SCA_210_041

Liquid Chromatography Mass Spectrometry


Sensitive method for the determination of Fipronil and its metabolite Fipronil Sulfone in egg using QuEChERS sample pretreatment and LC-MS/MS detection [LCMS-8060]

■ Introduction

Fipronil concerns a broad-spectrum insecticide from the group of phenylpyrazoles used in many countries as a biocide and plant protection product against fleas, lice, ticks, cockroaches, mites and other insects. Fipronil is an active compound in veterinary products fighting tick and flea infestations in dogs and cats. The use as plant protection product is restricted to seed treatment in the European Union since 2007. However, due to the illegal use as addition to the cleaning supplies used in chicken coops the eggs and meat might get contaminated as well.

The MRL (maximum residue levels) for Fipronil and its metabolite Fipronil sulfone (which is classified as having similar toxicity) in eggs is set to 0.005 mg/kg by the EU (by definition the sum of fipronil and fipronil-sulfone expressed as fipronil) [1], so that there is an actual requirement for the determination of both compounds in egg matrix at a relatively low level.

This application news presents a simple method using a standard QuEChERS extraction protocol followed by LC-MS/MS detection.



Fipronil

MF $C_{12}H_4Cl_2F_6N_4OS$
MW 437,1 g/mol

■ Sample preparation

Compound extraction was performed using a simple QuEChERS (Quick, Easy, Cheap, Effective, Rugged and Safe) method.

5 g of egg (egg white and egg yolk) were weighted into a 50 mL polypropylene tube, diluted with 5 mL

of water and spiked with a respective amount of Fipronil, Fipronil sulfone and in addition Fipronil-desulfinyl and Fipronil-sulfide (neochema, Germany).

10 mL of acetonitrile was added and the samples were mixed vigorously. After that ready to use QuEChERS extraction salts (Q-sep™ Q110, Pouch and tubes – cat. #26235, Restek) were added for sample drying and buffering. Samples were mixed again and centrifuged at 4500 rpm for 5 minutes. 1 mL of the supernatant was transferred into a dSPE tube (Q-sep™ QuEChERS dSPE – cat. #26217, Restek), shaken for 2 minutes, centrifuged, the supernatant was transferred into a glass vial and the pH was adjusted with 5% formic acid solution in acetonitrile (10 µL/mL supernatant).

■ Materials and methods

Extracts were analyzed using a method set up with Shimadzu's LC/MS/MS Method Package for Residual Pesticides Version 2 and a Nexera X2 UHPLC system coupled to a LCMS-8060 mass spectrometer. Analysis was carried out using MRM (Multi Reaction Monitoring) mode.



LC system	Nexera X2 (Shimadzu, Japan)
Analytical column	Raptor Biphenyl™ 100 x 2.1 mm, 2.7 µm (RESTEK)
Column oven temperature	35 °C
Injection volume	2 µl
Mobile Phase A	2 mM ammonium formate + 0.002% formic acid - Water
Mobile Phase B	2 mM ammonium formate + 0.002% formic acid - Methanol

Mass spectrometer	LCMS-8060 (Shimadzu, Japan)
Interface voltage	-3 kV
Q1 resolution	Unit (0.7 Da FWHM)
Q3 resolution	Unit (0.7 Da FWHM)
Nebulizing gas flow	3 L/min
Drying gas flow	10 L/min
Heating gas flow	10 L/min
DL temperature	150 °C
Heat block temperature	300 °C
Interface Temperature	350 °C

In addition, the so-called "MRM spectrum mode" was used for analysis. Here, not only the fragments of the quantifier and the qualifiers are determined, but also a higher number (typically 6-10) of MRM fragment ions. Using this MRM spectrum mode, conventional MRM quantification is combined with a high-quality MRM product ion spectrum, which can be used in a library search routine, thus increasing the specificity and verification of results (Figures 1 and 2).




Figure 2: The figure shows MRM chromatograms for Fipronil, one recorded with the usual 2 fragment ions, and compared with a method with higher number (12) of fragment ions which, despite this fact, have the same sensitivity.




Figure 1: Result of the library search, presented with LabSolutions Insight Screening software

Calibration

The matrix matched calibration curve (Figure 3-6) was prepared according to the method described before ranging from 0.0005 mg/kg to 0.05 mg/kg. Control samples at 0.001 mg/kg and 0.01 mg/kg correspond to the calibration curve.




Figure 3: Calibration curve of Fipronil in egg ranging from 0.0005 mg/kg to 0.05 mg/kg




Figure 4: Calibration curve of Fipronil-sulfide in egg ranging from 0.0005 mg/kg to 0.05 mg/kg

[1] EU Commission Regulation No 1127/2014 of 20 October 2014 Amending Annexes II and III to Regulation (EC) No 396/2005 of the European Parliament and of the Council in regards to maximum residual levels for amitrole, dinocap, fipronil, flufenacet, pendimethalin, propyzamide and pyridate in or on certain products.




Figure 5: Calibration curve of Fipronil-sulfone in egg ranging from 0.0005 mg/kg to 0.05 mg/kg




Figure 6: Calibration curve of Fipronil-desulfinyl in egg ranging from 0.0005 mg/kg to 0.05 mg/kg

Conclusion

By using the LC/MS/MS method package for residual pesticides V2 and a QuEChERS sample preparation a method for the determination of Fipronil and Fipronil-sulfone in eggs below the requested MRL of 0.005 mg/kg could be set up rapidly without further method development.



Shimadzu Europa GmbH

www.shimadzu.eu

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

Application News

No.C141

Liquid Chromatography Mass Spectrometry

High Sensitivity Analysis of Peanut Allergen in Cumin and Spice Mix [LCMS-8060]

Food allergens are a major public health concern. Among them, peanut allergy is one of the common food allergies. To avoid unexpected contact with food allergens, food labels are strictly used to indicate the presence of specific allergens. With the increasing awareness of food allergies, the presence of undeclared peanut in cumin lead to huge recalls in recent years. Although ELISA is the most commonly used technique to detect allergens, its false-positive rate is a major concern due to its cross-reactivity. We developed a method with high specificity and sensitivity to overcome this issue by using a high sensitivity triple quadrupole mass spectrometer to detect peanut allergen Ara h1 (Fig.1) in commercially available spices and seasonings.




Fig. 1 Structure of Ara h1 [357I] (68kDa) Vicilin Like Protein

■ Sample Preparation

Commercially available defatted peanut flour was purchased and used for the initial development work. The test samples were ground and protein content was enriched by liquid-liquid extraction. Extracted proteins were denatured, reduced and alkylated before subjecting to tryptic digestion to obtain peptides that were quantitated as proxies of original protein abundance.

Cinnamon, cumin, chilli pepper, ginger, garlic, mustard seed, nutmeg, oregano, rosemary, sage, turmeric and thyme were selected as test food samples for evaluating cross-reactivity and sensitivity of the developed method. Food samples were pretreated as above with or without 2 ppm peanut powder.

■ Selection of MRM Transitions Using Skyline

Ara h1 is known as is known as the sensitizing allergen in 95 % of peanut allergy. Tryptic digest of protein extracted from peanuts were analyzed by monitoring theoretically calculated transitions of peptides based on amino acid sequences of two clones P17 and P41B of Ara h1.

MRM transitions for each clone was determined by using Skyline (MacCoss Lab Software). The transition list, which contained more than ten peptides for each clone, was reviewed by removing several peptides that could be susceptible by post translational modification and Maillard reaction during food processing. Finally, nine peptides including three common peptides to both clones were selected based on sensitivity. Three transitions were set for each peptide.




Fig. 2 Workflow of MRM Transition Optimization Using Skyline




Fig. 3 AA Sequences of P17/P41B and Nine MRM Chromatograms

Table 1 Analytical Conditions

System	: Nexera X2	System	: LCMS-8060
Column	: Shim-pack XR-ODS II (50 mm L. x 2 mm I.D., 1.6 µm)	Ionization	: Heated ESI
Column Temperature	: 40 °C	Probe Voltage	: +1 kV (positive ionization)
Mobile Phases	: A: Water + 0.1 % formic acid B: Acetonitrile	Temperature	: Interface: 250 °C Desolvation Line: 150 °C Heater Block: 200 °C
Flowrate	: 500 µL/min	Gas Flow	: Nebulizing Gas: 3 L/min Heating Gas: 20 L/min Drying Gas: 5 L/min
Gradient	: 2 %B (0.00 min) > 25 %B (7.00 min) > 95 %B (7.10-8.00 min) > 2 %B (8.10-10.00 min)		
Injection Volume	: 10 µL		

Table 2 MS/MS Acquisition Parameters

MRM Transitions	Name	Polarity	Quan	Qual1	Qual2
	EGEQEWGTPGSEVR	+	780.85 > 802.40	780.85 > 644.35	780.85 > 316.10
	NNPFYFPSR	+	571.25 > 669.35	571.25 > 506.25	571.25 > 229.10
	IPSGFISYILNR	+	690.40 > 765.45	690.40 > 211.15	690.40 > 502.25
	SSDNEGVIVK	+	524.25 > 515.35	524.25 > 359.25	524.25 > 175.05
	GSEEDITNPINLR	+	793.90 > 726.45	793.90 > 612.40	793.90 > 402.25
	GTGNLELVAVR	+	564.80 > 686.40	564.80 > 557.40	564.80 > 444.30
	EGEQEWGTPGSHVR	+	784.85 > 652.35	784.85 > 555.30	784.85 > 316.10
	SSENNEGVIVK	+	588.30 > 515.35	588.30 > 359.25	588.30 > 246.20
	GSEEEGDITNPINLR	+	822.40 > 726.45	822.40 > 612.40	822.40 > 402.25
Dwell Time	: 41 to 130 msec depending upon the number of concomitant transitions to ensure to have at least 15 points per peak (max total loop time 400 msec).				
Pause Time	: 3 msec				
CID Pressure	: 300 kPa				
Quadrupole Resolution	: Q1: Unit Q3: Unit				

Interface Optimization

Ionization parameters optimization was performed using companion software ISSS (Interface Setting Support Software, Shimadzu Corp.). As a result, sensitivity was improved more than twofold compared to default values.




Fig. 4 Interface Optimization Results

Effect of Surfactant During Digestion

A higher intensity of peptides by addition of a surfactant during tryptic digestion was expected due to improved digestion efficiency. However, the intensity of peptides were relatively worse by adding surfactant. Thus, no surfactant was used for tryptic digestion.




Fig. 5 Difference of the Chromatograms of Peptide GTG... by Addition of Surfactant

■ Peanut Allergen in Other Nuts

Walnuts, cashew nuts, and almonds were analyzed to test specificity. These nuts were spiked with 2 ppm (2 mg/kg) of peanut before sample preparation. The spiked peanut peptides were successfully detected and any obvious peak was detected in blank samples.




Fig. 6 Chromatograms of Peptide GTG... in Other Kind of Nuts With or Without Spiking with Peanuts

■ Detection of ARA h1 in Spice Mixes and Seasonings

Several spice mixes and seasonings were analyzed using sample preparation and analytical conditions described here. Peaks of tryptic peptides of Ara h1 from samples without spiking of peanut peptides were detected.




Fig. 7 Detected Peaks of Peptide GTG... in Chili Mix and Seasoning

■ Peanut Allergen in Spices

Contaminated spice samples were prepared and analyzed to confirm that the low amount of peanuts added into the various spices can be detected. Peptides of Ara h1 were successfully observed from the spice samples spiked with 2 ppm of peanuts. It was also confirmed that there are no obvious false-positive peaks from the blank samples.




Fig. 8 Chromatograms of Peptide GTG... in Spices With or Without Spiking with Peanuts

■ Conclusion

A method for the analysis of Ara h1 in spices and seasonings was successfully developed. The combination of the developed method and a high sensitivity triple quadrupole mass spectrometer enabled the detection of 2 ppm or lower of peanut allergen Ara h1 in spices and seasonings.

First Edition: Dec. 2016



For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "®". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2016

Shimadzu Corporation
www.shimadzu.com/an/

Application News

No. AD-0095

Food Additives / Nexera X2

Fast and High Sensitivity Analysis of Six Preservatives in Beverages by UHPLC with Photodiode Array Detection

■ Introduction

Food preservatives are additives to inhibit, retard or prevent mould, acidification or other deterioration of foodstuffs caused by microbial contamination. The most commonly used preservatives in beverages are benzoic acid, sorbic acid and four parahydroxybenzoic acid esters (parabens). However, excess amounts of these additives can be harmful to consumer health. In this regard, the minimum permissible concentrations of preservatives are regulated in most countries to ensure safety for consumer *¹. Therefore, quantitative analysis of these preservatives in food is not only required for food quality assurance but also important for consumer interest and protection. High performance liquid chromatography (HPLC) has been used for analysis of the preservatives in beverage *^{2, *3, *4}. In this Application News, a new rapid and high sensitivity UHPLC method for simultaneous determination of the six preservatives in beverages is described. A gradient elution was optimized for separation and quantitation of the six preservatives with a photodiode array detector. A capillary flow cell with extra long optical path of 85 mm was employed to achieve high sensitivity for a very small injection amount of sample (1 µL) which was not cleaned up except filtration.

■ Experimental

Preparation of Standards and Samples

Benzoic acid, sorbic acid and parabens were obtained from chemicals suppliers. A mixed stock solution of 1.0 g/L of benzoic acid, sorbic acid and methyl, ethyl, propyl, butyl parabens were prepared with ethanol/water (70/30) solvent as the diluent. A set of nine working standards was prepared from the stock solution using the same diluent at the concentrations shown in Table 1. Soft drink, mango juice and cocoa drink were purchased at the local supermarket. The soft drink and mango juice were diluted 20 times and 2 times with diluent respectively while cocoa drink was not diluted. All the samples were filtered through a 0.45 µm syringe filter prior to injection to UHPLC.

Table 1 Concentrations of Working Standards of Six Preservatives for Setting Calibration Curves

No.	Working Standard	Benzoic acid (mg/L)	Sorbic acid (mg/L)	Parabens (mg/L)
1	S1	0.2	0.008	0.01
2	S2	2.0	0.08	0.1
3	S3	4.0	0.16	0.2
4	S4	20.0	0.8	1.0
5	S5	60.0	2.4	3.0
6	S6	80.0	3.2	4.0
7	S7	100.0	4.0	5.0
8	S8	150.0	6.0	7.5
9	S9	200.0	8.0	10.0

Instrumental and Analytical Conditions

A Nexera X2 UHPLC system (Shimadzu Corporation, Japan) was used in this work. The system is consisted of a high pressure binary gradient solvent delivery unit (LC-30AD pumps) and an UHPLC autosampler (SIL-30AC) coupled to a photodiode array detector (SPD-M30A) with a high sensitivity capillary flow cell (85 mm optical path length) featured as total reflection and low dispersion. A YMC Triart C18 column of 1.9 µm particle size (150 mmL. × 2.0 mm I.D.) was used for the separation of preservatives (benzoic acid, sorbic acid and methyl, ethyl, propyl, butyl parabens) with an optimized linear gradient program developed. The details of the LC conditions are shown in Table 2.

Table 2 Analytical Conditions of Preservatives in Beverages on Nexera X2 UHPLC

Column	: YMC Triart C18 1.9 µm 150 × 2.0 mm I.D.
Flow Rate	: 0.45 mL/min
Mobile Phase	: A: 1.5 % acetic acid + 1.5 % ammonium acetate in H ₂ O B: 1.5 % ammonium acetate in MeOH
Elution Mode	: Gradient elution: 40 % B (0.01 to 4.0 min) → 80 % B (4.01 to 5.5 min) → 40 % B (5.51 to 8.5 min)
Oven Temp.	: 45 °C
Injection Volume	: 1 µL
Detection (PDA)	: Wavelength 240 to 600 nm; Ref: 720 nm Quant, 240 nm for benzoic acid, 260 nm for other compounds

Results and Discussion

Method Development

The six preservatives were well-separated as sharp peaks between 1.7 min and 5.1 min as shown in Fig. 1. The total run time of the UHPLC method is 8.5 mins, which is several times faster than the HPLC method reported *2, *3, *4. It is worth to note that two wavelengths were selected for quantitative data processing, i.e., 240 nm for benzoic acid and 260 nm for the rest five compounds *4.




Fig. 1 Chromatograms of Mixed Standard (S3) with 1 µL Injection Volume on Nexera X2




Fig. 2 Calibration Curves of Six Preservatives at Concentration S1 to S9 (see Table 1)

Fig. 2 shows the calibration curves of the six compounds established with 1 µL injection volume. The linearity with correlation coefficient (R^2) greater than 0.999 across the wide calibration range of 0.008 to 200 mg/L was obtained for the six compounds.

The repeatability of the method was evaluated at the levels S2 and S5. The peak area %RSD for the six compounds were lower than 5.1 % and 0.3 % respectively (Table 3).

Table 3 Results of Repeatability Evaluation Using Working Standard S2 and S5 (n=6, 1µL Injection)

Compound	Conc. (mg/L)	RSD%	Conc. (mg/L)	RSD%
Benzoic acid	2.0	1.1	60.0	0.2
Sorbic acid	0.08	1.5	60.0	0.2
Methyl paraben	0.1	1.2	2.4	0.2
Ethyl paraben	0.1	3.8	3.0	0.2
Propyl paraben	0.1	3.2	3.0	0.2
Butyl paraben	0.1	5.1	3.0	0.3

The LOD and LOQ of the method, and peak identification criteria (RT & λ_{Max}) are summarized in Table 4. The results were obtained from the mixed standard S1 (Fig. 3). The high sensitivity achieved, i.e., LOQs ranging at 8 to 10 µg/L of the compounds except benzoic acid (280 µg/L), is attributed partially to the use of a high sensitivity SPDM30A detector with using a capillary cell of 85 mm optical path.

Table 4 LOD (S/N=3), LOQ (S/N=10) and Peak Identification Criteria of UHPLC Method Obtained from S1 Chromatogram

Compound	Conc. (µg/L)	RT	λ_{Max}	LOD (µg/L)	LOQ (µg/L)
Benzoic acid	200	1.702	238	90	280
Sorbic acid	8	2.183	257	2.7	8
Methyl paraben	10	2.866	258	3	10
Ethyl paraben	10	3.687	257	3	10
Propyl paraben	10	4.445	258	3	10
Butyl paraben	10	5.091	256	3	10




Fig. 3 Chromatogram of Mixed Standard S1 (1 µL)

Analysis of Beverage Samples

The UHPLC method established was applied for quantitation of preservatives in three kinds of beverages: soft drink B1, fruit juice B2 and cocoa drink B3. The chromatograms of the samples are shown in Fig. 4 and the results are summarized into Table 5. No preservatives was detected in cocoa drink. Benzoic acid and sorbic acid were detected in the soft drink and fruit juice. The identification of both benzoic acid and sorbic acid peaks were confirmed by UV spectra.




Fig. 4 Chromatogram of Beverage Samples (Injection: 1 mL):
Soft Drink with 20 Times Dilution (top);
Fruit Drink with 2 Times Dilutions;
Cocoa Drink without Dilution (Bottom)

Table 5 Quantitative Results of Six Preservatives in Three Beverages, Each with Duplicate Injections

Sample Name	Benzoicacid		Sorbicacid		Parabens	
	RT (min)	Conc. (mg/L)	RT (min)	Conc. (mg/L)	RT (min)	Conc. (mg/L)
B1	1.7	82.4	2.18	137.4	ND	
B2	1.7	142.4	2.18	13.8		
B3	ND					

Conclusions

A rapid and high sensitivity UHPLC method for quantitation of six preservatives, benzoic acid, sorbic acid and four para-hydroxybenzoic acid esters (parabens), in beverages was established using a reversed phase UHPLC column (1.9 μ m particle size). A capillary flow cell with extra long optical path of 85 mm was employed in the photodiode array detector. The method achieves LOQs ranging 8 to 10 μ g/L for the compounds except benzoic acid (280 μ g/L), with 1 μ L injection volume. The very small injection volume minimizes the contamination of beverage samples to the column and system, as such suitable for direct analysis of beverage samples without need for clean up procedure.

<References>

- *1: Singapore Agri-Food and Veterinary Services website: Food Regulations – AVA (2014)
- *2: United States Department of Agriculture Food Safety and Inspection Service, Office of Public Health Science: *Determination of Benzoic acid, Sorbic acid and Methyl, Ethyl, Propyl and Butyl Parabens by HPLC* (2004)
- *3: Tzu-Yun Chu, Chine-Lin Chen and Hsueh-Fang Wang, *A Rapid Method for the Simultaneous Determination of Preservatives in Soy Sauce*, Journal of Food and Drug Analysis, Vol. 11, No. 3, Pages 246-250 (2003)
- *4: Shu-Xiang Sun, Shu-Bai Li and Pei Yao, *Simultaneous Determination of Six Preservatives in Beverage by HPLC*, 3rd International Conference on Material, Mechanical and Manufacturing Engineering (2015)

Application News

No. AD-0096

Food Safety Analysis / Nexera X2

Development of an UHPLC Method for Simultaneous Determination of Thirteen Bisphenols in Milk Samples

■ Introduction

Bisphenol A (BPA) and Bisphenol F (BPF) are the monomers to make polycarbonate plastics and epoxy resins. Their diglycidyl esters, i.e., BADGE and BFDGE (Fig. 1), are also present in the polymeric products. These materials are made into a variety of consumer products or used as inner coatings for baby milk bottles and reusable food containers which are allowed to use in refrigerator and microwave for food storage and heating. It has been reported that polycarbonate plastics and epoxy-based coatings can release BPA, BPF, BADGE and BFDGE as well as their reaction products as illustrated in Fig. 1^{*1}. These leached chemicals can migrate into food and become contaminants consumed by consumers eventually. Although researches indicate that the migration of these chemicals is normally extremely low^{*2}, the specific migration limits (SMLs) of bisphenols were

listed in the EU legislation No 1895/2005 on the restriction of use of certain epoxy derivatives in materials contacted with food. BPA has estrogenic effect and can disrupt normal hormone levels and development in fetuses and babies. In U.S., FDA has published food additive regulations prohibiting the use of BPA-based epoxy resins as inner coatings of containers for infant formula packaging^{*3}. We describe in this Application News a new UHPLC method for simultaneous determination of thirteen concerned bisphenols including BPA, BPF, BADGE, BFDGE and some structural analogues. An UHPLC system (Nexera X2, Shimadzu Corporation) with a high sensitivity fluorescence detector^{*4} was adopted to develop a fast and high sensitivity method to meet the requirements of regulations.




Fig. 1 Chemical Structures and Abbreviation Names of Bisphenol A, Bisphenol F, Their Diglycidyl Esters and Derivatives^{*1}

*: positional isomers of BFDGE-2 and BFDGE-2HCl-2 are not shown.

Experimental

Instrumental and Analytical Conditions

An UHPLC system - Nexera X2 (Shimadzu Corporation) equipped with a fluorescence detector (RF-20Axs) was employed in this work. Separation of bisphenol A (BPA), bisphenol F (BPF) and other 11 derivatives are performed using a Shim-pack HR-ODS column (250 × 3.0 mm, 3 μm) with an optimized gradient elution program. Pure water and acetonitrile (ACN) were used as UHPLC mobile phases without any additive. The detailed analytical conditions of the UHPLC method are shown in Table 1.

Standards and Spiked Milk Samples

A mixed standard stock solution of thirteen bisphenols (Refer to Table 2) containing BPA, BPF and other derivative compounds were prepared in ACN/H₂O (30:70). A serial of calibration standards of concentrations from 5 μg/L to 2,000 μg/L were prepared from the stock to set up multi-point calibration curves. Two blank milk matrix spiked with known concentrations of standards (100 and 1,000 μg/L) obtained from a third party laboratory were used for evaluation of the method performance.

Table 1 UHPLC Conditions of Bisphenols and Derivatives

Column	: Shim-pack HR-ODS (250 × 3.0 mm, 3 μm)
Mobile phase	: A: Water B: Acetonitrile
Elution program	: 0.1 min, 30 % B; 13 min, 45 % B; 37 min, 70 % B; 38 to 43 min, 85 % B; 43.1 min, 30 % B.
Flow rate	: 0.40 mL/min
Detection	: Ex 235 nm, Em 317 nm
Oven temp.	: 30 °C
Injection	: 10 μL

Results and Discussion

Development of Fast UHPLC Method

For well separation of the thirteen bisphenols studied, a reference HPLC method has a long running time of 95 mins. As shown in Fig. 3, the current UHPLC method was optimized to achieve fast elution for every compounds with sufficient separation resolution, especially for the separation of BADGE-H₂O and BADGE-H₂O-HCl at 19.4 and 19.9 mins. Due to the similarity in compound structure and chemical properties, separation of these two peaks was a main obstacle to achieve fast analysis speed. The results obtained show clearly the advantages of an UHPLC column with small particle size (3 μm) of the C18 stationary phase.

The main targets BPA (18.1 min) and BPF (13.4 min), BADGE (34.3 min) are separated completely without any inference. BFDGE has two positional isomers, which appeared as a pair at 28.8 min and 30.0 min, respectively.

Noted that, another pair of positional isomers BFDGE-2HCl appeared just before the BFDGE peaks at 26.9 min and 27.9 min, respectively.




Fig. 2 UHPLC-RF Chromatogram of Mixed Standards of Thirteen Bisphenols at Concentration of 100 μg/L Each

Table 2 Summary of UHPLC Method and Performance Evaluation Results for Analysis of Thirteen Bisphenols

ID #	Name	Ret Time (min)	Calibration range: 5-2000 μg/L		RSD (%), n=6		Sensitivity **	
			R ²	Accuracy (%) *	5 μg/L	100 μg/L	LOD (μg/L)	LOQ (μg/L)
1	BFDGE-2H ₂ O	7.9	0.9998	100.9	1.0	0.9	0.7	2.0
2	BADGE-2H ₂ O	9.5	0.9991	98.7	0.4	0.4	0.4	1.2
3	BPF	13.4	0.9997	101.7	0.7	0.3	0.8	2.5
4	BPA	18.1	0.9997	101.2	0.8	0.3	1.0	3.1
5	BADGE-H ₂ O-HCL	19.4	0.9998	101.5	0.5	0.3	0.5	1.5
6	BADGE-H ₂ O	19.9	0.9997	103.2	0.3	0.3	0.5	1.5
7	BFDGE-2HCL-1	26.9	0.9996	98.8	0.6	0.1	1.5	4.6
8	BFDGE-2HCL-2	27.8	0.9997	100.8	2.3	0.2	1.2	3.6
9	BFDGE-1	28.9	0.9997	99.3	0.9	0.3	1.4	4.2
10	BFDGE-2	30.0	0.9997	101.3	1.2	0.4	1.6	4.7
11	BADGE-2HCL	31.8	0.9997	98.9	1.8	0.3	0.6	1.7
12	BADGE-HCL	33.0	0.9997	101.7	0.5	0.4	0.7	2.2
13	BADGE	34.3	0.9997	100.7	1.0	0.5	0.5	1.6

* : Average of 12 concentration levels 5-2000 μg/L

** : Estimated using 5 μg/L mixed stds data based on S/N=3 for LOD and S/N=10 for LOQ

Calibration Curves, Range and Linearity

Linear calibration curves of the thirteen bisphenols are established using mixed standards samples for concentrations ranging from 5 µg/L to 2000 µg/L as shown in Fig. 3. A total of 12 concentration levels were used with each compound in the mixture being 5, 10, 20, 50, 100, 200, 400, 600, 800, 1200, 1500 and 2000 µg/L. All of the thirteen bisphenols peaks give excellent linearity with R^2 greater than 0.999 as tabulated in Table 2.




Fig. 3 Calibration Curves of the Thirteen Bisphenols with Concentration Range from 5 µg/L to 2,000 µg/L

Evaluation of Method Performance

The accuracy of the method at every calibration levels were calculated and the average accuracy values for every compounds are presented in Table 2. To evaluate repeatability of the method, six consecutive runs of the lowest concentration mixed standard sample (5 µg/L) and a mixed standard sample of 100 µg/L were performed.

The RSD values for the 5 µg/L mixed standards are less than 2.3 %, while RSD values for 100 µg/L concentration level are less than 1 %, as can be seen in Table 2. The limit of detection (LODs) and limit of quantification (LOQs) were determined from the chromatogram of the lowest concentration mixed standards (5 µg/L) as shown in Fig. 4, following the rule of $S/N=3$ for LOD and $S/N=10$ for LOQ. The obtained LODs and LOQs are at 0.4 to 1.6 µg/L, and 1.5 to 4.7 µg/L for the thirteen bisphenols (Table 2).




Fig. 4 UHPLC-RF Chromatogram of Mixed Standards of Thirteen Bisphenols, 5 µg/L each Compound

Spiked Milk Samples and Recovery

A blank milk sample and two spiked samples of 13 mixed standards in the same blank were obtained from an analytical laboratory, labeled as Blank, S1 (spiked 100 µg/L) and S2 (spiked 1,000 µg/L). The blank milk matrix was analyzed first and the result showed no any detection of the 13 bisphenols studied.

The chromatogram of spiked sample S1 is shown in Fig. 5. The quantitative results and recovery data of the 13 bisphenols in both samples are tabulated in Table 3. In both samples, the measured concentrations of BADGE-H₂O-HCl are higher than the expected levels with a recovery around 130 %. On the other hand, three compounds with longer retentions (peaks 11 to 13) exhibit much lower concentrations as expected and low recovery of about 40 % and 60 %.




Fig. 5 UHPLC-RF Chromatogram of Spiked Milk Sample S1 (100 µg/L). The Sample was Diluted with Water for Two Time Prior to Injection (10 µL)

Table 3 Analysis Results of Spiked Milk Samples for 13 Bisphenols Determined by the UHPLC-RF Method Established

ID #	Name	Ret Time (min)	Spiked S1		Spiked S2	
			Conc. (µg/L)	Recovery %	Conc. (µg/L)	Recovery %
1	BFDGE-2H ₂ O	7.9	115.0	115	1140.8	114.1
2	BADGE-2H ₂ O	9.5	114.5	114.5	1123.6	112.4
3	BPF	13.4	114.9	114.9	1152.4	115.2
4	BPA	18.1	114.1	114.1	1162.9	116.3
5	BADGE-H ₂ O-HCL	19.4	131.3	131.3	1284.9	128.5
6	BADGE-H ₂ O	19.9	91.1	91.1	983.1	98.3
7	BFDGE-2HCL-1	26.9	79.3	79.3	947.6	94.8
8	BFDGE-2HCL-2	27.8	90.5	90.5	983.0	98.3
9	BFDGE-1	28.9	68.6	68.6	823.8	82.4
10	BFDGE-2	30.0	77.2	77.2	824.2	82.4
11	BADGE-2HCL	31.8	47.3	47.3	668.9	66.9
12	BADGE-HCL	33.0	43.3	43.3	635.7	63.6
13	BADGE	34.3	37.4	37.4	563.7	56.4

Conclusions

An UHPLC method with using a high sensitivity fluorescence detector was developed for fast, well-separation and high sensitivity analysis of thirteen bisphenols, including the most concerned BPA, BPF, BADGE and BFDGE, in milk samples.

This new method shows high sensitivity to low µg/L levels, high accuracy and excellent repeatability. The method was applied to spiked milk samples and the results indicated the good feasibility, high sensitivity and reliability in simultaneous determination of thirteen bisphenols in milk samples.

<References>

- *1 : Jan Poustka, Lenka Dunovská, Jana Hajšlová, Kateřina Holadová and Ivana Poustková, Czech J. Food Sci., Vol. 25, No. 4: 221–229 (2007)
- *2 : Jana Poustková, Jaroslav Dobiš, Ingrid Steiner, Jan Poustka, Michal Voldřich, Eur Food Res Technol (2004) 219:534–539.
- *3 : United States Food and Drug Administration web site: Bisphenol A (BPA) for Use in Food Contact Applications 4. Application News No. L408, Shimadzu Corporation.

Find us on 



Linked 



First Edition: November, 2017



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.

© Shimadzu Corporation, 2017